COVIDIEN

Introducer Products

Multi-language Instructions For Use (IFU)

TABLE OF CONTENTS

Instructions for use:

English	1
French	4
Portuguese	7
Spanish	10

COVIDIEN

ENGLISH

Introducer Products

INSTRUCTIONS FOR USE

DESCRIPTION

Introducer Products are stand-alone components, standard kits, or convenience kits built to user specifications. These products facilitate access to a vein or artery, allowing various catheters to be passed through for local or central diagnostic or therapeutic interventions.

These instructions for use cover all families of Introducer products and associated accessories.

Contents of these sets have been sterilized with ethylene oxide gas. They are only for single patient, single-use applications.

INDICATIONS FOR USE

The introducers are indicated for use in arterial and venous procedures requiring percutaneous introduction of intravascular devices.

WARNINGS

1. Carefully read these instructions before using these products. When using sets incorporating other manufacturers' components, also read their current Instructions For Use.
2. Physician should be familiar with insertion techniques of introducers using a guidewire (Seldinger/Modified Seldinger).
3. **CARE SHOULD BE EXERCISED DURING INSERTION, USE OR REMOVAL OF THE DEVICE TO PREVENT ASPIRATION OF AIR INTO THE VASCULATURE.**
4. Withdrawal, pullback, or manipulation of the guidewire when resistance is met may cause guide wire damage, breakage, or embolization.
5. Rapid withdrawal of the catheter or dilator through the hemostasis valve may cause misalignment of the valve gasket assembly, causing bleedback through the valve. **SHOULD THIS OCCUR, RESETTING OF THE VALVE MAY BE ACCOMPLISHED BY GENTLE PRESSURE INTO THE VALVE WITH THE TIP OF A DILATOR OR CATHETER.**
6. Inspect the packaging and product for damage prior to use. Do not use if the package is damaged as sterility may be compromised.
7. Care should be taken that acetone solutions or isopropyl alcohol do not come in contact with introducer hub. Such substances could weaken the hub and cause possible leakage.
8. Tighten all connections prior to use **WITHOUT OVERTIGHTENING**. Overtightening can cause damage to the components.
9. Periodically check all connections for tightness.

10. This device is not intended to be used as a catheter for performing angiographic studies.
 11. Injection of dye through any component should be via hand injection only.
 12. This Device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and/or reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. This device is extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- 
13. For products containing the DEHP symbol on the label: This product contains Di(2-ethylhexyl)phthalate (DEHP) which has been shown to cause reproductive harm in male neonates, pregnant women carrying male fetuses, and peripubertal males. The following procedures have been identified as posing the greatest risk for DEHP exposure: exchange transfusion in neonates, total parenteral nutrition (TPN) in neonates (with lipids in polyvinylchloride (PVC) bag), multiple procedures in sick neonates (high cumulative exposure), heart transplantation or coronary artery bypass graft surgery (aggregate dose), and massive infusion of blood into trauma patient. It is recommended that DEHP-free medical products be considered when these procedures are to be performed on male neonates, pregnant women who are carrying male fetuses, and peripubertal males.

CAUTIONS

1. Federal Law (U.S.A.) restricts these devices to sale by or on the order of a physician.
2. Percutaneous introducer sheaths should not remain indwelling without internal support of the sheath wall.
3. Care should be exercised during insertion, manipulation and withdrawal of the catheter through the hemostasis valve. Do not insert blunt rigid objects through the hemostasis valve. The integrity of the hemostasis valve may be compromised with the use of excessive force, rapid withdrawal or rigid blunt accessories. If resistance is encountered, **DO NOT FORCE**.
4. If using the Tuohy-Borst Introducer, **DO NOT** tighten the adapter without internal support.
5. If repositioning of the sheath suture tab is necessary, ensure that internal support is provided by a vessel dilator prior to insertion to avoid twisting the sheath.
6. Upon removal of a catheter or other medical device, it is recommended that an obturator be immediately placed in the sheath.
7. Upon withdrawal of a device, the sideport should be aspirated in order to remove any clotted material that may have accumulated. This should be accomplished before the insertion of another medical device.
8. Depth markings on outside diameter of introducer tubing are references only for visual confirmation of introducer insertion distance. NOTE: Depth markings are not intended for use as a finite measurement.

PRECAUTIONS

1. Simultaneous advancement of sheath and dilator with a rotating motion is recommended to prevent sheath and vessel damage.
2. Maintain a firm grip on the guidewire at all times.
3. At no time should the introducer or guidewire be advanced or withdrawn when resistance is met without determining the cause by fluoroscopic examination or suitable imaging method.

PROCEDURES FOR USE

Percutaneous Introducers

The following instructions are supplied for information about the techniques of percutaneous catheter introduction. Procedural details may be altered at the discretion of the physician.

1. Prepare and drape the skin at the intended puncture site; achieve local anesthesia as necessary.

2. Prepare the introducer by assembling the dilator(s) into the introducer and locking the luer assemblies. Ensure all connections are tight without overtightening. Overtightening may result in damage to components.
3. If using the detachable Hemostasis Valve, attach proximal end of valve to sheath by twisting on the luer lock connector.
4. Locate vessel with an appropriate sized needle and syringe.
5. Enter vessel using thin walled needle or O-T-N catheter. Remove locator needle. If O-T-N catheter is used, advance catheter over needle into vessel and remove needle when good flow verifies placement of catheter tip within the vessel. **IF INSERTION IS UNSUCCESSFUL, WITHDRAW COMPLETE ASSEMBLY AS A UNIT. DO NOT ADVANCE NEEDLE INTO CATHETER OR ATTEMPT TO WITHDRAW CATHETER BACKWARDS OVER THE NEEDLE WHILE CATHETER IS IN THE PATIENT.**
6. Straighten "J" tip of guidewire with plastic insertion sleeve.
7. Insert guidewire into the O-T-N catheter or thin walled needle and gently advance it to the desired length. It may be necessary to gently rotate the "J" tip for successful advance of the guidewire. Avoid vigorous manipulation to prevent damage to vessel or shearing of the guidewire tip. **MAINTAIN A FIRM GRIP ON GUIDEWIRE AT ALL TIMES.**
8. If an obstruction is met that cannot be passed, remove needle and guidewire together and select another introduction site. **DO NOT ATTEMPT TO WITHDRAW GUIDEWIRE BACKWARDS THROUGH NEEDLE OR CATHETER AS THIS MAY RESULT IN SHEARING GUIDEWIRE OR DAMAGING CATHETER.**
9. When guidewire is advanced to desired location remove O-T-N catheter or needle proximally.
10. Thread assembled sheath/dilator until over guidewire and advance to puncture site. Enlarge puncture site with a small knick if needed, insert sheath/dilator unit into vessel with a slight rotary motion.
11. Remove the dilator and guidewire together leaving sheath in vessel.
12. Suture sheath in place using suture tab or hub. **DO NOT PLACE SUTURE ON SHEATH TUBING AS THIS MAY RESTRICT FLOW OR DAMAGE TUBING.**
13. Connect hemostasis valve securely to sheath when utilizing an introducer with removable valve.
14. A catheter one-half to one French size smaller than introducer sheath is recommended to allow adequate flow through introducer side-port.
15. Tuohy-Borst adapters should be gently hand-tightened to prevent blood reflux and/or catheter migration. Overtightening may compromise indwelling catheter lumen.
16. Care should be exercised not to pull the guidewire or catheter to extreme angles while the device is through the hemostasis valve, as it will distort the valve leaflets and promote leakage.
17. Follow hospital protocol for puncture site dressing and maintenance.

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Manufactured for Covidien, LLC