**ARGON MEDICAL DEVICES, INC.**  
 1445 Flat Creek Road, Athens, Texas 75751 USA  
 Tel: +1 800-927-4669;  
 Tel: +1 903-675-9321  
 www.argonmedical.com

EC	REP	<b>EMERGO EUROPE</b> Prinsessegracht 20 2514 AP The Hague The Netherlands +31 70 345 8570
----	-----	---


Figure 1


Figure 2

**Pro-Mag™ Ultra 2.2 Automatic Biopsy Instrument**

**Description:** The Pro-Mag™ Ultra 2.2 Automatic Biopsy Instrument is an automatic, reusable instrument for histological core biopsies. The Pro-Mag™ Ultra 2.2 Automatic Biopsy instrument has a throw (advancement) of 22mm and is designed to be used with disposable ACN™ needles.

**Indications for Use:** The Pro-Mag™ Ultra 2.2 Automatic Biopsy Instrument is used to obtain tissue core biopsy in the following clinical areas:

- **Urology** – Transrectal or transperineal biopsy of the prostate.
- **Radiology** – Percutaneous kidney, liver, lung and breast biopsy.

**Contraindications:** For use only for core biopsies of soft tissue as determined by a licensed physician. The instrument should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of core needle biopsy. Physician judgment is required when considering biopsy on patients with bleeding disorder, or receiving anti-coagulant medications.

**Note:** These instructions for the Pro-Mag™ Ultra 2.2 Automatic Biopsy Instrument are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

**Caution:** The Pro-Mag™ Ultra 2.2 Automatic Biopsy Instrument should be used by or under the supervision of a physician who is completely familiar with the use and limitations of this device and biopsy sampling.

**Caution:** The Pro-Mag™ Ultra 2.2 Instrument is designed to be compatible with ACN™ needles, sold separately. Compatibility of the instrument with other manufacturer's needle has not been tested.

Percutaneous introduction of the automatic cutting needle (ACN™ or SACN™) should be performed under ultrasound, computerized tomography (CT) or fluoroscopy in order to maximize needle placement.

Potential complications of biopsy may include, but not limited to: hemorrhage, hematoma, infection, hematuria, urinary retention, dysuria, pain, and perforation of the anterior wall of the bladder using transperineal approach.

**Directions For Use:**

After following the directions for reprocessing:

1. Inspect the disposable needle set package for damage. If undamaged, open the package using proper aseptic technique.
2. Remove the needle from the pouch. Inspect the stylet for damage to the point, the cannula for damage to the cutting edge, or other imperfections that would prevent proper operation of the needle assembly.
3. Open the lid of the biopsy instrument, load needle set into the biopsy instrument, remove the spacer, close the lid, and remove the protective sheath on the needle.
4. Safety button must be set in the position indicated by the arrow on the lid before cocking.
5. Cock the actuator handle twice in order to activate the device for firing (this automatically sets the safety.)
6. Release the safety by pushing the "Safety Release Button" located on the side of the instrument.
7. Firing the device can be achieved by pushing in on either the front trigger button or the rear trigger button. This action will trigger the spring mechanism, propelling the needle forward for a core biopsy sample.  
**NOTE:** The needle need not be removed from the biopsy instrument to retrieve the sample.
8. After removing the needle from the patient, reset the safety button and cock the actuator handle once to expose the biopsy notch located at the distal end of the needle. This action exposes the specimen for removal. Cocking the actuator handle a second time will activate the device for a second biopsy.  
**NOTE:** The disposable ACN™ needle set was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.

**REPROCESSING INSTRUCTIONS FOR THE PRO-MAG™ ULTRA 2.2 AUTOMATIC BIOPSY INSTRUMENT (Per ISO 17664-2004):**

<b>Warnings / Precautions</b>	<ul style="list-style-type: none"> <li>• Never place the device in bleach or bleach solution.</li> <li>• Lubrication on the sides of the carrier blocks will slow the action of the device</li> </ul>																
<b>Limitations</b>	After cleaning and sterilization, verify functionality prior to re-use.																
<b>INSTRUCTIONS</b>																	
<b>Point of Use</b>	This product is provided non-sterile and must be cleaned and sterilized before using.																
<b>Containment and Transportation</b>	No particular requirements																
<b>Preparation for cleaning</b>	No particular requirements. Ultrasonic treatment of the device has not been tested.																
<b>Cleaning: Automated</b>	<p>The Pro-Mag™ Automatic Biopsy instrument has been validated for effective cleaning using an automated washer/disinfector cycle of a minimum of 30 minutes total time, including a pre-wash, enzyme wash, main wash, rinse and dry.</p> <p>Immerse instrument, with lid open, in the washer-disinfector and refer to the manufacturer's recommendations and instructions. Examine device for any residual soil. Thoroughly dry the device or it will not function properly.</p>																
<b>Cleaning: Manual</b>	<p>The Pro-Mag™ Automatic Biopsy Instrument should be cleaned in a mild detergent (e.g. ENZOL Enzymatic Detergent) and lukewarm water solution, per manufacturer's recommendations.</p> <p>Immerse instrument, with lid open, into solution and use a soft brush to clean any residue left from the procedure. After washing completely, rinse the device in hot water. Thoroughly dry the device or it will not function properly.</p>																
<b>Disinfection</b>	A high level of disinfection can be achieved with CIDEX or ProCide NS; follow the manufacturer's recommended procedure. Thoroughly dry the device or it will not function properly.																
<b>Packaging</b>	A standard, sterilization wrap may be used. In the US, an FDA approved surgical wrap must be used.																
<b>Sterilization (Temperatures are minimum required, times are minimum required)</b>	The instrument should first be cleaned, lubricated and packaged prior to sterilization.																
	<table border="1"> <thead> <tr> <th>Cycle:</th> <th>Gravity</th> <th>Gravity</th> <th>Pre-vac</th> </tr> </thead> <tbody> <tr> <td>Temp.:</td> <td>250° F (121°C)</td> <td>270° F (132°C)</td> <td>270° F (132°C)</td> </tr> <tr> <td>Time:</td> <td>40 - 50 min</td> <td>10 - 25 min</td> <td>4 min</td> </tr> <tr> <td>Drying:</td> <td colspan="3">20 min, or until visibly dry</td> </tr> </tbody> </table>	Cycle:	Gravity	Gravity	Pre-vac	Temp.:	250° F (121°C)	270° F (132°C)	270° F (132°C)	Time:	40 - 50 min	10 - 25 min	4 min	Drying:	20 min, or until visibly dry		
	Cycle:	Gravity	Gravity	Pre-vac													
	Temp.:	250° F (121°C)	270° F (132°C)	270° F (132°C)													
Time:	40 - 50 min	10 - 25 min	4 min														
Drying:	20 min, or until visibly dry																
<b>Maintenance, Inspection and Testing</b>	<p>Washing detergents may eventually remove the lubricant originally applied. Lubricate device with silicone free, steam sterilization compatible lubrication (e.g., Aesculap® Sterili® 1 Oilspray JG 600) prior to sterilization. Lubricate at the following points as displayed on Figure 2 :</p> <ul style="list-style-type: none"> <li>1.1 Trigger Button (Must move freely)</li> <li>1.2 Top End of Actuator Pin</li> <li>1.3 Safety Button (Must move freely; lube both sides)</li> <li>1.4 Lever Pivot</li> <li>1.5 Top Edge of Ramps on Latch</li> <li>1.6 Sides of Catch on Stylet Holder Block</li> <li>1.7 Cam Surfaces</li> <li>1.8 Link Pivots</li> </ul>																
	After cleaning and sterilization, verify functionality prior to re-use.																
<b>Storage</b>	No specific requirements																
<b>Additional Information</b>	None																

Note: The instructions provided above have been validated by the manufacturer of the medical device as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using

equipment, materials and personnel in the reprocessing facility achieves the desired result. This normally requires validation and routine monitoring of the process.

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