**ARGON MEDICAL DEVICES, INC.**  
 1445 Flat Creek Road Athens, Texas 75751 USA  
 Tel: +1 800-927-4669;  
 Tel: +1 903-675-9321  
 www.argonmedical.com


EC REP **EMERGO EUROPE**  
 Prinsessegracht 20  
 2514 AP The Hague  
 The Netherlands  
 +31 70 345 8570

 The symbols glossary is located electronically at [www.argonmedical.com/symbols](http://www.argonmedical.com/symbols)

RxOnly        

**ENGLISH**

**Homer BLN  
Breast Localization Needle**

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

**Indication for Use:**

The Homer Breast Localization Needle (BLN) is intended to be used to localize nonpalpable breast lesions.

**Description:**

The Homer Breast Localization Needle (BLN) consists of a flexible nitinol "J" wire that can be retracted and repositioned to achieve accurate placement. The nitinol wire is loaded in a sharp 20ga stainless steel needle. A two-part stabilizer on the wire helps to advance and retract the wire.

| Catalogue # | Description   | Qty. per Box |
|-------------|---------------|--------------|
| 231030G     | 20ga x 3cm    | 10           |
| 231050G     | 20ga x 5cm    | 10           |
| 231075G     | 20ga x 7.5cm  | 10           |
| 231100G     | 20ga x 10cm   | 10           |
| 231125G     | 20ga x 12.5cm | 10           |

**Warnings:**

- Do not use the Homer BLN with an electrocautery device.
- The Homer BLN is not designed for use as a retractor. Do not use excessive force on any component.
- Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.
- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.

**How Supplied**

The Homer BLN is sterilized by ethylene oxide gas. Content is sterile in unopened and undamaged package. Do not use device if package has been damaged or has been opened.

**Storage**

Store in a cool, dry place.

**Disposal**

After use the Homer BLN may be potential biohazard. Handle and dispose of in accordance with accepted medical practice and all applicable laws and regulations.

**Instruction for Use**

- Select the correct length Homer BLN according to lesion depth. With "J" curve in extended position, tighten stabilizer securely to the wire against the needle hub. (Figure 1).
- Pulling on the stabilizer, retract the "J" curve into the needle. (Figure 2).
- After prepping and draping breast, and numbing the skin at the point of needle entry, insert the needle into the breast to the desired depth (Figure 3). In especially thick skin, a tiny pinpoint scalpel nick might be helpful.
- By pushing on the stabilizer while holding the hub, advance wire into breast until stabilizer butts against the hub. This ensures formation of the "J" curve in the parenchyma. The "J" tip will form in the parenchyma in line with the extension of the stabilizer. The extension should be aligned with the arrow engraved in the needle hub. Check proper needle position with craniocaudal and lateral mammograms (Figure 4).
- If initial wire placement is not satisfactory, reposition needle by pulling on the stabilizer while holding the hub, and retracting the wire to initial configuration (Figure 2). Redirect to proper position. It may not be necessary to remove the entire needle to reposition. Advance wire into breast until stabilizer butts against the needle hub. Repeat procedure as necessary until post-localization mammograms indicate proper needle/wire position (Figure5).
- When proper needle position is achieved, patient is sent to surgery with both needle and wire in vivo.
- If dye is to be used to mark a lesion, the stabilizer may be removed and dye can be injected through any standard IV tubing by attaching it to the hub (over the exposed end of the sterile wire). The dye will pass through the needle with the wire inside and the curve end anchored within the breast parenchyma. Reattach the stabilizer and send the patient to surgery with two markers in place, the dye and "J" wire with the needle.

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