ARGON MEDICAL DEVICES, INC.
 1445 Flat Creek Road
 Athens, Texas 75751 USA
 Tel: +1 800-927-4669;
 Tel: +1 903-675-9321
 www.argonmedical.com

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands
 +31 70 345 8570

Argon Medical Devices UK Ltd
 Eastgate Business Centre
 Eastern Avenue
 Burton-on-Trent
 Staffordshire
 DE13 0AT

X9585303/0321A

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
German	Das Symbolglossar ist elektronisch verfügbar unter www.argonmedical.com/symbols
French	Le glossaire des symboles est disponible sous forme électronique à l'adresse www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en la dirección www.argonmedical.com/symbols
Swedish	Symbolordlistan finns elektroniskt på www.argonmedical.com/symbols
Portuguese (EU)	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Finnish	Le glossaire des symboles est disponible sous forme électronique à l'adresse www.argonmedical.com/symbols
Norwegian	Symbolforklaringen finnes elektronisk på www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

ENGLISH

Hawkins, Accura and D Wire BLN Breast Localization Needle

Intended Use/Purpose

The BLN (Breast Localization Needle) are intended to be used with imaging techniques for accurate placement. Craniocaudal and lateral mammograms.

Indication for Use

Hawkins BLN are intended to be used to localize non-palpable breast lesions.

Accura and D Wire BLN devices are intended for use as a needle localization of breast lesions.

Device Description

Hawkins BLN offers with or without echogenic tip, and a barb retention design to provide strong anchoring strength in the breast. The device offers centimeter markings that provides reference for depth placement.

Accura BLN offers with and without stiffening cannula. The wire is marked in two places: I indicates wire at distal tip of needle; II indicates barb deployment.

D Wire BLN features a welded hookwire that help increase anchoring strength. It also offers a centimeter markings that provide reference for depth placement.

Warnings

- Rapid advancement of this needle into dense breast tissue may result in needle bending at the eyelet; advance carefully if resistance is experienced.¹
- Do not cut wire. Cutting this strandwire will cause unraveling and possible loose strand in operative site.² Under no circumstances should a hookwire engaged in tissue be removed without surgical intervention.

¹ = HAWKINS I & HAWKINS II
² = HAWKINS II FLEXSTRAND™ & HAWKINS III FLEXSTRAND™

- Do not use the BLN with an electrocautery device.
- BLN is not designed for use as a retractor. Do not use excessive force on any component.
- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged and expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- This device is intended for use by trained, qualified physician.

Possible complications

- Pain
- Bleeding at insertion site
- Infection at insertion site
- Wire migration
- Wire breakage

Direction for Use

Hawkins I BLN (Repositionable)

- Clean and prep area in a sterile fashion.
- Locate the lesion using appropriate radiographic technique.
- Insert needle into the breast and direct needle towards the lesion. Utilize the reference depth markers on the needle for accurate depth placement. (depth reference markers are 1cm apart)
- Confirm needle placement with radiographic technique. Adjust needle placement accordingly to assure adequate position within the lesion.
- Pull back on disk handle assembly to expose wire side barb to lock the needle in place. (Advance disk handle assembly to disengage barb and reposition if necessary.)
- Unscrew white needle grip from lock-down disk and remove to expose the distal end of needle.
- Loosen black screw on the side of lock-down disc and slide disk down to the skin entry site. Tighten black screw to secure disk in place.
- Cover/secure proximal end of needle if necessary, before transferring patient to surgical area.

Hawkins II BLN (Repositionable)

- Clean and prep area in a sterile fashion.
- Locate the lesion using appropriate radiographic technique.
- Insert needle into the breast and direct needle towards the lesion. Utilize the reference depth markers on the needle for accurate depth placement. (depth reference markers are 1cm apart)
- Confirm needle placement with radiographic technique. Adjust needle placement accordingly to assure adequate position within the lesion.
- With the needle tip in the desired position temporary side barb deployment can be checked by holding the wire in place and advancing the needle to the distal marker on the wire. This will expose the side barb through the temporary deployment window. (If any adjustments need to be made, retract needle and side barb back into the wire and adjust needle position.)
- With desired needle placement, deploy the wire by retracting the needle to the first distal wire marker at end of the needle hub. (This informs the user that wire is at the distal tip of needle.)
- Retract the needle to the second wire marker. (This indicates side barb is deployed.)
- Retract needle until third wire marker is at distal hub of needle (This indicates distal barb is deployed.)
- Needle may be removed, and imaging may be taken to confirm placement of the barbs within the lesion.
- Apply the skin retention clip at the skin entry site to prevent forward migration of the wire during patient transport.
- Cover/secure proximal end of wire if necessary, before transferring patient to surgical area.

Hawkins III (Non-repositionable needle)

- Clean and prep area in a sterile fashion.
- Locate the lesion using appropriate radiographic technique.
- Insert needle into the breast and direct needle towards the lesion. Utilize the reference depth markers on the needle for accurate depth placement. (depth reference markers are 1cm apart)
- Confirm needle placement with radiographic technique. Adjust needle placement accordingly to assure adequate position within the lesion.
- With the needle tip in the desired position advance the wire through the needle until the first wire marker is at the distal end of the needle hub. (This informs the user that wire is at the distal tip of needle.) *Be sure to insert wire with barb aligned with arrow on needle hub to assure proper deployment of barb.
- Retract the needle until second wire marker is at distal hub of needle (This indicates distal barb deployment)
- Needle may be removed, and imaging may be taken to confirm placement of the barbs within the lesion.
- Apply the skin retention clip at the skin entry site to prevent forward migration of the wire during patient transport.
- Cover/secure proximal end of wire if necessary, before transferring patient to surgical area.

Accura & Accura II BLN (Non-repositionable)

- Clean and prep area in a sterile fashion.
- Locate the lesion using appropriate radiographic technique.

- Insert needle into the breast and direct needle towards the lesion. Utilize the reference depth markers on the needle for accurate depth placement. (depth reference markers are 1cm apart)
- Confirm needle placement with radiographic technique. Adjust needle placement accordingly to assure adequate position within the lesion.
- With the needle tip in the desired position advance the wire through the needle until the first wire marker is at the distal end of the needle hub. (This informs the user that wire is at the distal tip of needle.) *Be sure to insert wire with barb aligned with arrow on needle hub to assure proper deployment of barb.
- Retract the needle until second wire marker is at distal hub of needle (This indicates distal barb deployment)
- Needle may be removed, and imaging may be taken to confirm placement of the barbs within the lesion.
- Apply the skin retention clip at the skin entry site to prevent forward migration of the wire during patient transport.
- Cover/secure proximal end of wire if necessary, before transferring patient to surgical area.

D Wire (Non-repositionable)

- Clean and prep area in a sterile fashion.
- Locate the lesion using appropriate radiographic technique.
- Insert needle into the breast and direct needle towards the lesion. Utilize the reference depth markers on the needle for accurate depth placement. (depth reference markers are 1cm apart)
- Confirm needle placement with radiographic technique. Adjust needle placement accordingly to assure adequate position within the lesion.
- With the needle tip in the desired position advance the wire through the needle until the first wire marker is at the distal end of the needle hub. (This informs the user that wire is at the distal tip of needle.) *Be sure to insert wire with barb aligned with arrow on needle hub to assure proper deployment of barb.
- Retract the needle until second wire marker is at distal hub of needle (This indicates distal barb deployment)
- Needle may be removed, and imaging may be taken to confirm placement of the barbs within the lesion.
- Apply the skin retention clip at the skin entry site to prevent forward migration of the wire during patient transport.
- Cover/secure proximal end of wire if necessary, before transferring patient to surgical area.

Disposal

After use the device may be potential biohazard. Handle and dispose of in accordance with accepted medical practice and all applicable laws and regulations.

Storage: Store at controlled room temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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