ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669
www.argonmedical.com


EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

Argon Medical Devices UK Ltd

Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT


X9585308/0321A


- 3**
EN: Medial Malleolus
BG: Среден малеол
CS: Malleolus medialis (vnitřní kotník)
DA: Malleolus medialis
DE: Malleolus medialis
EL: Έσω σφυρό
ES: Maléolo interno
ET: Mediaalne malleolus
FI: Sisäkehräs
FR: Malléole médiale
HR: Medijalni maleol
HU: Belső boka
IT: Malleolo Mediale
LT: Vidinė kulkšnis
LV: Mediālā potīte
NL: Mediale malleolus
NO: Medial malleolus
PL: Kość przysródkowa
PT: Maléolo Medial
RO: Maleolă medială
RU: Медиальная лодыжка
SK: Vnútorý členok (malleolus medialis)
SL: Medialni malleolus
SV: Mediala Malleolen
TR: Medial Malleol


<https://www.argonmedical.com/resources/product-information>

Pediatric Bone Marrow Access Needles

- 1**
EN: Epiphysal Plate
BG: Епифизна плоча
CS: Epifyzární ploténka
DA: Cartilago epiphysialis
DE: Epiphysenfuge
EL: Επιφυσαϊκή πλάκα
ES: Placa epifisaria
ET: Kasvuplaat
FI: Kasvulevy
FR: Plaque épiphysaire
HR: Epifizna ploča
HU: Növekedési porc
IT: Cartilagine Di Coniugazione
LT: Augimo plokštelė
LV: Kaulu epifizes plātne
NL: Epifysairschijf
NO: Epifysealplate
PL: Chrząstka nasadowa
PT: Cartilagem Epifisária
RO: Placă epifizară
RU: Эпифизарная пластинка
SK: Epifyzová platnička
SL: Epifizna plošča
SV: Cartilago Epiphysialis
TR: Epifiz Plağı

- 2**
EN: Tibial Tuberosity
BG: Издатина на тибиялната кост
CS: Tuberositas tibiae (drsnatina na přední proximální ploše tibie)
DA: Tuberositas tibiae
DE: Tuberositas tibiae
EL: Κνημιαίο κύρτωμα
ES: Tuberosidad tibial
ET: Sääreluu kõber
FI: Sääriluun kyhmy
FR: Tubérosité tibiale
HR: Tuberoza tibije
HU: Sípcontní tuberozítás (csontdudor)
IT: Tuberosità Tibiale
LT: Blauzdikaulio gumburai
LV: Tibiālā tuberozītāte
NL: Tuberositas tibiae
NO: Tibial tuberøsitet
PL: Guzowatość kości piszczelowej
PT: Tuberosidade Tibial
RO: Tuberozitate tibială
RU: Бугристость большеберцовой кости
SK: Drsnatina tibie (tuberositas tibiae)
SL: Tuberoznost golenice
SV: Tuberositas Tibiae
TR: Tibial Tüberozite

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Latvian	Simboliju žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Lithuanian	Simbolų glosarijs ir atrodams elektroniską veidą vietnė www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Russian	Электронный глоссарий символов расположен по адресу www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

Intended Use/Purpose: The Pediatric Bone Marrow Access Needle is ideal for intraosseous therapy. The needle can be used to access medullary cavities when initiating resuscitative infusion.

Device Description: The Pediatric Bone Marrow Access Needle offers 18ga regular wall (thicker wall) to provide extra support and each stylet has a trocar point for easy insertion. The luer lock connection is provided for secure syringe attachment.

Indications for Use: For obtaining access into medullary cavities for the purposes of initiating resuscitative infusion or for aspirating marrow in pediatric patients.

Contraindications: For use only for biopsies of bone/bone marrow as determined by a licensed physician. These needles should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of bone biopsy. Physician judgment is required when considering biopsy on patients with bleeding disorder or receiving anti-coagulant medications. This device should not be used in patients who have a compromised extremity, skin infection or comminuted fracture at insertion site. This device also should not be used in those patients who have osteogenesis imperfecta or osteopetrosis.

Possible Complications:

- Bruising
- Prolonged bleeding
- Soft tissue Damage
- Bone Fractures or injury to growth plates
- Infection

Warning

- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expir date has been exceeded.
- Do not continue to use if any of the component are damaged during procedure.
- Do not bend the needle excessively, it may result damage/breakage to the device.
- This device was designed, tested and manufactured for single patient use only.
- When used repeatedly in the same patient, the device should be inspected for damage or wear after each core sample is taken. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness and/or injury. DO NOT reuse, reprocess or re-sterilize this device.
- Reconditioning, refurbishing, repair, modification, of the device to enable further use is expressly prohibited.

Precaution

- This device should be manipulated only while under fluoroscopic observation with radiographic equipment that provides high quality images.
- Access to the vertebral body via pedicle requires a minimum pedicle width of 5mm.
- Knowledge of the dimensions of the site of insertion as assessed by MRI, CT, or other imaging method.

Note: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

Recommended Insertion Sites:

The preferred site is in the tibial tuberosity of the proximal tibia. The actual site is located 2cm below (inferior) the tuberosity.
A secondary site is the medial malleolus of the distal tibia.

Suggested Instructions for Insertion (for Proximal Tibial Tuberosity):

1. Stabilize limb in a semi-externally rotated position.
2. Prep area with iodine and/or alcohol.
3. Hold needle in place at right angle to the plate of the bone. The needle should be pointing anteroposterior with the leg semi-externally rotated. Ensure the needle is angled away from the knee joint. Special care should be taken to avoid the epiphyseal plate.
4. The insertion is done in two movements:
 - a) The first includes penetration of the skin and subcutaneous tissues.
 - b) The second includes penetration of the periosteum and the bone itself.
5. Using a firm twisting motion, pass the needle through the skin and cortex of the bone until an obvious "pop" or give is felt. There is usually no doubt as to the exact sensation.
6. Withdraw the stylet. It is common for a small amount of marrow and/or blood to seep out of the top of the needle after withdrawing the stylet.
7. Proper placement is confirmed by loss of resistance.
8. The needle will feel firmly fixed by the bone.
9. Aspiration should occur without resistance. It is important to stress that no material should be injected or infused until a "blood-marrow" mixture has been easily aspirated.
10. If necessary to further confirm location, inject 5cc of sterile saline and withdraw blood samples. Proper position can be assured by the free flow of blood during aspiration. Aspirate no more than 1cc of marrow.
11. Stabilize needle in a standard fashion at skin level and initiate infusion.

In the Event of Obstruction or Failure:

- a. Reassess insertion site.
- b. Flush cannula, re-insert stylet and repeat above procedure.
- c. If it is not possible to obtain aspirate or begin infusion, consider an alternate insertion site.

Needle Removal:

To remove the needle, apply gentle traction and rotate the needle cannula back and forth to disengage from bone cortex. Dress insertion site using accepted practices.

Special Considerations:

Most complications from intraosseous infusions occur after 24 hours of infusion. Since this procedure is intended as a "last-resort", it should only be used until conventional venous infusion can be initiated. Its use for longer than 24 hours is not recommended.

Disposal: After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

Storage: Store at controlled room temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

References:

J. Friery, K. Weiner: "Start an IV in that Bone," *Emergency*, November (1987), 28-33.

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