en	Instruction For Use
bg	Инструкции за употреба
cs	Pokyny pro použití
da	Brugsanvisning
de	Gebrauchsanweisung
el	Οδηγίες χρήσης
es	Instrucciones de uso
et	Kasutusjuhised
fi	Käyttöohjeet
fr	Mode d'emploi/Notice d'utilisation
hu	Használati utasítás
it	Istruzioni per l'uso
lt	Naudojimo instrukcijos
lv	Lietošanas instrukcija
nl	Gebbruiksaanwijzing
no	Bruksanvisning
pl	Instrukcja obsługi
pt	Instruções de utilização
ro	Instrucțiuni de utilizare
sk	Návod na použitie
sv	Bruksanvisning
tr	Kullanım talimatları

**ARGON MEDICAL DEVICES, INC.**

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669
www.argonmedical.com

**EMERGO EUROPE**

Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

The symbols glossary is located electronically at www.argonmedical.com/symbols

IFU7144/1019C

en - English**Indication for use**

The product is for nephrostomies

Duration

Up to 12 weeks

Warnings

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Do not use if package appears open or if the expiry date has been exceeded.
- This product contains the phthalate DEHP. When using the device, it should be taken into consideration that possible health affects when used for long durations in pregnant women, nursing mothers and children cannot be entirely discounted.
- The product must only be used by qualified personnel who are familiar with the technique.

Preparation

- Loosen the straightener together with the guidewire.
- Activate the hydrophilic coating by flushing the dispenser with saline.
- Advance the straightener carefully over the curvature of the catheter, while straightening the curvature with your fingers.
- The straightener now keeps the curvature straight allowing testing of the guidewire through the catheter before use.

Procedure

- Perform the skin incision under local anaesthesia.
- Introduce the needle as far as renal pelvis using ultrasound guidance, CT or fluoroscopy.
- Verify correct position by aspiration or by injection of contrast medium and fluoroscopy.
- When correctly positioned introduce the guidewire through the needle.
- Withdraw the needle over the guidewire.
- Dilate the canal.
- Introduce the catheter (with the straightener) over the guidewire.
- Remove the straightener
- When the catheter is positioned in renal pelvis, remove the guidewire.
- Connect the catheter to a drainage bag using a suitable connecting tube.

It is recommended affixing the catheter in a straight line and that any curvature be applied to the connecting tube.

It is recommended that the physician follow the hospital's standard of care procedures for drainage catheters.

Removal of catheter

- Disconnect the drainage bag connector tube from catheter.
- Pull the catheter out gently. If access is to be maintained, a straight floppy tip guidewire passed through the catheter will facilitate removal while maintaining access.

Storage

Store in a dark, dry, cool place. Avoid extended exposure to light.

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