en	Instructions For Use	تعليمات الاستعمال
ar		
cs	Pokyny pro použití	
da	Brugsanvisning	
de	Gebrauchsanweisung	
el	Οδηγίες χρήσης	
es	Instrucciones de uso	
et	Kasutusjuhised	
fi	Käyttöohjeet	
fr	Mode d'emploi	
hu	Használati utasítás	
it	Istruzioni per l'uso	
lt	Naudojimo instrukcijos	
nl	Gebbruiksaanwijzing	
no	Bruksanvisning	
pl	Instrukcja obsługi	
pt-br	Instruções de uso	
pt-eu	Instruções de utilização	
ro	Instrucțiuni de utilizare	
sk	Návod na použitie	
sl	Navodila za uporabo	
sv	Bruksanvisning	
tr	Kullanım talimatları	
Ru	Инструкция по применению	

**en - English**

**Caution:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

**Device description**

The SKATER® Introducer System is intended for fine needle puncture and tissue dilation with subsequent positioning of a guidewire for catheter placement in non-vascular procedures. The SKATER® Introducer consists of a 4F inner dilator with stiffening cannula that accepts a 0.018" guidewire, and a 6F outer sheath that accepts both 0.018" and 0.035" or 0.038" guidewires.

**Indication for use**

SKATER® Introducer Systems with radiopaque marker facilitate introduction and placement of a guidewire.

**Contraindications**

Physician judgement should be used in determining for which patients percutaneous interventional procedures are contraindicated.

**Duration**

Transient, less than 60 minutes.

**Warnings**

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Do not use if package appear open or if the expiry date has been exceeded.
- The product must only be used by qualified personnel who are familiar with the technique.

**Preparation of guidewire (if applicable)**

- Loosen the straightener together with the guidewire.
- Flush the dispenser with saline.
- Form the tip by running the spring tip over the edge of the colored straightener.

**Procedure**


- Perform the skin incision under local anaesthesia.
- Perform a fine needle puncture using ultrasound or X-ray.
- When the needle is correctly positioned, remove the stylet and introduce the 0.018" guidewire.
- Remove the needle and dilate with the assembled coaxial introducer system.
- Loosen the inner dilator, checking that the tip of the outer sheath has been properly placed inside the cavity (note the radiopaque ring marker).
- Remove the stiffening cannula and inner dilator.
- Introduce the 0.038" guidewire (not included).
- Remove the outer sheath and withdraw the 0.018" guidewire.
- Insert and advance the required catheter (not included) over the guidewire. Dilate further if necessary.

**Storage**

Store in a dark, dry, cool place. Avoid extended exposure to light.

**ARGON MEDICAL DEVICES, INC.**

1445 Flat Creek Road  
Athens, Texas 75751 USA  
Tel: +1 (903) 675 9321  
Tel: +1 (800) 927 4669  
www.argonmedical.com

**EMERGO EUROPE**

Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
+31 70 345 8570

The symbols glossary is located electronically at [www.argonmedical.com/symbols](http://www.argonmedical.com/symbols)

**Disposal**

After use, this product may be a potential biohazard. Dispose in accordance with applicable laws and regulations.

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