**Intended Use/Purpose**

The SKATER™ Nephrostomy Catheter is used for fluid drainage procedures.

**Device Description**

SKATER™ Nephrostomy Catheter offer locking or non-locking pigtail, with metal or flexible stiffening cannula. The product features SLIP-COAT™ hydrophilic coating at the distal tip for easy insertion, large lumen and large oval drainage holes designed for maximal drainage capacity.

**Indication for use**

The product is for nephrostomies.

**Duration**

Up to 12 weeks

**Warnings**

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package appear open or if the expiry date has been exceeded.
- During insertion, avoid contact with bone, cartilage, and scar tissue which can damage the catheter tip.

**Caution (Locking Pigtail):** If the catheter is to be removed in another department, we recommend that these guidance notes accompany the patient's case notes to ensure that the relevant personnel are aware of the presence of a locking catheter. It is also advisable to inform the patient.

**Precaution**

- The product must only be used by licensed, knowledgeable physician and qualified personnel who are familiar with the technique and who can ensure that patients will be monitored and appropriately treated for expected complications.
- Ensure the drain is secured and the system is intact to prevent dislodgement. Secure with a catheter fixation device, suture, or tape.
- Assess drain insertion site for signs of leakage, redness, or oozing. These signs may indicate an infection or irritation of the surrounding skin.
- Monitor changes in character or volume of fluid or bleeding.
- It is recommended affixing the catheter in a straight line and that any curvature be applied to the connecting tube.
- The drain also provides a pathway for bacteria to get into the wound. The risk of infection can be due to ascending bacterial invasion, foreign body reaction, decreased local tissue resistance, bacterial hiding places, poor placement due to kinked drain, and poor postoperative management.
- It is recommended that the physician follow the hospital's standard of care procedures for drainage catheters.
- As with any draining procedure, there is risk of infection, pain, discomfort, inadvertent tissue damage, and fluid loss. When used as intended by a physician, the benefits of using the SKATER™ Nephrostomy Catheter outweigh the risks associated with the use of the device.

**Preparation**

- Advance the straightener carefully over the curvature of the catheter, while straightening the curvature with your fingers.
- **Locking Pigtail** - Gently pull the thread to avoid unintentional looping of the thread.
- Introduce the chosen stiffener (if plastic: activate the hydrophilic coating with saline) and fasten the stiffener onto the catheter (Luer Lock).
- Remove the straightener
- Activate the coating by saline.

**Procedure**

- Select and prepare the drainage site using standard technique.
- Insert an appropriate guidewire using standard technique.
- Dilate the tract, if necessary.
- Advance the catheter over the guidewire in fluoroscopy.
- When the catheter tip is inside renal pelvis, loosen the stiffener.
- Hold the stiffener tight and advance the catheter over the stiffener until the complete length of the pigtail is in the cavity.
- **Non-Locking Pigtail**
  1. Remove the stiffener.
  2. Verify correct position of the catheter using fluoroscopy.
  3. Remove the guidewire.

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Български	Речникът на символите е наличен в електронен вид на адрес <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Česky	Slovníček symbolů v elektronické formě je umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Dansk	Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Deutsch	Das Symbolglossar ist elektronisch verfügbar unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Ελληνικά	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Español	El glosario de símbolos puede consultarse en línea en la dirección <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Eesti	Sümbolite sõnastik on elektrooniliselt saadaval aadressil <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Suomi	Symbolisanasto on sähköisessä muodossa osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .
Français	Le glossaire des symboles est disponible sous forme électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Magyar	A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Italiano	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Lietuvių k.	Simbolių žodynas elektronine forma pateiktas <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Latviešu val.	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Nederlands	Uitleg over de symbolen is elektronisch beschikbaar op <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Norsk	Symbolforklaringen befinner seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Polski	Wersja elektroniczna słownika symboli znajduje się na stronie <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Português	O glossário de símbolos está localizado eletronicamente em <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Română	Glosarul simbolurilor este stocat electronic la <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Svenska	Symbolordlistan finns elektroniskt på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Türkçe	Semboller sözlüğü elektronik ortamda <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> adresinde yer almaktadır

4. Connect the catheter to a drainage bag using a suitable connecting tube.
- **Locking Pigtail**
  1. Remove the stiffener and the guidewire.
  2. Gently pull the thread to ensure correct position and to complete curvature of the pigtail.
  3. Fix the thread tightly by winding it around the slot in the hub and finish by pressing the clip onto the slot and twist off the handle from the clip.
  4. Connect the catheter to a drainage bag using a suitable connecting tube.
  5. Verify correct position of the catheter using fluoroscopy.

#### Removal of catheter

- Disconnect the drainage bag connector tube from catheter.
- **Locking Pigtail** - Remove the clip and unwind the suture. Check that both threads are loose and cut one thread in order to loosen the pigtail.
- Pull the catheter out gently. If access is to be maintained, a straight floppy tip guidewire passed through the catheter will facilitate removal while maintaining access.

#### Disposal

After use, handle and dispose in accordance with hospitals policies and procedures concerning biohazard materials and waste.

#### Storage

Store in controlled room temperature.

**NOTE:** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at [quality.regulatory@argonmedical.com](mailto:quality.regulatory@argonmedical.com) as well as to the competent health authority where the user/patient resides.

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