## ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road  
Athens, Texas 75751 USA  
Tel: +1 (903) 675 9321  
Tel: +1 (800) 927 4669  
www.argonmedical.com

## Argon Medical Devices UK Ltd

Eastgate Business Centre  
Eastern Avenue  
Burton-on-Trent  
Staffordshire  
DE13 0AT

## EMERGO EUROPE

Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
+31 70 345 8570

35-2000-10C-0221/0221A

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Bulgarian	Речникът на символите е наличен в електронен вид на адрес <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Czech	Slovník symbolů v elektronické formě je umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Danish	Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
German	Das Glossar der Symbole kann unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Spanish	El glosario de símbolos puede consultarse en línea en la dirección <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Finnish	Elektroninen symbolisanasto on osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .
French	Le glossaire des symboles est disponible au format électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Lithuanian	Simbolių žodyno elektroninė forma pateikta <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Latvian	Elektroniskā simbolu vārdnīca ir pieejama <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Norwegian	Symbolforklaringen befinner seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Slovenian	Slovar simbolov se nahaja v elektronski obliki na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .
Turkish	Semboller sözlüğü elektronik ortamda <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> adresinde yer almaktadır

**Intended Use/Purpose**

Guidewires are intended for use in angiographic procedure to introduce and position catheters and interventional devices within the peripheral vasculature.

**Device Description**

Guidewires are designed to facilitate the passage of catheters and sheath introducers into the vasculature utilizing the Seldinger or Modified Seldinger techniques. The device is made with stainless steel and available with or without PTFE coating.

**Indications for Use**

Guidewires are indicated to facilitate the percutaneous placement of intravascular devices during diagnostic and intervention procedures.

**Contraindications for Use**

There are no known contraindications. Refer to the device label for any additional product-specific contraindication that may apply.

**Possible Complications/Adverse Events**

The following complications associated with Guidewire may include, but are not limited to:

- Air embolism
- Hematoma at the puncture site
- Infection
- Perforation of vessel walls
- Perforation of adjacent organs and/or tissues.

**Warnings**

- This device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and/or reuse.
- Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning.
- Guidewires are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- Withdrawal, pullback, or manipulation of the guidewire when resistance is met may cause guidewire damage, breakage, or embolization.
- Avoid withdrawing guidewire through metal needles, Guidewire may shear, or PTFE coating may scrape off against the needle bevel.
- Do not excessively torque the guidewire.
- Do not entrap or over-rotate the distal tip of the guidewire, it can lead to wire rupture.
- Do not excessively bend the guidewire, it may result in wire fracture.

**Precautions**

- Contents sterile unless package is damaged or opened.
- Inspect the package integrity before use.
- Guidewires are intended to be used by trained medical personnel in a clinical setting.
- The provided dispenser is best means of storage and handling of the guidewire, avoid coiling guidewires in less than 8" diameter during use.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- Use prior to the "Exp" date.
- If strong resistance is met during manipulation, discontinue the procedure and determine the cause for the resistance before proceeding.
- Avoid bending, kinking or modifying the shape of the wire.
- Guidewire placement should be monitored during placement and manipulation using fluoroscopic or suitable imaging method.
- Sufficient Guidewire length must remain exposed to maintain suitable grip proximally at all times.
- PTFE guidewires should not be sterilized by irradiation sterilization.
- Twisting or entanglement of the guidewires may occur when more than one wire is being used simultaneously. This can be prevented by carefully isolating and marking the proximal end of the wires.

**CAUTIONS for Movable Core Guidewires**

While the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

**Preparation**

1. Make sure the product is compatible with other equipment.  
**Note:** If unable to remove guide wire from dispenser, separate straightener from hoop dispenser.
2. Flush the guidewire with sterile heparinized normal saline or a similar isotonic solution by connecting a syringe to the luer hub of the guidewire dispenser.
3. Inject saline with syringe in a vertical position.

#### Direction for Use

4. Insert the flexible end of the guidewire into the needle. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device. If "J" wire is used, prepare "J" tip by sliding J-Straightener over tip to straighten. Insert J-tip inside J-Straightener into the needle hub gently.
5. Advance the guidewire through the needle.
6. Position the guidewire.
7. Hold the guidewire in place and withdraw the needle.
8. Advance the intravascular device or catheter over the guidewire.

#### Movable Guidewires

To lengthen the flexible segment of a movable-core guidewire, gently withdraw the core. Do not advance movable-core within vasculature.

#### Disposal

After use, handle and dispose in accordance with end-user facility's policies and procedures concerning biohazardous materials and waste.

#### Storage

Store at controlled room temperature.

**NOTE:** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at [quality.regulatory@argonmedical.com](mailto:quality.regulatory@argonmedical.com) as well as to the competent health authority where the user/patient resides.

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