
Endomyocardial Biopsy (EMB) Forceps - JAWZ™

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669
www.argonmedical.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Portuguese	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Russian	Электронный глоссарий символов расположен по адресу www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolförklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Chinese	符號詞彙表以電子方式在以下網站提供： www.argonmedical.com/symbols
Arabic	www.argonmedical.com/symbols يوفر مسرد الرموز إلكترونياً على

Made in USA
65-2000-01/ Rev. 0420U

EN - ENGLISH

Intended Use/Purpose

Endomyocardial Biopsy (EMB) Forceps are manually operated medical devices used to collect endomyocardial tissue samples from the inner wall of the heart for histopathological examination.

Device Description

The Endomyocardial Biopsy Forceps are radiopaque, handle-held and manually operated devices. It is actuated using a three-ring plastic handle. The spring in the handle maintains the jaws in a normally closed position. The thumb ring, at the proximal end, rotates to accommodate any thumb position. It has a flexible metal coil or plastic tube with a pair of scissors-like cups with cutting edges around their rims at the distal end which are operated by a handle attached at proximal end.

Indication for Use

The 50cm and 70cm Disposable Radiopaque Endomyocardial Biopsy Forceps are designed for Right Ventricular Biopsies using the Jugular Approach, and for Right or Left Ventricular Pediatric Biopsies using the Femoral Approach.

The 105cm Disposable Radiopaque Endomyocardial Biopsy Forceps are designed for Right or Left Ventricular Biopsies using the Femoral Approach.

Contra-indications for Use

There are no medical conditions (contraindications) associated with the Endomyocardial Biopsy Forceps. Assessment of patient overall physical state by the physician is required prior to use of the device.

Side Effect/Risks:

Endomyocardial biopsy (EMB) complication may include access site hematoma, transient right bundle branch block, transient arrhythmias, tricuspid regurgitation, and occult pulmonary embolism. When used as intended by a physician, the benefits of using the Endomyocardial Biopsy Forceps outweigh the risks associated with the use of the device.

Warnings

- Inspect the packaging and product for damage prior to use.
- Do not use if the package is open or damaged and if the expiry date has exceeded.
- The Biopsy Forceps are designed and intended for single, one patient use. DO NOT REUSE. Structural or functional integrity of the device may be compromised through cleaning, re-sterilization, and subsequent reuse.

Precautions

- This device is intended for use by trained medical personnel in interventional cardiology or heart transplantation.
- The Heart should be routinely monitored by EKG during the procedure.
- All mechanical forces are manually applied to the handle interface of the EMB forceps.
- The EMB forceps are subjected to forward axial motion via handle to activate the opening of the Jawz assembly, reverse motion via handle to activate the closing of the jaw, and shaft and distal tip rotation to traverse device into ventricle(s). Distal curve angulation is applied by the user to assist in the passage of the device when traversing the anastomotic juncture and passage across the tricuspid valve to the ventricle apex.
- The duration of the procedure is typically less than 1 hour.

SUGGESTED DIRECTIONS FOR USE

The Biopsy Forceps should be advanced through a Sheath Introducer or Guiding Catheter and into the right or left ventricle using fluoroscopic or echocardiographic guidance,

NOTE: Confirm the position of the Long Sheath Introducer or Guiding Catheter in the ventricle prior to and after each specimen collection.

- The jaws of the Biopsy Forceps should be opened as the forceps emerge from the distal end of the Long Sheath Introducer or Guiding Catheter.
- The opened jaws of the Biopsy Forceps should be positioned at the heart wall, closed firmly, and sufficient pressure maintained on the handle to assure retention of the specimen during withdrawal through the Sheath Introducer or Guiding Catheter.
- Continuous flushing of the Introducer Sheath or Guiding Catheter with Heparinized Saline should be conducted while slowly withdrawing the Biopsy Forceps while maintaining pressure to retain the specimen.
- Rinse and/or wipe the inside of the jaws of the biopsy forceps between specimen acquisition to remove biological debris. Ensure the forceps open and close prior to each insertion of the device into the introducer sheath or guide catheter.

Caution: The Forceps should be thoroughly rinsed with Heparinized Saline prior to each specimen collection during the procedure.

NOTE: Biopsy Forceps that are labeled as formable can be shaped at the distal tip. To form or modify the angle of the curve, grasp the distal tip between the thumb and index finger. Slide the fingers along tip body while bending to the desired angle. Do not kink or form the Biopsy Forceps to more than 90 degrees.

Disposal

After use, handle and dispose in accordance with hospital policies and procedures concerning biohazard materials and waste.

Storage

Store at controlled room temperature.

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