Pro-Mag™ Ultra Automatic Biopsy Instrument

Figure 1



Figure 2

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road,
Athens, Texas 75751 USA
Tel: +1 (903) 675-9321;
Tel: +1 (800) 927-4669
www.argonmedical.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

Argon Medical Devices UK Ltd
Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT

X9585760/1020A

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Czech	Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolu glosarijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols

Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar s simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolförklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Arabic	www.argonmedical.com/symbols يتوفر مسرد الرموز الإلكتروني على

EN - English

Pro-Mag™ Ultra Automatic Biopsy Instrument

Intended Use/Purpose

Pro-Mag™ Ultra Automatic Biopsy Instrument is a reusable device that utilizes disposable needles to harvest high quality histological specimens.

Device Description

Pro-Mag™ Ultra Automatic Biopsy Instrument features a built-in safety button to prevent accidental firing, easy to trigger from the front or rear of the instrument and easy specimen retrieval without removing the needle from the instrument. This device also provides automatic sequential firing of a stylet and subsequent firing of a cannula after cocking the handle twice and releasing the safety mechanism.

The Pro-Mag™ Ultra Automatic Biopsy Instrument are available in three (3) models:

- **Pro-Mag™ Ultra Automatic Biopsy Instrument** has a throw (advancement) of 25mm and used in conjunction with Pro-Mag™ Biopsy Needle Sets.
- **Pro-Mag™ Ultra ST Automatic Biopsy Instrument** has a throw (advancement) of 14mm and used in conjunction with Pro-Mag™ "Short Throw" Biopsy Needle.
- **Pro-Mag™ Ultra 2.2 Automatic Biopsy Instrument** has a throw (advancement) of 22mm and used in conjunction with ACN™ Biopsy Needle.

Note: The Pro-Mag Ultra Automatic Biopsy Instrument is designed to be compatible with listed Biopsy Needles, sold separately. Compatibility of the instrument with other manufacturer's needles has not been tested.

Indication for Use

The Pro-Mag™ Ultra Automatic Biopsy Instrument is indicated to obtain core samples from soft tissue such as the liver, kidney, prostate, breast and various soft tissue lesions for examinations and diagnosis.

Caution: Changes in the indication for use may result a new intended uses if new questions of safety and effectiveness are raised.

Contra-indications for Use

For use only for core biopsies of soft tissue as determined by a licensed physician. The instrument should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of core needle biopsy. Physician judgement is required when considering biopsy on patients with bleeding disorder or receiving anti-coagulant medications.

Warnings

- The instrument is reusable and is supplied non-sterile.
- The instrument is intended to be cleaned, disinfected and sterilized by the user as defined in the reprocessing instruction.
- Do not continue to use if any of the component are damaged during the procedure.
- Hard tissue biopsy is considered an off-label misuse of the instrument.

Precautions

- The product must only be used by licensed, knowledgeable physician and qualified personnel who are familiar with possible side effects, typical findings, limitations, indications and contraindications of this device and biopsy sampling.
- These instructions for the Pro-Mag™ Ultra Automatic Biopsy Instrument are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with the device.
- The device will typically be utilized in an outpatient clinical environment and in the following clinical areas:
 - > Urology – Transrectal or transperineal biopsy of the prostate.
 - > Radiology – Percutaneous kidney, liver and breast biopsy

- Before loading needles into the instruments, slide cannula over the stylet several times noting the proper free movement when necessary.
- DO NOT use if needles specimen notch does not pass into the outer cannula freely. (Note: A slight drag may be felt in the final 6mm of needle set closure.)
- Each device can be reused in the same procedure on the same patient for multiple biopsies.
- When used repeatedly in the same patient, the devices should be inspected for damage or wear after each sample is taken.
- Never test fire automatic reusable biopsy instrument while loaded with a needle set, damage to needle set may occur.
- Patient movement may affect the outcome of the procedure relating to quality of sample.

Side Effect

- It is not uncommon for patients to experience bleeding, infection at insertion site, fever, pain at the biopsy site, swelling at the biopsy site, adjacent tissue injury and hematoma after the biopsy. When used as intended by a physician, the benefits of using the Biopsy Instrument Needle outweigh the risks associated with the use of the device.

Reprocessing Instructions for Pro-Mag Ultra Automatic Biopsy Instrument (Per ISO 17664:2017/AAMI TIR 30:2011(R)2016)

Warnings / Precautions	<ul style="list-style-type: none"> • Never place the device in bleach or bleach solution. • Lubrication on the sides of the carrier blocks will slow the action of the device 			
Limitations	After cleaning and sterilization, verify functionality prior to re-use.			
INSTRUCTIONS				
Point of Use	This product is provided non-sterile and must be cleaned and sterilized before using.			
Containment and Transportation	No particular requirements			
Preparation for cleaning	No particular requirements. Ultrasonic treatment of the device has not been tested.			
Cleaning: Automated	<p>The Pro-Mag™ Automatic Biopsy instrument has been validated for effective cleaning using an automated washer/disinfector cycle of a minimum of 30 minutes total time, including a pre-wash, enzyme wash, main wash, rinse and dry.</p> <p>Immerse instrument, with lid open, in the washer-disinfector and refer to the manufacturer's recommendations and instructions. Examine device for any residual soil. Thoroughly dry the device or it will not function properly.</p>			
Cleaning: Manual	<p>The Pro-Mag™ Automatic Biopsy Instrument should be cleaned in a mild detergent (e.g. ENZOL Enzymatic Detergent) and lukewarm water solution, per manufacturer's recommendations.</p> <p>Immerse instrument, with lid open, into solution and use a soft brush to clean any residue left from the procedure. After washing completely, rinse the device in hot water. Thoroughly dry the device or it will not function properly.</p>			
Disinfection	A high level of disinfection can be achieved with CIDEX or ProCide NS; follow the manufacturer's recommended procedure. Thoroughly dry the device or it will not function properly.			
Packaging	A standard, sterilization wrap may be used. In the US, an FDA approved surgical wrap must be used.			
Sterilization (Temperatures are minimum required, times are minimum required)	The instrument should first be cleaned, lubricated and packaged prior to sterilization.			
	Cycle:	Gravity	Gravity	Pre-vac
	Temp.:	250° F (121°C)	270° F (132°C)	270° F (132°C)
	Time:	40 - 50 min	10 - 25 min	4 min
	Drying:	20 min, or until visibly dry		
Maintenance, Inspection and Testing	<p>Washing detergents may eventually remove the lubricant originally applied. Lubricate device with silicone free, steam sterilization compatible lubrication (e.g., Aesculap® Sterilit® 1 Oilspray JG 600) prior to sterilization. Lubricate at the following points as displayed on Figure 2:</p> <ol style="list-style-type: none"> 1.1 Trigger Button (Must move freely) 1.2 Top End of Actuator Pin 1.3 Safety Button (Must move freely; lube both sides) 1.4 Lever Pivot 1.5 Top Edge of Ramps on Latch 1.6 Sides of Catch on Stylet Holder Block 1.7 Cam Surfaces 1.8 Link Pivots <p>After cleaning and sterilization, verify functionality prior to re-use.</p>			
Storage	No specific requirements			
Additional Information	None			

Note: The instructions provided above have been validated by the manufacturer of the medical device as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieves the desired result. This normally requires validation and routine monitoring of the process.

Product Lifetime

- The life of the instrument is determined by the frequency of use by the user. The Pro-Mag™ Biopsy Instrument is qualified for 5,000 cocking and firing cycles.

Pre-Procedure Preparation

- Inspect the needle product package for damage.
- If undamaged, open the needle package using proper aseptic technique.
- Remove the needle from the pouch. Inspect the stylet for damage to the point, the cannula for damage to the cutting edge, or other imperfections that would prevent proper operation of the needle assembly.
- Open the lid of the biopsy instrument device, load the needle into the biopsy instrument device, close the lid and remove protective sheath on the needle.
- Safety button must be set in the position indicated by the arrow on the lid before cocking.
- Percutaneous introduction of the automatic cutting needle should be performed under ultrasound, computerized tomography (CT) or fluoroscopy in order to maximize needle placement.

Procedure

- Cock the actuator handle twice in order to activate the device for firing (this automatically sets the safety).
- Release the safety by pushing the "Safety Release Button" located on the side of the instrument device.
- Firing the device can be achieved by pushing in on either the front trigger button or the rear trigger button. This action will trigger the spring mechanism, propelling the needle forward for a core biopsy sample.
Note: The needle need not be removed from the biopsy device to retrieve the sample.
- After removing the needle from the patient, reset the safety button.
- Cock the actuator handle once to expose the biopsy notch located at the distal end of the needle. This action exposes the specimen for removal.
- Cocking the actuator handle a second time will activate the device for a second biopsy.

Disposal

After use, dispose the needle in accordance with hospitals policies and procedures concerning biohazard materials and waste. The Pro-Mag™ Ultra Automatic Biopsy Instrument is intended to be cleaned, disinfected and sterilized by the user as defined in the reprocessing instruction. However, if the device instrument reached to its end life cycle, then dispose in accordance with all federal, state and/or local regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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