Homer Breast Localization Needle (BLN)

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road Athens, Texas 75751 USA
Tel: +1 800-927-4669
Tel: +1 903-675-9321
www.argonmedical.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
+31 70 345 8570

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
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

©2013-2017 Argon Medical Devices, Inc. All rights reserved.
X9585295-0321
Date of Issue: 05/15/2023 Rev. B

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en la dirección www.argonmedical.com/symbols
French	Le glossaire des symboles est disponible sous forme électronique à l'adresse www.argonmedical.com/symbols
Portuguese (BR)	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Portuguese (EU)	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
German	Das Symbolglossar ist elektronisch verfügbar unter www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Swedish	Symbolordlistan finns elektroniskt på www.argonmedical.com/symbols
Greek	Το γλῶσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Finnish	Le glossaire des symboles est disponible sous forme électronique à l'adresse www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Polish	Wersja elektroniczna słownika symboli znajduje się na stronie www.argonmedical.com/symbols
Czech	Slovníček symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols

Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Lithuanian	Simbolių žodynas elektronine forma pateiktas www.argonmedical.com/symbols
Estonian	Sümbolite sõnastik on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Arabic	يُوفّر مصدر الرموز إلكترونيًا على www.argonmedical.com/symbols

ENGLISH

Homer BLN Breast Localization Needle

Intended Use/Purpose :

The Homer BLN contains a strong, flexible J-curve wire that may be retracted and repositioned as often as necessary to achieve accurate placement.

Indication for Use:

The Homer BLN is intended to be used to localize nonpalpable breast lesions.

Device Description:

The Homer BLN consists of a flexible nitinol "J" wire that can be retracted and repositioned to achieve accurate placement. The nitinol wire is loaded in a sharp 20ga stainless steel needle. A two-part stabilizer on the wire helps to advance and retract the wire.

Catalogue #	Description	Qty. per Box
231030G	20ga x 3cm	10
231050G	20ga x 5cm	10
231075G	20ga x 7.5cm	10
231100G	20ga x 10cm	10
231125G	20ga x 12.5cm	10

Contra-indications for Use

There are no contra-indications.

Side Effects

There are no known side effects.

Possible complications

- Pain
- Bleeding at insertion site
- Infection at insertion site
- Wire migration
- Wire breakage

Warnings:

1. Do not use the Homer BLN with an electrocautery device.

- The Homer BLN is not designed for use as a retractor. Do not use excessive force on any component.
- Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.
- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged and expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- This device is intended for use by trained, qualified physician.

Precautions:

- Rapid advancement of this needle into dense breast tissue may result in needle bending at the eyelet; advance carefully if resistance is experienced.¹
- Do not cut wire. Cutting this standwire will cause unraveling and possible loose strand in operative site.²
- Under no circumstances should a hookwire engaged in tissue be removed without surgical intervention.

Pre-Procedural Preparation

- Select the correct length Homer BLN according to lesion depth. With "J" curve in extended position, tighten stabilizer securely to the wire against the needle hub. (Figure 1).
- Pulling on the stabilizer, retract the "J" curve into the needle. (Figure 2).

Directions for Use

- After prepping and draping breast and numbing the skin at the point of needle entry, insert the needle into the breast to the desired depth (Figure 3). In especially thick skin, a tiny pinpoint scalpel nick might be helpful.
- By pushing on the stabilizer while holding the hub, advance wire into breast until stabilizer butts against the hub. This ensures formation of the "J" curve in the parenchyma. The "J" tip will form in the parenchyma in line with the extension of the stabilizer. The extension should be aligned with the arrow engraved in the needle hub. Check proper needle position with craniocaudal and lateral mammograms (Figure 4).
- If initial wire placement is not satisfactory, reposition needle by pulling on the stabilizer while holding the hub, and retracting the wire to initial configuration (Figure 2). Redirect to proper position. It may not be necessary to remove the entire needle to reposition. Advance wire into breast until stabilizer butts against the needle hub. Repeat procedure as necessary until post-localization mammograms indicate proper needle/wire position (Figure 5).
- When proper needle position is achieved, patient is sent to surgery with both needle and wire in vivo.
- If dye is to be used to mark a lesion, the stabilizer may be removed and dye can be injected through any standard IV tubing by attaching it to the hub (over the exposed end of the sterile wire). The dye will pass through the needle with the wire inside and the curve end anchored within the breast parenchyma. Reattach the stabilizer and send the patient to surgery with two markers in place, the dye and "J" wire with the needle.

Disposal

After use the Homer BLN may be potential biohazard. Handle and dispose of in accordance with accepted medical practice and all applicable laws and regulations.

Storage: Store at controlled room temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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