ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 800-927-4669;
Tel: +1 903-675-9321
www.argonmedical.com


EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
+31 70 345 8570

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Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλῶσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Latvian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Lithuanian	Simbolų glosarijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen finnes seg elektronisk på www.argonmedical.com/symbols
Polish	Głosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Russian	Электронный глоссарий символов расположен по адресу www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

X9585309-0321

Date of Issue: 07/31/2023 Rev. B

ENGLISH

"I" Type Bone Marrow Aspiration Needle

Intended Use/Purpose: The "I" Type Bone Marrow Aspiration Needle is used to control depth penetration for aspiration at the sternal or iliac crest.

Devices Description: The device features a luer lock connector on the handle for secure syringe attachment. It has an adjustable, removable needle stop that safely controls depth of penetration.

Indications for Use: For Sternal or Iliac Crest Aspiration.

Contraindications: For use only for biopsies of bone/bone marrow as determined by a licensed physician. These needles should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of bone biopsy. Physician judgment is required when considering biopsy on patients with bleeding disorder or receiving anti-coagulant medications.

Possible complications:

Potential Complications/Risks: It is not uncommon for patients to experience the following :

- Nerve injury including puncture of the spinal cord or nerve roots potentially resulting in radiculopathy, paresis or paralysis.
- Bleeding or hemorrhage
- Hematoma
- Pain
- Allergic Reaction
- Discomfort

Warning

- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expir date has been exceeded.
- Do not continue to use if any of the component are damaged during procedure.
- Do not bend the needle excessively, it may result damage/breakage to the device.
- This device was designed, tested and manufactured for single patient use only.
- When used repeatedly in the same patient, the device should be inspected for damage or wear after each core sample is taken. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness and/or injury. DO NOT reuse, reprocess or re-sterilize this device.
- Reconditioning, refurbishing, repair, modification, of the device to enable further use is expressly prohibited.

Precaution

- This device should be manipulated only while under fluoroscopic observation with radiographic equipment that provides high quality images.
- Access to the vertebral body via pedicle requires a minimum pedicle width of 5mm.
- Knowledge of the dimensions of the site of insertion as assessed by MRI, CT, or other imaging method.

Note: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

Instructions for Use:

1. Using sterile technique, prepare the skin with antiseptic and drape.
2. Infiltrate the selected area with local anesthetic.

For Sternal Aspiration (Steps 3-4):

3. Adjust the needle stop for correct needle length. The adjustable screw-stop allows for a minimum needle length of 3/8"; when the screw-stop is removed, the needle length is 1-7/8".

Note: Do not use for sternal puncture with stop assembly removed. (See Iliac Crest instructions).

4. Insert needle into the sternum in usual fashion. Proceed with Step 7.

For Iliac Crest Aspiration (Steps 5-6):

5. The screw-stop can be removed from the screw thread for a needle length of 1-7/8". If more length is desired, screw stop all the way up (with needle guard in place), grab stop assembly and pull briskly straight along needle. The entire needle stop assembly will now be removed and the needle length will be 2-5/8".

Note: Do not use if needle is bent.

6. Insert needle into Iliac crest in usual fashion.
7. Once needle tip is positioned in marrow cavity, remove stylet by rotating upper section of handle 90°, and pulling straight out.
8. A luer-lock connection is provided for syringe connection.
Note: Either luer-lock or luer-slip syringe may be used.
9. Aspirate marrow and process samples in accordance with hospital procedure.
10. Dispose of all used needles in a correct and safe manner.

Disposal: After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

Storage: Store at controlled room temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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