COVIDIEN

Guidewires

Multi-language Instructions For Use
(IFU)

TABLE OF CONTENTS


Instructions for use:

English..... 1

French..... 3

Portuguese..... 5

Spanish..... 7

COVIDIEN

ENGLISH

GUIDEWIRES INSTRUCTIONS FOR USE

DESCRIPTION

Stainless steel Guidewire with or without PTFE coating and with or without heparin bonding.

INDICATIONS FOR USE

To facilitate the percutaneous placement of intravascular devices during diagnostic and intervention procedures. Refer to the device label for any additional product-specific indications that may apply. Prior to use please see the complete 'Instructions For Use' for more information on Indications, Contraindications, Warnings, Cautions, Adverse Events and Guidelines for Use.

CONTRAINDICATIONS

There are no known Contraindications. Refer to the device label for any additional product-specific contraindications that may apply.

WARNINGS

1. This Device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and/or reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Guidewires are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
2. Withdrawal, pullback, or manipulation of the guidewire when resistance is met may cause guidewire damage, breakage, or embolization.
3. Avoid withdrawing guidewire through metal needles; guidewires may shear or PTFE coating may scrape off against the needle bevel.

CAUTIONS

1. Store in a cool, dark, dry place.
2. The provided dispenser is the best means of storage and handling of the guidewire, avoid coiling guidewires in less than 8" diameter during use. Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire. Use prior to the "Exp" date.
3. If strong resistance is met during manipulation, discontinue the procedure and determine the cause for the resistance before proceeding. Avoid bending, kinking or modifying the shape of the wire.

4. Guidewire placement should be monitored during placement and manipulation using fluoroscopic, ultrasound, or suitable imaging method.
5. Sufficient Guidewire length must remain exposed to maintain suitable grip proximally at all times
6. PTFE guidewires should not be sterilized by irradiation sterilization.
7. Contents sterile unless package is damaged or opened.
8. Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

ADVERSE EVENTS

The following complications may occur: air embolism, hematoma at the puncture site, infection, or perforation of the vessel wall.

GUIDE LINES FOR USE

1. Flush the guidewire with sterile heparinized normal saline or a similar isotonic solution by connecting a syringe to the luer hub of the guidewire dispenser.
2. Insert the flexible end of the guidewire into the needle. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device. If “J” wire is used, prepare “J” tip by sliding J-Straightener over tip to straighten. Insert J-tip inside J-Straightener into needle hub gently.
3. Advance the guidewire through the needle.
4. Position the guidewire.
5. Hold the guidewire in place and withdraw the needle.
6. Advance the intravascular device or catheter over the guidewire.

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Manufactured for Covidien, LLC.