Traveler™21 Portal Vein Access Set

ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669

EMERGO EUROPE

Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Emergo Consulting (UK) Limited

c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

MedEnvoy Switzerland

Golthardstrasse 28
6302 Zug
Switzerland

<https://www.argonmedical.com/resources/product-information>

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EN - ENGLISH

Device Description


The device is for single-use, supplied sterile and sterilized by Ethylene Oxide.

The Portal Vein Access Set contains a 5F MPA catheter, a 14ga Stiffening Cannula with a 7F Cannula Sheath, a 21ga Needle with a 5F Needle Catheter (separated with a removable spacer clip), and a bending tool. The 14ga stiffening cannula with cannula sheath has a directional handle that indicates the direction of the curve.

These components are used through any compatible introducer sheath to create a pathway through the liver parenchyma through which an endoprosthesis can be delivered. The Portal Vein Access Set is used to gain access to the hepatic vein and guide a sharp puncture tool toward the parenchyma. The puncture tool is used to make a pathway from the hepatic vein to the portal vein, and then the pathway is dilated to provide access for a larger sheath. The shunt is inserted through the sheath and deployed through the pathway. Then, all of the Portal Vein Access Set components are removed.

The Traveler™21 Portal Vein Access Set is compatible with a 10 F introducer sheath no longer than 45.5cm in overall length.

Set Components

- A. 5F x 80cm MPA Catheter
- B. 14ga Stiffening Cannula / 7F x 52cm Cannula Sheath
- C. 21ga x 72cm Coaxial Needle
- D. 5F x 67cm Needle Catheter
- E. Bending Tool

Intended Purpose/ Intended Use

The Traveler™21 Portal Vein Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

Intended Patient Population

Patients may be adults of any gender, any race or ethnicity, undergoing transjugular intrahepatic portosystemic shunt (TIPS) as treatment for the following:

- Variceal bleeding, secondary prevention, and acute bleeding refractory to medical and endoscopic treatments
- Refractory ascites
- Hepatorenal syndrome (types 1 and 2)
- Budd–Chiari syndrome
- Hepatic veno-occlusive disease
- Hepatic hydrothorax
- Portal hypertensive gastropathy
- Hepatopulmonary syndrome
- Portal vein thrombosis

Indication for Use

The Traveler™21 Portal Vein Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

Intended Users

The product is intended for use by medical professionals and must only be used by qualified personnel who are familiar with the technique.

Duration/Lifetime

In use in procedures up to 4 hours

20-2007-45

Date of Issue: 2025/04/23 Rev F

Clinical Benefits

The Portal Vein Access sets provide an indirect clinical benefit by facilitating creation of the parenchymal tract between the hepatic vein and the portal vein for delivery of a transjugular intrahepatic portosystemic shunt used to treat the complications of portal hypertension.

Risk and Side-Effects

- Intraoperative hemorrhage
- Puncture site hematoma
- Gall bladder puncture
- Cardiac arrhythmia
- Arteriovenous fistula
- Arterio-biliary fistula
- Death

Contraindications

None known

Warnings

- Contents are supplied sterile and are intended for single use only. Do not re-sterilize.
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection, or other injury.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure
- Do not use excessive force when advancing, retracting, or manipulating kit components and accessories.
- Use of power injection of contrast media to verify access to the portal vein may cause a loss of access.
- A patient may experience a negative reaction to angiography if the patient is allergic to contrast media or has compromised renal function.
- Risks and benefits should be carefully assessed in pregnant patients as radiation from fluoroscopic imaging may endanger the fetus.
- Do not bend the needle/catheter assembly more than 40° at any single point, as this might cause kinking of the needle/catheter.

Precautions

- Inspect package integrity before use.
- Inspect all set components prior to use
- For use by physicians trained in vascular diagnostic and interventional techniques only
- Do not overtighten hub connections.
- Use fluoroscopic guidance when advancing, retracting, or manipulating kit components and accessories in the vasculature.
- The patient should be monitored (EKG, oximetry, blood pressure, pulse) during the procedure.
- If the Portal Vein Access set is to be used for shunt delivery, the introducer sheath selected for portal vein access must also be compatible with the shunt, per the manufacturer recommendation.

Instructions for Use

1. Remove all components from packaging using sterile technique.
2. Flush components with saline (not supplied).
3. Assemble the 10F Introducer Sheath (not supplied) with the 10F Dilator (not supplied).
4. Assemble the 14ga Stiffening Cannula with the 7F Cannula Sheath
5. If necessary, apply the Bending Tool over the 14ga Stiffening Cannula and 7Fr Cannula Sheath and adjust the angle of the cannula/sheath assembly.
6. Obtain jugular access and advance a $\leq 0.038''$ guidewire through the jugular vein, into the inferior vena cava.
7. Advance a compatible 10F Introducer Sheath/Dilator assembly over the guidewire into the IVC.
8. Remove the 10F Dilator.
9. Insert the 5F MPA Catheter over the guidewire into the hepatic vein. Advance the guidewire into the hepatic vein.

10. Advance the 10F Introducer Sheath into the hepatic vein and maintain access without occluding the vessel.
11. Remove the 5F MPA Catheter, leaving the guidewire in place.
12. Introduce the 14ga Stiffening Cannula/7F Cannula Sheath assembly into the hepatic vein and remove the guidewire from the 14ga Stiffening Cannula. Be careful to maintain hepatic access with the 14ga Stiffening Cannula and Introducer Sheath while removing the guidewire.
13. Insert 21ga needle/5F Needle Catheter assembly into 14ga Stiffening Cannula until the black line on the Needle Catheter aligns with the end of the safety funnel on the 14ga Stiffening Cannula. The 21ga needle tip should not be advanced past the catheter or cannula tip during insertion.
14. Remove the spacer clip from the needle/catheter assembly and retract the 5F Needle Catheter until the hub is locked with the needle handle to expose the needle tip.
15. Expose the needle tip by connecting the hub of the 21ga needle hub to the 5F Needle Catheter hub. Ensure the 21ga needle tip does not advance past the distal tip of the 14ga Stiffening Cannula.
16. Orient the 14ga Stiffening Cannula toward the portal vein (directional handles on needle and stiffening cannula indicate the direction of the needle bevel and cannula curve).
17. Wedging the 14ga Stiffening Cannula against the vein wall, advance the 21ga needle forward through the hepatic parenchyma and toward the portal vein. If the 21ga Needle is advanced beyond the portal vein, retract it until the tip of the catheter is within the portal vein.
18. Remove the stylet from the needle and confirm access to the portal vein.
19. Advance the 5F Needle Catheter over the 21ga needle and into the portal vein.
20. Remove the 21ga needle from the 5F Needle Catheter.
21. Introduce guidewire (not supplied) through the 5F Needle Catheter and into the portal vein.
22. If necessary, dilate the parenchymal tract to accommodate the 10F Introducer Sheath.
23. Advance the 10F Introducer Sheath across the parenchymal tract until the tip of the Introducer Sheath is in the portal vein. Be careful to avoid occluding the portal vein with the Introducer Sheath.
24. Remove the 5F Needle Catheter, 14ga Stiffening Cannula, and 7F Cannula Sheath, leaving the Introducer Sheath in place.
25. Proceed with any indicated diagnostic or interventional procedures through the 10F sheath.

Storage

Store at controlled room temperature.

Disposal

After use, handle and dispose in accordance with hospital's policies and procedures concerning biohazard materials and waste.

Note: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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