
**First PICC™
First Midcath™
L-Cath™
PICC and Midline Catheter™**

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669
www.argonmedical.com

Brasil Detentor de Registro: CM HOSPITALAR S.A.
Endereço: Avenida Luiz Maggioni, 2727 – Distrito Empresarial - CEP: 14072-055 - Ribeirão Preto/SP
Registro nº: 80743990006, 80743990008, 80743990009, 80743990010, 80743990011, 80743990013

 <https://www.argonmedical.com/resources/product-information>

PN 20-2007-07

Date of Issue: 12/07/2023 Rev. H

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Português do Brasil	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Français	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
繁體中文	符號詞彙表以電子方式在以下網站提供： www.argonmedical.com/symbols

Catheter	Size (GA)	Catheter Length (CM)	Priming Volume (mL)	Water Flow Rate - Gravity (mL/hr)	Water Flow Rate - 8 psi. (mL/hr)
L-Cath	28Ga	25	0.06	7	30
L-Cath	26GA	30	0.1	50	210
L-Cath	20GA 2 lumen	60	22GA 0.16 23GA 0.16	15 7	80 40
L-Cath	20GA	20	0.13	400	>1000
L-Cath	20GA	60	0.27	260	>1000
L-Cath	18GA	20	0.13	600	>1000
L-Cath	18GA	60	0.28	400	>1000
L-Cath	18GA 2 lumen	60	20GA 0.31 21GA 0.26	80 40	465 175
L-Cath	18GA	60	0.72	>1000	>1000
L-Cath	16GA 2 lumen	60	18GA 0.60 19GA 0.48	520 250	>1000 >1000
Y-Site			0.33		
First PICC	26GA (1.9F)	50	0.13	20	90
First PICC	22GA PED (2.8F)	50	0.19	125	600
First Midcath	20GA (3F)	20	0.16	260	>1000
First PICC	20GA (3F)	65	0.26	110	530
First Midcath	18GA (4F)	20	0.19	525	>1000
First PICC	18GA (4F)	65	0.33	225	>1000
First PICC	16GA (5F)	65	0.41	400	>1000
First PICC	16GA (5F) 2 lumen	65	20GA 0.27 20GA 0.27	150	720
T-Port			0.32		

Flow rates and priming volumes indicated are for the untrimmed catheter without T-Port or Y-Site.

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ENG-English

INTENDED USE/PURPOSE

Peripherally inserted central catheters (PICC) are intended to provide short or long term venous access to the large central veins and are generally used to give medications or liquid nutrition.

DEVICE DESCRIPTION

A family of peripherally inserted central catheters made of polyurethane. Catheters are sold as standalone devices which include the following accessories: Trim tool, tape measure, and stylet. Catheters are also sold as a component in a basic kit (catheter and introducer packaged together) which includes the following: BD Introsyte Autoguard Introducer or Splittable Introducer, stylet, trim tool, tourniquet.

INDICATION FOR USE

The FIRST PICC™, First Midcath™, L-CATH Single, Dual Lumen and Midline Catheters are indicated for short or long term peripheral access to the central venous system for the administration of fluids, medications, and nutrients; the sampling of blood; and monitoring blood pressure and temperature intravenously.

CONTRAINDICATIONS FOR USE

- The devices are contraindicated whenever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway.
- The devices are contraindicated when the patient's body size is insufficient to accommodate the size of the device.
- The patient is known or is suggested to be allergic to materials contained in the device.
- Clinical assessment of the patient must be completed to ensure no contraindications exist.

SIDE EFFECTS / ADVERSE EVENTS

- Air Embolism
- Bleeding
- Cardiac Arrhythmia / Tamponade
- Catheter Embolism
- Catheter Occlusion
- Exit Site Infection / Necrosis
- Hematoma
- Intolerance Reaction to Implanted Device
- Laceration of Vessels of Viscus
- Thromboembolism

WARNINGS

- Contents are supplied sterile and are intended for single use only. Do not re-sterilize
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- Do not use as a right atrial or umbilical catheter. PICC catheters shall have the distal tip dwelling in the lower one third of the superior vena cava to the junction of the superior vena cava and the right atrium.
- Lower extremity placement shall have the distal tip dwelling in the thoracic inferior vena cava above the level of the diaphragm.
- Midline catheters shall have the distal tip residing in the basilic, cephalic, or brachial vein at or below the axillary level and distal to the shoulder. Infusates should be at or near normal serum osmolality.
- Do not use catheter for high-pressure injection. Syringes smaller than 10mL, and mechanical higher-pressure injectors can generate pressures capable of rupturing the catheter (greater than 40psi).
- Clinicians must be familiar with and trained in the placement, maintenance and use of PICC and /or midline catheters.
- Advancement of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardial erosion, or cardiac tamponade.
- Solutions containing high concentrations of alcohol can temporarily weaken the structure of polyurethane material. For polyurethane catheters, minimize exposure to alcohol solutions.
- Do not expose the catheter to acetone or acetone/alcohol solutions.
- Do not use force to advance the catheter. Resistance could indicate a vein obstruction or malposition of the catheter.
- Do not use force to remove the stylet. Resistance can damage the integrity of the catheter and the stylet.
- Do not use force to flush the catheter if resistance is met.
- Do not hold the catheter with forceps while removing the stylet.
- Do not use hemostats or clamps on catheter or hub connection.
- Do not place tape strips over the catheter tubing. This will compromise the strength and integrity of the tubing.
- Single lumen polyurethane and all silicone catheters should be secured with tape or a manufactured catheter stabilization device. These devices should not be sutured.

PRECAUTIONS

- This device is intended for use by trained medical personnel and must only be used by qualified personnel who are familiar with the technique.
- Sterilized by ethylene oxide. Do not re-sterilize.
- Carefully read and follow all instructions prior to use.
- Do not cut stylet. If cut, the stylet can potentially damage the catheter and harm the patient.
- Follow sterile technique per institutional policy. Sterile technique and proper skin preparation are essential for proper use. The use of maximal barrier precautions during catheter insertion is recommended to help reduce the risk of catheter-related infection, including the use of sterile gown, mask, cap, sterile gloves, and large sterile drapes to cover the body.
- Observe Standard/Universal precautions for all patients to reduce the potential exposure to bloodborne pathogens.
- If suturing is an option available through institutional guidelines, suture only through eyelets provided on dual lumen polyurethane catheters. Never suture around catheter tubing.

DIRECTION FOR USE

1. Identify Appropriate Vein
2. Position and Measure Patient
 - For PICC Placement Upper Extremity: Position the arm at a 90° angle, and measure from the planned insertion site to the right clavicular head, then down to the third palpated intercostal space for superior vena cava (SVC) placement.
 - For PICC Placement Lower Extremity: Measure from planned insertion site to the right of the umbilicus and up to the xiphoid process for inferior vena cava (IVC) placement.
 - For MIDLINE Placement: Measure from the planned insertion site to the appropriate tip position.
3. Setup Sterile Field
 - Open package to expose the contents of the tray.
 - Setup and prepare the necessary equipment using aseptic technique.
4. Prep Site and Drape Patient
 - Prep the anticipated puncture site according to institutional policy.
 - Drape patient.
5. Prime and Trim Catheter
 - Attached syringe filled with sterile normal saline, and prime catheter.
 - Peel open catheter's protective sleeve to the length measured in step 2.
 - If present, withdraw the stylet one-half centimeter shorter than the length measured in step 2.
 - If necessary, trim the distal end of the catheter squarely to the length measured in step 2.
6. Apply Tourniquet
7. Perform Venipuncture and Withdraw Introducer
 - Use introducer per Instructions for Use enclosed with introducer.
8. Advance Catheter
 - Peel-away-cannula may be removed at any point in the insertion process based on clinician preference.
 - Continue to advance the catheter slowly in small increments.
 - For SVC: As the catheter tip passes the shoulder, have the patient turn head toward the insertion site with chin on shoulder to help prevent possible malposition into the jugular vein.
 - Complete catheter advancement to the desired tip position.
9. Remove Stylet (if present)
 - Extend the catheter tubing and extension in a straight line with the vein.
 - While stabilizing the catheter by holding the oval junction or Y-site, grasp the stylet handle and withdraw gently and slowly. Catheter bunching indicates stylet resistance.
 - Do not use the catheter tubing for stabilization during stylet removal.
 - If resistance or bunching of the catheter is encountered, stop stylet withdrawal and allow the catheter to return to normal shape.
 - Then withdraw both the catheter and stylet together about one inch, repeat priming with saline, and try removing the stylet again.
 - Repeat this procedure until the stylet is easily removed.
 - Once the stylet is out, advance the catheter into the desired position without the stylet.
10. Aspirate and Flush
 - Use a 10 mL syringe design to flush the catheter. Do not use force to flush the catheter as a 10 mL syringe design has the potential to generate high pressure (greater than 40 psi) capable of rupturing the catheter.
 - Aspirate to visualize blood return and flush to ensure patency of all lumens.
 - Attach the accessories and flush per institutional policy.
 - Remove Y-site or T-connector and replace with Luer lock extension set or injection cap.
11. Clean Insertion Site
 - Remove drape.
 - Clean insertion site.
 - Apply skin protectant if desired.
12. Secure Catheter and Apply Sterile Dressing
 - The external portion of the catheter must be adequately secured. Any change in the length of the external portion of the catheter at the insertion site indicates that the tip location has been altered as well.
 - Dress according to institution policy and procedures.
13. Confirm catheter tip placement radiographically if indicated. Radiographic confirmation of tip placement is required for all central tip placements.

CATHETER REMOVAL

1. Remove dressing.
2. Grasp catheter near insertion site.
3. Remove slowly. NOTE: Do not use excessive force.
4. If resistance is felt, stop removal. Apply warm compress and wait 20-30 minutes.
5. Resume removal procedure.

DISPOSAL

After use, handle and dispose in accordance with hospitals policies and procedures concerning biohazard materials and waste.

STORAGE: Store at controlled room temperature.

Note: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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