ARGON MEDICAL DEVICES, INC.
 1445 Flat Creek Road
 Athens, Texas 75751 USA
 Tel: +1 (903) 675 9321
 Tel: +1 (800) 927 4669
 www.argonmedical.com

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands
 +31 70 345 8570

Argon Medical Devices UK Ltd
 Eastgate Business Centre
 Eastern Avenue
 Burton-on-Trent
 Staffordshire
 DE13 0AT

X9585764/Rev.0321A

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgaria	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols

Swedish	Symbolförklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

English

Hawkins™ Blunt Needle

Intended Use/Purpose

Hawkins™ Blunt Needle is intended for gaining percutaneous access.

Device Description

The Hawkins Blunt Needle System permit safer biopsies under CT, Ultrasound and Fluoroscopy by reducing the possibility of arterial penetration and allowing multiple passes by using a co-axial needle. It has an adjustable needle stop in the shape of a spring clip which allows the user to prevent from forward migration, effectively localizing the needle tip to the biopsy site. The outer cannula has an echogenic tip, promoting accurate placement under ultrasound guidance. It has centimeter markings on the outer cannula that provide reference for depth placement.

Indications for Use

Hawkins™ Blunt Needle is used for gaining percutaneous access for interventional radiology procedures.

Contraindications for Use

There are no contraindications; however, physicians should assess and determine if percutaneous interventional procedures are suitable for each patient.

Warnings

- Contents are supplied sterile and are intended for single use only. Frequency of use is at the discretion of the user.
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expiry date has been exceeded.

Precautions

- This device is intended for use by qualified physician with experienced in interventional radiology procedures.
- The risks associated using the Hawkins™ Blunt Needle are minimal and acceptable when used as intended by a physician.
- No limitations known to the lifetime when used in a single procedure.
- The device is intended to be mobile or portable. The Cannula could stay in the patient during CT scan, if required.

Directions for Use

- Perform a needle puncture using ultrasound or x-ray.
- When the needle is correctly positioned, remove the needle and introduce the blunt stylet and advance to intended area.
- Remove the blunt stylet.
- Introduce the 0.038" guidewire (not included).
- Remove the outer cannula.
- Insert and advance the required catheter (not included) over the guidewire. Dilate further if necessary.

Possible Complications:

- Bleeding
- Thrombus
- Venous puncture

Disposal:

After use, the devices should be placed in a sharp's container for disposal. Handle and dispose in accordance with hospitals policies and procedures concerning biohazard materials and waste.

Storage: Store at controlled temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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