Intended Use/Purpose

The Lunderquist Steel Guidewire is extremely stiff guidewire use to provide extra support during non-vascular procedure.

Device Description

The device has a flexible tip for access to delicate or difficult areas. The core material is stainless steel with no coating and available in variety of lengths and tip configurations.

Indication for use

The guidewire is only to be used for non-vascular procedures.

Duration

Transient, less than 60 minutes.

Warnings

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package appear open or if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged prior and during the procedure.
- Do not advance the wire against resistance until the cause of the resistance has been determined by fluoroscopy. Excess force against resistance may result in damage to guidewire, or catheter, or vessel perforation.
- Do not withdraw a guidewire through a needle. Straighten the guidewire in order to withdraw the needle.
- The product must only be used by qualified personnel who are familiar with the technique.

Potential Complications

Potential complications include, but are not limited to the following:

- Puncture site hematoma
- Infection
- Perforation of adjacent vessel wall, tissues/organs
- Minor pain/discomfort
- Minor bleeding

Preparation

- Make sure that the guidewire is undamaged.
- Loosen the straightener together with the guidewire.
- Inject saline with syringe in a vertical position.
- Make sure the product is compatible with other equipment.

Directions for Use

- Insert the flexible end of the guidewire into the device. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device.
- Advance the guidewire through the device.
- Position the guidewire.
- Hold the guidewire in place and withdraw the device.
- Advance the catheter over the guidewire.

Storage: Store at controlled room temperature.

Disposal: After use, this product may be a potential biohazard. Handle in a manner, which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675-9321
Tel: +1 (800) 927-4669
www.argonmedical.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
+31 70 345 8570

Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

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