|                     |  |
|---------------------|--|
| English             | The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                           |
| Português do Brasil | O glossário dos símbolos está localizado eletronicamente em <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                 |
| Czech               | Slovníček symbolů je elektronicky umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                |
| Danish              | Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                    |
| German              | Das Glossar der Symbole kann unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> abgerufen werden                         |
| Greek               | Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                |
| Spanish             | El glosario de símbolos puede consultarse electrónicamente en <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>               |
| Estonian            | Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>             |
| Finnish             | Elektroninen symbolisanasto on osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .                                |
| French              | Le glossaire des symboles est disponible au format électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> |
| Hungarian           | A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                  |
| Italian             | Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> |
| Lithuanian          | Simbolių žodyno elektroninė forma pateikta <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                  |
| Dutch               | Uitleg over de symbolen is elektronisch beschikbaar op <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                      |
| Norwegian           | Symbolforklaringen befinner seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                              |
| Polish              | Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>       |
| Portuguese          | O glossário de símbolos está localizado eletronicamente no endereço <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>         |
| Romanian            | Glosarul simbolurilor este stocat electronic la <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                             |
| Slovak              | Glosár symbolov je v elektronickej podobe k dispozícii na stránkach <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>         |
| Swedish             | Symbolforklaringen hittas elektroniskt på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                   |
| Turkish             | Semboller sözlüğü elektronik ortamda <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> adresinde yer almaktadır               |
| Arabic              | يُوفّر مسرد الرموز الإلكترونيًا على <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>   |

IFU1300

Date of Issue: 10/10/2023 Rev. D

**Intended Use/Purpose**

The device is used to facilitate the placement of catheters.

**Device Description**

AMPLATZ Guidewires are made of stainless steel core with a polytetrafluoroethylene (PTFE) coated stainless steel spring wire that helps reduce friction during procedure. It offer a safe and flexible guidewire.

**Indications for use**

AMPLATZ Guidewires are for use in facilitating the placement of catheters within the coronary and peripheral vasculature. They are also for use in the placement of non-vascular catheters.

**Duration**

Transient, less than 60 minutes.

**Warnings**

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package appears open or if the expiry date has been exceeded.
- Do not advance the wire against resistance until the cause of the resistance has been determined by fluoroscopy. Excess force against resistance may result in damage to guidewire, or catheter, or vessel perforation.
- Do not excessively torque the guidewire.
- Do not entrap or over-rotate the distal tip of the guidewire, it can lead to wire rupture.
- Do not excessively bend the guidewire, it may result in wire fracture.
- Guidewire placement should be monitored during placement and manipulation using fluoroscopic or suitable imaging method.
- Do not withdraw a guidewire through a needle. Straighten the guidewire in order to withdraw the needle.
- The product must only be used by qualified personnel who are familiar with the technique.

**Precaution**

- Twisting or entanglement of the guidewires may occur when more than one wire is being used simultaneously. This can be prevented by carefully isolating and marking the proximal end of the wires.

**Potential Complications**

Potential complications include, but are not limited to the following:

- Perforation of vessel or arterial wall
- Thrombus formation
- Puncture site hematoma
- Vasospasm
- Infection
- Vessel dissection

**Directions for Use**

- Insert the flexible end of the guidewire into the device. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device.
  - If "J" wire is used, prepare "J" tip by sliding J-straightener over tip to straighten. Insert J-tip inside J-Straightener into device hub gently.
- Advance the guidewire through the device.
- Position the guidewire.
- Hold the guidewire in place and withdraw the device.
- Advance the catheter over the guidewire.

**Preparation**

- Make sure that the guidewire is undamaged.
- Loosen the straightener together with the guidewire.
- Inject saline with syringe in a vertical position.
- Make sure the product is compatible with other equipment.

**Storage**

Store at controlled room temperature. Avoid extended exposure to light.

**Disposal**

After use, handle and dispose in accordance with end-user facility's policies and procedures concerning biohazardous materials and waste.

**NOTE:** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at [quality.regulatory@argonmedical.com](mailto:quality.regulatory@argonmedical.com) as well as to the competent health authority where the user/patient resides.

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