ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675-9321
Tel: +1 (800) 927-4669
www.argonmedical.com

EC REP **EMERGO EUROPE**
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
+31 70 345 8570

UK REP **Emergo Consulting (UK) Limited**
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

IFU9014
Date of Issue: 08/18/2023 Rev. B

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Arabic	بنوفر مصدر الرموز إلكترونياً على www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Portuguese-Brazil	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

Intended Use/Purpose

The Access Guidewire is intended for use in the placement of non-vascular percutaneous drainage catheters.

Device Description

The Access Guidewire features a high torque-control and smooth transition from rigid core to radiopaque coil tip that provides excellent manoeuvrability.

Indication for use

The product is only for non-vascular procedures.

Duration

Transient, less than 60 minutes.

Warnings

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package appear open or if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- Do not advance the wire against resistance until the cause of the resistance has been determined by fluoroscopy. Excess force against resistance may result in damage to guidewire, or catheter, or vessel perforation.
- Do not withdraw a guidewire through a needle. Straighten the guidewire in order to withdraw the needle.
- Make sure the product is compatible with other equipment.
- The product must only be used by qualified personnel who are familiar with the technique.

Potential Complications

Potential complications include, but are not limited to the following:

- Puncture site hematoma
- Infection
- Perforation of the adjacent vessel wall, tissues/organs
- Minor pain/discomfort
- Minor bleeding

Preparation

- Make sure that the guidewire is undamaged.
- Loosen the straightener together with the guidewire.
- Inject saline with syringe in a vertical position.
- Stainless Steel Wire** - Shape the tip into a J-curve using your fingers if necessary.
- Nitinol Wire**
 - Form the tip by running the coil over the edge of the coloured plastic straightener.
 - Activate the coating by saline.

Directions for Use

- Insert the flexible end of the guidewire into the device. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device.
- Advance the guidewire through the device.
- Position the guidewire.
- Hold the guidewire in place and withdraw the device.
- Advance the catheter over the guidewire.

Storage: Store at controlled room temperature.

Disposal: After use, this product may be a potential biohazard. Handle in a manner, which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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