ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669
www.argonmedical.com

EMERGO EUROPE

Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

Argon Medical Devices UK Ltd

Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT

© 2015 Argon Medical Devices, Inc.
X9585745/0221A

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolų glosarijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Latvian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Dutch	Symbolforklaringen befinnen seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols

English

SEMI-AUTOMATED / SINGLE ACTION BIOPSY DEVICE (SABD™)

Bevel and Trocar Tip Design

Intended Use/Purpose:

The Semi-Automated/Single Action Biopsy Device (SABD™) is semi-automatic biopsy instrument intended to be used to obtain soft tissue biopsies.

Device Description

Semi-Automated/Single Action Biopsy Devices (SABD™) is a lightweight, adjustable and easily maneuverable biopsy instrument. It composed of an echogenic spring-loaded biopsy needle fitted into a plastic handle to permit single handed specimen collection. The position of the needle may be visualized by X-Ray, CT or Ultrasound.

Indication for Use

To be used for soft tissue and tumor biopsy of such organs as the liver, spleen, kidney, prostate, lung, breast and lymph nodes. When used for breast biopsy, product is for diagnosis only.

Contraindications

No absolute contraindication, however, sound professional judgment should be used when performing biopsies on patients who are using anticoagulant therapy and/or patients with known hemophilia problems.

Warnings

- The SABD™ device is supplied sterile by ethylene oxide gas as a single use product.
- The device is packaged in a plastic pouch. Inspect the package integrity before use. Do not use the device if package is open or appears to be damaged or defective.
- Do not re-use, re-process, or re-sterilize the device, tray or any items provided in the tray, it has not been evaluated and it may compromise the structural integrity and subsequent patient illness, infection and/or other injury.
- Do not use if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- Verify integrity of the needle before loading instrument and after firing(s). If needle is damaged, replace the whole system (needle set and instrument).
- Functional tests are performed during manufacture. Do not try to "dry fire" the device.
- Not for use with MRI

Precautions

- Device should only be used by or under supervision of a physician who is completely familiar with the use and limitations of biopsy sampling.
- Use sound professional judgment when deciding whether it is advisable to use the product with the aid of visualization modalities such as ultrasound, computerized tomography (CT), fluoroscopy or direct visualization.
- Ensure that no critical structures are beyond tip of device within needle stroke length (throw).


NOTE: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.


Potential Complications/Risk:

Possible complications may include but are not limited to the following. When used as intended by a physician the benefits of using the SABD Instrument outweigh the risks associated with the use of the device.

- Hemorrhage
- Infection
- Adjacent Tissue Injury
- Pain
- Hemoptysis
- Hemothorax
- Organ or Vessel Perforation
- Pneumothorax
- Air Embolism

Preparation and Instructions for Use:

- I. Ready the device by pulling back on the plunger (1) until a firm click is felt, indicating that the cutting cannula is locked into "Ready" position.
- II. With the plunger (1) fully retracted, so that specimen notch is completely covered by the cutting cannula, advance Device Tip (2) proximal to the area to be biopsied. **Note: Do not advance Plunger (1) until stylet is in position.**
- III. While maintaining the Device Tip (2) position, advance the Plunger (1) with thumb to expose specimen notch within area to be biopsied. Activate cutting cannula by fully depressing Plunger (1) with thumb to capture tissue within specimen notch.
- IV. Remove the device from the biopsy area. To retrieve tissue specimen, pull back Plunger (1) until a firm click is felt, indicating cutting cannula is locked into position. Advance Plunger (1) forward to expose tissue specimen within notch. Retrieve tissue from specimen notch.

Handling and Disposal:

Handle in a manner which will prevent accidental puncture, percutaneous injuries and exposure to bloodborne pathogens. After use, this product may be a potential biohazard. Dispose in accordance with applicable laws and regulations.

Storage:

Store at standard ambient temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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