**INSTRUCTIONS FOR USE**  
**11 Gauge Bone Access Needle**  
**11 Gauge Bone Access Needle, Trocar**  
**11 Gauge Bone Access Needle, Bevel**

<https://www.argonmedical.com/resources/product-information>

**DISTRIBUTED BY:**

Medtronic Sofamor Danek USA, Inc.  
1800 Pyramid Place.  
Memphis, TN 38132  
U.S.A.  
Telephone 800 933 2635 (In U.S.A.)  
901 396 3133 (Outside of U.S.A.)  
Fax 901 396 0356

**MANUFACTURED BY:**

**Argon Medical Devices Inc.**  
1445 Flat Creek Road  
Athens, TX 75751 USA  
Tel: +1 (903) 675 9321  
Tel: +1 (800) 927 4669  
[www.argonmedical.com](http://www.argonmedical.com)

**EC REP**  
EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands  
+31 70 345 8570

**UK REP**  
**Emergo Consulting (UK) Limited**  
c/o Cr 360 - UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
United Kingdom

MEDTRONIC is a registered trademark of Medtronic Inc.


Figure 1, Figura 1, Figure 1, Figure 1, Kuva 1, Figura 1, Figure 1, Figure 1, Figure 1, Figure 1, Obrázek 1, Figure 1, Obr. 1, Slika 1, 1. ábra, 1. attéls, Joonis 1, Εικόνα 1, 1 pav., Abbildung 1, Figura 1, Rysunek 1, Figura 1, Şekil 1

|   |  |  |
|---|--|--|
|   |  |  |
| Diamond; Diamante, Forme de diamant, Diamant, Timanti, Diamante, Diamant, Ruitvormig, Kosoštercová, Rombisk, Kosoštercová, Diamant, Gyémánt, Rombveida gals, Romb, Ρομβοειδής, Rombo forma, Diamant, Diamante, Diamant, Romb, Elmas | Trocar; Trocar, Trocart, Trokar, Trokaari, Trocarte, Trokar, Trocart, Trokarová, Troakar, Trokarová, Trokar, Trokár, Trísstürveida gals, Trokaar, Trokáp, Troakaro forma, Trokar, Trocar, Trokar, Trocar, Trokar | Bevel; Bisel, Biseau, Skrá, Viistokarki, Bisel, Affasning, Afgeschuind, Zkosená, Fasad, Šikmá, Poševni rob, Ferde, Silpais gals, Kaldots, Λοβήτητο, Kúginé forma, Schrägung, Smussata, Skos, Bizou, Eğimli |

Figure 2; Figura 2, Figure 2, Figure 2, Kuva 2, Figura 2, Figure 2, Figure 2, Obrázek 2, Figure 2, Obr. 2, Slika 2, 2. ábra, 2. attéls, Joonis 2, Εικόνα 2, 2 pav., Abbildung 2, Figura 2, Rysunek 2, Figura 2, Şekil 2

**Tip Configuration; Configuración de punta, Configuration de l'embout, Spisskonfigurasjon, Kärjen kokoonpano, Configuração da ponta, Spidskonfiguration, Tipconfiguratie, Konfigurace hrotu, Spetskonfiguration, Konfigurácia špičky, Konfiguracija konice, Hegykonfiguráció, Galu konfigurácia, Otsa konfiguratsioon, Διαμόρφωση άκρου, Galuiko konfiguracija, Spitzkonfiguration, Configurazione della punta, Konfiguracja końcówki, Configurația vârfului, Uç Yapılandırması**

|            |  |
|------------|--|
| English    | The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                           |
| Spanish    | El glosario de símbolos puede consultarse electrónicamente en <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>               |
| French     | Le glossaire des symboles est disponible au format électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> |
| Norwegian  | Symbolforklaringen finnes seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                |
| Finnish    | Elektroninen symbolisanasto on osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .                                |
| Portuguese | O glossário de símbolos está localizado eletronicamente no endereço <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>         |
| Danish     | Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                    |
| Dutch      | Uitleg over de symbolen is elektronisch beschikbaar op <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                      |
| Czech      | Slovníček symbolů je elektronicky umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                |
| Swedish    | Symbolförklaringen hittas elektroniskt på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                   |
| Slovak     | Glosár symbolov je v elektronickej podobe k dispozícii na stránkach <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>         |
| Slovenian  | Slovar simbolov se nahaja v elektronski obliki na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .                         |
| Hungarian  | A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                  |
| Latvian    | Simbolu glosārijs ir atrodams elektroniskā veidā tīmekļa vietnē <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>             |
| Estonian   | Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>             |
| Greek      | Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                |
| Lithuanian | Simbolių žodyno elektroninė forma pateikta <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                  |
| German     | Das Glossar der Symbole kann unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> abgerufen werden                         |
| Italian    | Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> |
| Polish     | Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>       |
| Romanian   | Glosarul simbolurilor este stocat electronic la <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                             |
| Turkish    | Semboller sözlüğü elektronik ortamda <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> adresinde yer almaktadır               |

**ENGLISH**

**INSTRUCTIONS FOR USE**

**11 Gauge Bone Access Needle**  
**11 Gauge Bone Access Needle, Trocar**  
**11 Gauge Bone Access Needle, Bevel**

**INTENDED USE/PURPOSE**

Bone Access Needles are designed to provide access to bone for biopsy procedure.

**INDICATIONS FOR USE**

The 11 Gauge Bone Access Needles are intended only for percutaneous access to bone.

**DEVICE DESCRIPTION**

The Bone Access Needle is available in size 11ga with the following feature designs.

**11 Gauge Bone Access Needle (1) (See Figure 1)**

- A. Stylus
- B. Stylet
- C. Direction Marker
- D. Tip

**Tip configuration (See Figure 2)**

- Diamond
- Trocar
- Bevel

**WARNINGS**

Breakage of the device may require intervention or retrieval.

This device was designed, tested and manufactured for single patient use only. When used repeatedly in the same patient, the device should be inspected for damage or wear after each core sample is taken. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness and/or injury. DO NOT reuse, reprocess or re-sterilize this device.

**PRECAUTIONS**

- It is important to read the Instructions For Use and these precautions carefully prior to device operation.
- Use the device prior to the use by date noted on the package.
- Do not use if package is opened or damaged because product integrity including sterility may be compromised.
- Do not use damaged products. Before use, inspect the device and packaging to verify that no damage has occurred.
- Do not use this product if you have not been properly trained in its use. Physicians using the device should be familiar with the physiology and pathology of the selected anatomy, and be trained in the performance of the chosen surgical technique.
- This device should be manipulated only while under fluoroscopic observation with radiographic equipment that provides high quality images.
- Access to the vertebral body via the pedicle requires a minimum pedicle width of 5 mm. Insertion of the instrument requires specific knowledge of the dimensions of the site of insertion as assessed by MRI, CT, or other imaging method.
- Never use electric power (or any other alternative power sources) in conjunction with the 11 Gauge Bone Access Needle. Use only manual power when using the 11 Gauge Bone Access Needle.
- Reconditioning, refurbishing, repair, modification, or resterilization of the device to enable further use is expressly prohibited.

**Potential Complication/Side Effects:**

Complications potentially associated with use of the device include:

- Nerve injury including puncture of the spinal cord or nerve roots potentially resulting in radiculopathy, paresis or paralysis
- Thrombus or other materials resulting in symptomatic pulmonary embolism or other clinical sequelae
- Hemothorax or pneumothorax

- Deep or superficial wound infection
- Unintended puncture wounds including vascular puncture and dural tear
- Bleeding or hemorrhage
- Hematoma
- Pain

#### DIRECTIONS FOR USE

1. Make a skin incision with a scalpel blade over the selected area for initial access.
2. Hold the 11 Gauge Bone Access Needle with the handle in the palm of the hand.
3. Using manual control and appropriate imaging guidance, advance the needle with gentle but firm pressure through the incision to the selected bone surface. Rotate the needle in an alternating clockwise and counterclockwise motion for insertion into bone to the desired depth.
4. To remove the stylet from the access needle stylus, hold the handle of the access needle stylus while rotating the stylet handle counterclockwise 90° and pull it straight out from the access needle.
5. Advance a guide pin (not included) within the lumen of the access needle stylus, to the desired depth using manual control and appropriate imaging guidance. Hold the guide pin in place and remove the access needle stylus using an alternating clockwise and counterclockwise rotation.
6. A working channel can now be created using the bone access tools.
7. To re-insert stylet into the access needle stylus, insert it into the access needle and rotate the stylet clockwise to lock the handles. The handles will be aligned when locked.

#### STERILIZATION

Sterilized with Ethylene Oxide.

#### HOW SUPPLIED

The 11 Gauge Bone Access Needle is supplied sterile in a peel-open package. In the event of damage to the sterile packaging, do not use and notify the manufacturer.

#### STORAGE

The instrument should be stored in their original shipping materials. Proper care should be taken to ensure that the instruments will not be damaged. Store in controlled room temperature.

#### DISPOSAL

After use, this product may be potential biohazard. Handle and dispose of in accordance with accepted medical practice and all applicable laws and regulations.

**NOTE:** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

#### LIMITATION OF LIABILITY

IN NO EVENT SHALL MEDTRONIC BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE 11 GAUGE BONE ACCESS NEEDLE, BASED UPON BREACH OF CONTRACT (INCLUDING BREACH OF WARRANTY).

#### REQUESTS FOR INFORMATION

For further information, please contact Customer Service, Medtronic Sofamor Danek USA, Inc. 1800 Pyramid Place, Memphis, TN 38132 U.S.A.; phone: 800 933 2635, toll free in the U.S.: 901 396 3133, fax: 901 396 0356.

© 2011 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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© 2011 Medtronic Sofamor Danek USA, Inc. Wszelkie prawa zastrzeżone.