ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road, Athens, Texas 75751 USA
Tel: +1 (903) 675-9321
Tel: +1 (800) 927-4669
www.argonmedical.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
+31 70 345 8570

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
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols
Portuguese-Brazil	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Swedish	Symbolförklaringen hittas elektroniskt på www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols
Norwegian	Symbolforklaringen finnes elektronisk på www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese – Portugal	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Czech	Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusán itt található: www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols
Arabic	www.argonmedical.com/symbols يوفر مسرد الرموز إلكترونياً على

ENGLISH

V•Stick™ Vascular Access Set

Intended Use: V•Stick™ Vascular Access Sets facilitate fine needle venous or arterial puncture and tissue dilation for subsequent positioning of a 0.035" or 0.038" guidewire.

Device Description: The V•Stick™ Vascular Access Set consists of a 4F or 5F coaxial introducer set (with a 3F standard or 3F stiff dilator), a 21ga entry needle (with an echogenic or non-echogenic tip), and a 0.018" nitinol or stainless steel guidewire with a floppy radiopaque tip.

Indication for Use: The V•Stick™ Vascular Access Sets is intended for use in the introduction and placement of guidewires and/or catheters.

Contraindications: This device is contraindicated for use on patients who cannot receive antiplatelet and/or anticoagulant therapy. Physician judgment should be used in determining when percutaneous interventional procedures are contraindicated.

Complications/Risk:

- Bleeding
- Hematoma
- Occlusion
- Infection
- Thrombosis
- Phlebitis

Warning/Precaution:

- The product must only be used by qualified personnel who are familiar with the technique.
- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- The V•Stick™ Vascular Access Set is sterilized by ethylene oxide gas. Content is sterile and non-pyrogenic in unopened and undamaged package. Do not use device if package has been damaged or has been opened.

Caution:

1. The guidewire should not be withdrawn through the needle. Damage or shearing of the guidewire may occur.
2. The 0.018" guidewire must be in appropriate position in the vessel before advancing the coaxial introducer system.

Direction for Use:

1. Perform a skin incision under local anesthesia.
2. Do a fine needle puncture of the vessel using the 21ga needle.
3. Introduce the 0.018" guidewire, making sure to loosen the straightener together with the guidewire on the circular dispenser.
4. Remove the needle, place the fully assembled coaxial introducer system over the guidewire, and advance it to dilate the vessel.
5. Loosen the inner dilator, making sure that the tip of the sheath introducer has been properly placed inside the vessel.
6. Remove the inner dilator and the 0.018" guidewire.
7. Introduce a 0.035" guidewire or a 0.038" guidewire (not included in the set) through the sheath introducer.
8. Remove the sheath introducer.
9. Advance the required catheter over the guidewire.

Storage: Store at controlled room temperature.

Disposal: After use, this product may be a potential biohazard. Handle in a manner, which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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