Fine Needle Aspiration (FNA) Biopsy Needles (Chiba, Franseen, Greene, MaxiCELL™, MCN Core, Spinal, Trocar, Westcott)

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road, Athens, Texas 75751 USA
Tel: +1 (903) 675-9321;
Tel: +1 (800) 927-4669
www.argonmedical.com

EC REP **EMERGO EUROPE**
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
+31 70 345 8570

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UK REP **Emergo Consulting (UK) Limited**
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

X9585392-0321
Date of Issue: 05/18/2023 Rev. B

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът за символите в електронен вариант се намира на www.argonmedical.com/symbols
Portuguese	O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλῶσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sõmbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Głosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Russian	Электронный глоссарий символов расположен по адресу www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols

Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Arabic	يتوفر مسرد الرموز الإلكتروني على www.argonmedical.com/symbols

EN - ENGLISH

Intended Use/Purpose

Fine Needle Aspiration (FNA) Biopsy Needles are intended to harvest aspirate for cytological evaluations or specimens for histological evaluations.

Device Description

The FNA Needles consist of a stainless-steel stylet inside of a stainless-steel cannula with molded plastic hubs. A locking spring is provided on the cannula and is used to secure the needle against the skin and prevent forward movement. Each needle has numerically ordered centimeter markings on the outer cannula to provide reference for depth placement. The Biopsy Needles are available in a wide variety of tip configurations, gauge sizes, and overall lengths.

Indications for Use:

The FNA Needles are indicated to be used in percutaneous soft tissue (pleura cavity and associated organs, abdominal cavity and associated organs, etc.) biopsies effected via the insertion of the device by a qualified physician, under guidance (CT, Ultrasound, MRI, Sequential Mammography, or Fluoroscopy), and manipulated in such a manner that a tissue specimen is retained within a collection orifice on the needle.

Contraindications for Use: Physician judgment is required when considering biopsy on patients with bleeding disorders or receiving anti-coagulant medications who may have the potential of an increased risk of bleeding complications.

Possible Complications:

Depending on the location of the biopsy, patient can exhibit the following rare complications:

- Minor Pain/discomfort
- Bleeding
- Hematoma
- Infection
- Post FNA infarction of the target tissue

Warnings

- The device is designed, tested and manufactured for single use only.
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Do not reuse, reprocess or re-sterilize the device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not attempt to use the device on hard tissue.
- Do not use excessive force, it may result in damage/breakage of the device.
- Do not bend the needle excessively, it may result in damage/breakage of the device.
- Verify the integrity of the needle before inserting. If needle is damaged, replace the needle.

Precautions:

- These devices should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of a FNA procedure.
- The device is intended for continuous use for less than 60 minutes.
- Device can be used during a single procedure and will allow multiple biopsies to be taken. No known limitations to the lifetime when used in a single procedure.
- Patient movement may affect outcome of procedure relating to quality of sample. Minimal risk to patient.
- The Needle can be used with a syringe to create negative pressure (aspiration).

NOTE: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

BIOPSY PROCEDURE:

1. Using sterile technique, prepare the biopsy puncture site.
2. Insert the needle through the skin with the stylet in place, and advance the needle, typically under imaging guidance, directly into the lesion.
For MaxiCELL™ - Insert the MaxiCELL needle through the skin with the stylet in place, and advance the needle, typically under imaging guidance, directly into the lesion with the scraping aperture at the 12 o'clock position. The scraping aperture is visible on direct inspection of the needle prior to insertion and there is a notch on the stylet which marks the position of the aperture.
3. Remove the stylet and attach a syringe.
4. While applying suction, move the needle up and down rapidly several times within the lesion.
For MaxiCELL™ - While applying suction, gently angle the needle shaft over forefinger to create a convex needle shape which opens the side port of the needle to engage the scraping aperture in the lesion and move the needle up and down rapidly several times within the lesion.
5. When the aspiration procedure is complete, maintain gentle suction and withdraw the needle.
For MaxiCELL™ - When the aspiration procedure is complete, maintain gentle suction and angle the MaxiCELL™ needle under the forefinger to create a concave shape relative to the scraping aperture and thus disengage the scraping aperture to allow for easy withdrawal.
6. The sample (aspirate or tissue fragments) are now contained within the needle and ready to be deposited on a slide or in formalin for cytological analysis.
7. Remove the needle.
8. Apply pressure to puncture site to control possible bleeding.

Disposal

After use, this product may be potential biohazard. Handle in a manner, which will prevent accidental sharps puncture. Dispose in accordance with hospital policies and procedures for waste and biohazard materials.

Storage

Store at controlled room temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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