**ARGON MEDICAL DEVICES, INC.**

1445 Flat Creek Road  
Athens, Texas 75751 USA  
Tel: +1 (903) 675 9321  
Tel: +1 (800) 927 4669  
www.argonmedical.com

**EMERGO EUROPE**  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands  
+31 70 345 8570

**Emergo Consulting (UK) Limited**  
c/o Cr 360 - UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
United Kingdom

© 2015 Argon Medical Devices, Inc.  
X9585745-0221

Date of Issue: 08/14/2023 Rev. B

<https://www.argonmedical.com/resources/product-information>

|                   |                                                                                                                                                            |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>English</b>    | The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                           |
| <b>Bulgarian</b>  | Речникът на символите е наличен в електронен вид на адрес <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                   |
| <b>Czech</b>      | Slovník symbolů v elektronické formě je umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                          |
| <b>Danish</b>     | Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                    |
| <b>German</b>     | Das Glossar der Symbole kann unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> abgerufen werden                         |
| <b>Greek</b>      | Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                |
| <b>Spanish</b>    | El glosario de símbolos puede consultarse en línea en la dirección <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>          |
| <b>Estonian</b>   | Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>             |
| <b>Finnish</b>    | Elektroninen symbolisanasto on osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                  |
| <b>French</b>     | Le glossaire des symboles est disponible au format électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> |
| <b>Croatian</b>   | Pojmovnik simbola nalazi se u elektroničkom obliku na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                       |
| <b>Hungarian</b>  | A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                  |
| <b>Italian</b>    | Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> |
| <b>Lithuanian</b> | Simbolų glosarijs ir atrodams elektroniską veidą vietnė <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                     |
| <b>Latvian</b>    | Simbolių žodyno elektroninė forma pateikta <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                  |
| <b>Dutch</b>      | Symbolforklaringen befinnen seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                              |
| <b>Polish</b>     | Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>       |
| <b>Portuguese</b> | O glossário de símbolos está localizado eletronicamente no endereço <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>         |
| <b>Romanian</b>   | Glosarul simbolurilor este stocat electronic la <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                             |
| <b>Slovak</b>     | Glosár symbolov je v elektronickej podobe k dispozícii na stránkach <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>         |
| <b>Slovenian</b>  | Slovar simbolov se nahaja v elektronski obliki na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                           |
| <b>Swedish</b>    | Symbolforklaringen hittas elektroniskt på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>                                   |

**English****SEMI-AUTOMATED / SINGLE ACTION BIOPSY DEVICE (SABD™)**

Bevel and Trocar Tip Design

**Intended Use/Purpose:**

The Semi-Automated/Single Action Biopsy Device (SABD™) is semi-automatic biopsy instrument intended to be used to obtain soft tissue biopsies.

**Device Description**

Semi-Automated/Single Action Biopsy Devices (SABD™) is a lightweight, adjustable and easily maneuverable biopsy instrument. It composed of an echogenic spring-loaded biopsy needle fitted into a plastic handle to permit single handed specimen collection. The position of the needle may be visualized by X-Ray, CT or Ultrasound.

**Indication for Use**

To be used for soft tissue and tumor biopsy of such organs as the liver, spleen, kidney, prostate, lung, breast and lymph nodes. When used for breast biopsy, product is for diagnosis only.

**Contraindications**

No absolute contraindication, however, sound professional judgment should be used when performing biopsies on patients who are using anticoagulant therapy and/or patients with known hemophilia problems.

**Warnings**

- The SABD™ device is supplied sterile by ethylene oxide gas as a single use product.
- The device is packaged in a plastic pouch. Inspect the package integrity before use. Do not use the device if package is open or appears to be damaged or defective.
- Do not re-use, re-process, or re-sterilize the device, tray or any items provided in the tray, it has not been evaluated and it may compromise the structural integrity and subsequent patient illness, infection and/or other injury.
- Do not use if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- Verify integrity of the needle before loading instrument and after firing(s). If needle is damaged, replace the whole system (needle set and instrument).
- Functional tests are performed during manufacture. Do not try to "dry fire" the device.
- Not for use with MRI

**Precautions**


- Device should only be used by or under supervision of a physician who is completely familiar with the use and limitations of biopsy sampling.
- Use sound professional judgment when deciding whether it is advisable to use the product with the aid of visualization modalities such as ultrasound, computerized tomography (CT), fluoroscopy or direct visualization.
- Ensure that no critical structures are beyond tip of device within needle stroke length (throw).


**NOTE:** These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

**Potential Complications/Risk:**

Possible complications may include but are not limited to the following. When used as intended by a physician the benefits of using the SABD Instrument outweigh the risks associated with the use of the device.

- Hemorrhage
- Infection
- Adjacent Tissue Injury
- Pain
- Hemoptysis
- Hemothorax
- Organ or Vessel Perforation
- Pneumothorax
- Air Embolism

**Preparation and Instructions for Use:**

- I. Ready the device by pulling back on the plunger (1) until a firm click is felt, indicating that the cutting cannula is locked into "Ready" position.
- II. With the plunger (1) fully retracted, so that specimen notch is completely covered by the cutting cannula, advance Device Tip (2) proximal to the area to be biopsied. **Note: Do not advance Plunger (1) until stylet is in position.**
- III. While maintaining the Device Tip (2) position, advance the Plunger (1) with thumb to expose specimen notch within area to be biopsied. Activate cutting cannula by fully depressing Plunger (1) with thumb to capture tissue within specimen notch.
- IV. Remove the device from the biopsy area. To retrieve tissue specimen, pull back Plunger (1) until a firm click is felt, indicating cutting cannula is locked into position. Advance Plunger (1) forward to expose tissue specimen within notch. Retrieve tissue from specimen notch.

#### Handling and Disposal:

Handle in a manner which will prevent accidental puncture, percutaneous injuries and exposure to bloodborne pathogens. After use, this product may be a potential biohazard. Dispose in accordance with applicable laws and regulations.

#### Storage:

Store at standard ambient temperature.

**NOTE:** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at [quality.regulatory@argonmedical.com](mailto:quality.regulatory@argonmedical.com) as well as to the competent health authority where the user/patient resides.

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