
Introducer Products

ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 800-927-4669;
Tel: +1 903-675-9321
www.argonmedical.com

EMERGO EUROPE

Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
+31 70 345 8570

Emergo Consulting (UK) Limited

c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

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INSTRUCTIONS FOR USE

INTENDED USE/PURPOSE

The Introducer devices are designed to provide a percutaneous conduit for the atraumatic passage of other intravascular devices.

DESCRIPTION

Introducer Products are stand-alone components, standard kits, or convenience kits built to user specifications. These products facilitate access to a vein or artery, allowing various catheters to be passed through for local or central diagnostic or therapeutic interventions.

These instructions for use cover all families of Introducer products and associated accessories including the Septishield Protective Catheter Sleeve.

INDICATIONS FOR USE

The introducers are indicated for use in arterial and venous procedures requiring percutaneous introduction of intravascular devices.

Radiopaque banded introducers identify location of Introducer distal tip for positioning.

The Argon Guiding Catheter Introducer is designed to provide a pathway through which a PTCA Guiding Catheter may be percutaneously introduced into the vasculature. These introducers also provide a mechanism through which catheter exchanges may be made.

WARNINGS

1. Contents of these sets have been sterilized with ethylene oxide gas. They are only for single patient, single-use applications.
2. Carefully read these instructions before using these products. When using sets incorporating other manufacturers' components, also read their current Instructions For Use.
3. Physician should be familiar with insertion techniques of introducers using a guide wire (Seldinger/Modified Seldinger).
4. **CARE SHOULD BE EXERCISED DURING INSERTION, USE OR REMOVAL OF THE DEVICE TO PREVENT ASPIRATION OF AIR INTO THE VASCULATURE.**
5. Withdrawal, pullback, or manipulation of the guide wire when resistance is met may cause guide wire damage, breakage, or embolization.
6. Rapid withdrawal of the catheter or dilator through the hemostasis valve may cause misalignment of the valve gasket assembly, causing bleedback through the valve. **SHOULD THIS OCCUR, RESETTING OF THE VALVE MAY BE ACCOMPLISHED BY GENTLE PRESSURE INTO THE VALVE WITH THE TIP OF A DILATOR OR CATHETER.**
7. Inspect the packaging and product for damage prior to use. Do not use if the package is damaged as sterility may be compromised.
8. Do not use if expiry date has been exceeded.
9. Do not continue to use if any of the component are damaged during the procedure.
10. Care should be taken that acetone solutions or isopropyl alcohol do not come in contact with introducer hub. Such substances could weaken the hub and cause possible leakage.
11. Tighten all connections prior to use **WITHOUT OVERTIGHTENING**. Overtightening can cause damage to the components.
12. Periodically check all connections for tightness.
13. This device is not intended to be used as a catheter for performing angiographic studies.
14. Injection of dye through any component should be via hand injection only.

This Device is designed and intended for ONE-TIME USE ONLY. Do not resterilize and/or reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. This device is extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

15. For products containing the  symbol on the label: This product contains Di(2-ethylhexyl)phthalate (DEHP) which has been shown to cause reproductive harm in male neonates, pregnant women carrying male fetuses, and peripubertal males. The following procedures have been identified as posing the greatest risk for DEHP exposure: exchange transfusion in neonates, total parenteral nutrition (TPN) in neonates (with lipids in polyvinylchloride (PVC) bag), multiple procedures in sick neonates (high cumulative exposure), heart transplantation or coronary artery bypass graft surgery (aggregate dose), and massive infusion of blood into trauma patient. It is recommended that DEHP-free medical products be considered when these procedures are to be performed on male neonates, pregnant women who are carrying male fetuses, and peripubertal males.

CAUTIONS

1. Percutaneous introducer sheaths should not remain indwelling without internal support of the sheath wall.
2. Care should be exercised during insertion, manipulation and withdrawal of the catheter through the hemostasis valve. Do not insert blunt rigid objects through the hemostasis valve. The integrity of the hemostasis valve may be compromised with the use of excessive force, rapid withdrawal or rigid blunt accessories. If resistance is encountered, **DO NOT FORCE**.
3. If using the Tuohy-Borst Introducer, **DO NOT** tighten the adapter without internal support.
4. If repositioning of the sheath suture tab is necessary, ensure that internal support is provided by a vessel dilator prior to insertion to avoid twisting the sheath.
5. Upon removal of a catheter or other medical device, it is recommended that an Obturator be immediately placed in the sheath.
6. Upon withdrawal of a device, the sideport should be aspirated in order to remove any clotted material that may have accumulated. This should be accomplished before the insertion of another medical device.
7. Depth markings on outside diameter of Introducer tubing are references only for visual confirmation of introducer insertion distance. NOTE: Depth markings are not intended for use as a finite measurement.
8. If using the Septishield II Protective Catheter Sleeve, care should be exercised not to overtighten the Tuohy-Borst mechanism as this may occlude the lumen(s) of the catheter.

PRECAUTIONS

1. Simultaneous advancement of sheath and dilator with a rotating motion is recommended to prevent sheath and vessel damage.
2. Maintain a firm grip on the guide wire at all times.
3. At no time should the introducer or guide wire be advanced or withdrawn when resistance is met without determining the cause by fluoroscopic examination.
4. Introducer devices are intended to be used by trained medical personnel in a clinical setting.

PREPARATION AND PROCEDURES FOR USE

Percutaneous Introducers

The following instructions are supplied for information about the techniques of percutaneous catheter introduction. Procedural details may be altered at the discretion of the physician.

1. Prepare and drape the skin at the intended puncture site; achieve local anesthesia as necessary.
2. Prepare the Introducer assembling the dilator(s) into the introducer and locking the luer assemblies. Ensure all connections are tight without overtightening. Overtightening may result in damage to components.
3. If using the detachable Hemostasis Valve, attach proximal end of valve to sheath by twisting on the luer lock connector.
4. Locate vessel with an appropriate sized needle and syringe.
5. Enter vessel using thin walled needle or O-T-N catheter. Remove locator needle. If O-T-N catheter is used, advance catheter over needle into vessel and remove needle when good flow verifies placement of catheter tip within the vessel. **IF INSERTION IS UNSUCCESSFUL, WITHDRAW COMPLETE ASSEMBLY AS A UNIT. DO NOT ADVANCE NEEDLE INTO CATHETER OR ATTEMPT TO WITHDRAW CATHETER BACKWARDS OVER THE NEEDLE WHILE CATHETER IS IN THE PATIENT.**
6. Straighten "J" tip of guide wire with plastic insertion sleeve.
7. Insert guide wire into the O-T-N catheter or thin walled needle and gently advance it to the desired length. It may be necessary to gently rotate the "J" tip for successful advance of the guide wire. Avoid

vigorous manipulation to prevent damage to vessel or shearing of the guide wire tip. **MAINTAIN A FIRM GRIP ON GUIDE WIRE AT ALL TIMES.**

8. If an obstruction is met that cannot be passed, remove needle and guide wire together and select another introduction site. **DO NOT ATTEMPT TO WITHDRAW GUIDE WIRE BACKWARDS THROUGH NEEDLE OR CATHETER AS THIS MAY RESULT IN SHEARING GUIDE WIRE OR DAMAGING CATHETER.**
9. When guide wire is advanced to desired location remove O-T-N catheter or needle proximally.
10. Thread assembled sheath/dilator until over guide wire and advance to puncture site. Enlarge puncture site with a small nick if needed, insert sheath/dilator unit into vessel with a slight rotary motion.
11. Remove the dilator and guide wire together leaving sheath in vessel.
12. If using an 8F or 9F Guiding Catheter Introducer, the guide wire and the Inner Dilator may be removed providing a pathway through which a .063 inch guide wire may be introduced. The Secondary Dilator may then be removed leaving the sheath and .063 inch guide wire in the vasculature. A PTCA Guiding Catheter may then be introduced over the .063 inch guide wire and maneuvered into the vasculature through the sheath. Alternately, the Inner Dilator, Secondary Dilator and guide wire can be removed and the large (.063 inch) guide wire and PTCA Guiding Catheter introduced simultaneously into the sheath.
13. Suture sheath in place using suture tab or hub. **DO NOT PLACE SUTURE ON SHEATH TUBING AS THIS MAY RESTRICT FLOW OR DAMAGE TUBING.**
14. Connect hemostasis valve securely to sheath when utilizing an introducer with removable valve.
15. A catheter one-half to one French size smaller than introducer sheath is recommended to allow adequate flow through introducer side-port.
16. Tuohy-Borst adapters should be gently hand-tightened to prevent blood reflux and/or catheter migration. Overtightening may compromise indwelling catheter lumen.
17. Care should be exercised not to pull the guide wire or catheter to extreme angles while the device is through the hemostasis valve, as it will distort the valve leaflets and promote leakage.
18. Follow hospital protocol for puncture site dressing and maintenance.

Septishield II Protective Catheter Sleeve

The following aseptic technique is supplied only to explain the use of the Septishield II Protect Catheter Sleeve. At the discretion of the physician the procedural details may be altered in accordance with a preferred technique or protocol.

1. Upon removal from pouch, please check the Tuohy-Borst mechanism of the proximal hub to insure that the lumen is sufficiently open to accept the catheter. Accepts maximum 9F (0.118" OD) catheter. Insert the catheter through the proximal hub (blue) of the Septishield and advance Septishield well onto the catheter.
2. After placement of Argon's sheath introducer and hemostasis valve has been completed, introduce the catheter to proper position.
3. Advance the Septishield on the catheter until the distal hub meets the hemostasis valve. The Septishield may now be connected to the hemostasis valve by simply luer locking both fittings together.
4. Extend protective inner sleeve to desired length.
5. Once inner sleeve is extended to desired length, secure onto catheter by adjusting the Tuohy-Borst mechanism of the proximal hub.

Possible Complications/Risk:

- Injury at site
- Bleeding, infection
- Pain & hematoma.

Storage: Store at controlled room temperature.

Disposal: After use, this product may be a potential biohazard. Dispose in accordance with applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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