SCORPION® Portal Vein Access Set

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
Device Description

The device is for single-use, supplied sterile and sterilized by Ethylene Oxide.

The SCORPION® Portal Vein Access Set contains a 5F MPA catheter, a 14ga Stiffening Cannula, a 0.040" Nitinol Stylet with a 5F PEEK Catheter, and a 7F Cannula Sheath. The 14ga stiffening cannula with cannula sheath has a directional handle that indicates the direction of the curve. These components are used through an compatible introducer sheath to create a pathway through the liver parenchyma through which an endoprosthesis can be delivered. The SCORPION® Portal Vein Access Set is used to gain access to the hepatic vein and guide a sharp puncture tool (0.040" Stylet) through the parenchyma. The puncture tool (Stylet) is used to make a pathway from the hepatic vein to the portal vein, and then the pathway is dilated to provide access for a larger sheath. The shunt is inserted through the sheath and deployed through the pathway. Then, all of the SCORPION® Portal Vein Access Set components are removed.

The SCORPION® Portal Vein Access Set is compatible with a 10 F introducer sheath no longer than 45.5cm in overall length

Intended Purpose/ Intended Use

The SCORPION® Portal Vein Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

Intended Patient Population

Patients may be adults of any gender, any race or ethnicity, undergoing transjugular intrahepatic portosystemic shunt (TIPS) as treatment for the following:

- Variceal bleeding, secondary prevention, and acute bleeding refractory to medical and endoscopic treatments
- Refractory ascites
- Hepatorenal syndrome (types 1 and 2)
- Budd–Chiari syndrome
- Hepatic veno-occlusive disease
- Hepatic hydrothorax
- Portal hypertensive gastropathy
- Hepatopulmonary syndrome
- Portal vein thrombosis

Indication for Use

The SCORPION® Portal Vein Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

Intended Users

The product is intended for use by medical professionals and must only be used by qualified personnel who are familiar with the technique.

Duration/Lifetime

In use in procedures up to 4 hours

Clinical Benefits

The Portal Vein Access sets provide an indirect clinical benefit by facilitating creation of the parenchymal tract between the hepatic vein and the portal vein for delivery of a transjugular intrahepatic portosystemic shunt used to treat the complications of portal hypertension.

Risks and Side effects

- Intraperitoneal hemorrhage
- Puncture site hematoma
- Gall bladder puncture
- Cardiac arrhythmia
- Arteriovenous fistula
- Arterio-biliary fistula
- Death

Contraindications

None known

Warnings

- Contents are supplied sterile and are intended for single use only. Do not re-sterilize
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during the procedure.
- Do not use excessive force when advancing, retracting, or manipulating kit components and accessories.
- Use of power injection of contrast media to verify access to the portal vein may cause a loss of access.
- A patient may experience a negative reaction to angiography if the patient is allergic to contrast media or has compromised renal function.
- Risks and benefits should be carefully assessed in pregnant patients as radiation from fluoroscopic imaging may endanger the fetus.
- Do not bend the needle/catheter assembly more than 40° at any single point, as this might cause kinking of the needle/catheter.

Precautions

- Inspect package integrity before use.
- Inspect all set components prior to use
- For use by physicians trained in vascular diagnostic and interventional techniques only
- Do not overtighten hub connections.
- Use fluoroscopic guidance when advancing, retracting, or manipulating kit components and accessories in the vasculature.
- The patient should be monitored (EKG, oximetry, blood pressure, pulse) during the procedure.
- If the Portal Vein Access set is to be used for shunt delivery, the introducer sheath selected for portal vein access must also be compatible with the shunt, per the manufacturer recommendation.

Instructions for Use

1. Remove all components from packaging using sterile technique.
2. Flush components with saline (not supplied).
3. Prepare the 10F Introducer Sheath with the 10F Dilator.
4. Assemble 7F Cannula Sheath over the 14ga Stiffening Cannula.
5. Assemble the PEEK Catheter over the 0.040" Nitinol Stylet, then place the spacer clip in between the PEEK Catheter and 0.040" Nitinol Stylet.

Figure 1: Spacer Clip in Place

6. Obtain jugular access and advance a ≤ 0.035 " guidewire (not supplied) through the jugular vein, into the inferior vena cava.
7. Advance a compatible 10F Introducer Sheath and 10F Dilator (not supplied) over the guidewire into the IVC.
8. Remove the 10F Dilator.
9. Insert the 5F MPA Catheter over the guidewire into the hepatic vein. Advance the guidewire into the hepatic vein.
10. Advance the 10F Introducer Sheath to the hepatic vein.
11. Remove the 5F MPA Catheter, leaving the guidewire in place.
12. Introduce the 7F Cannula Sheath /14ga Stiffening Cannula assembly through the Introducer Sheath into the hepatic vein. Then remove guidewire. Be careful to maintain hepatic access with the 14ga Stiffening cannula and Introducer Sheath while removing the guidewire.
13. Insert the 0.040" Nitinol Stylet / PEEK Catheter assembly until the black line on the PEEK Catheter aligns with the end of the safety funnel on the 14ga Stiffening Cannula. The 0.040" Stylet tip should not be advanced past the PEEK Catheter.
14. Orient the 14ga Stiffening Cannula toward the portal vein.
15. Remove the spacer clip from the 0.040" Nitinol Stylet / PEEK Catheter assembly and retract the PEEK Catheter until the hub is locked with the handle on the 0.040" Nitinol Stylet to expose the stylet tip.
16. Wedging the 14ga Stiffening Cannula against the vein wall and orienting the 0.040" Nitinol Stylet towards the portal vein, advance the 0.040" Stylet forward through the parenchyma and toward the portal vein.
17. Remove the 0.040" Nitinol Stylet from the PEEK Catheter.
18. Confirm access to the portal vein.
19. Introduce a guidewire through the PEEK Catheter into the portal vein.
20. If necessary, dilate the parenchymal tract to accommodate the 10F Introducer Sheath.
21. Advance the 10F Introducer Sheath across the parenchymal tract and into the portal vein.
22. Remove the PEEK Catheter, 14ga Stiffening Cannula, and SCORPION® Cannula Sheath leaving the 10F Introducer Sheath in place.
23. Proceed with any indicated interventional procedures.

Storage

Store at controlled room temperature.

Disposal

After use, handle and dispose in accordance with hospitals policies and procedures concerning biohazard materials and waste

CAUTION

Federal Law (USA) restricts this device to sale by or on the order of a physician.

Note: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/ patient resides.

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SCORPION® is a registered trademark of Hatch Medical, L.L.C.