SCORPION®X Portal Vein Access Set

ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669

EMERGO EUROPE

Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Emergo Consulting (UK) Limited

c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

MedEnvoy Switzerland

Golthardstrasse 28
6302 Zug
Switzerland


<https://www.argonmedical.com/resources/product-information>

The symbols glossary is located electronically at
Речникът на символите е наличен в електронен вид на адрес
O glossário dos símbolos está localizado eletronicamente em
Slovníček symbolů v elektronické formě je umístěn na
Symbolforklaringen findes elektronisk på
Das Symbolglossar ist elektronisch verfügbar unter
Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση
El glosario de símbolos puede consultarse en línea en la dirección
Sõmbolite sõnastik on elektrooniliselt saadaval aadressil
Symbolisanasto on sähköisessä muodossa osoitteessa
Le glossaire des symboles est disponible sous forme électronique à l'adresse
Pojmovnik simbola nalazi se u elektroničkom obliku na
A szimbólumok összefoglalása elektronikusan itt található:
Il glossario dei simboli è disponibile in formato elettronico all'indirizzo
Simbolių žodynas elektronine forma pateiktas
Simbolu glosārijs ir atrodams elektroniskā veidā vietnē
Wersja elektroniczna słownika symboli znajduje się na stronie
O glossário de símbolos está localizado eletronicamente em
Glosarul simbolurilor este stocat electronic la
Glosár symbolov je v elektronickej podobe k dispozícii na stránkach
Slovar simbolov se nahaja v elektronski obliki na
Symbolordlistan finns elektroniskt på
www.argonmedical.com/symbols

EN - ENGLISH

Set Components

- A. 5F x 80cm MPA Catheter
- B. 13ga x 52cm Scorpion Stiffening Cannula
- C. 8F x 52cm Cannula Sheath
- D. 17ga x 73cm Scorpion Nitinol Needle
- E. 6.2F x 71cm Scorpion Needle PEEK Catheter
- F. Spacer Clip

Device Description

The device is for single-use, supplied sterile and sterilized by Ethylene Oxide.

The SCORPION®X Portal Vein Access Set contains a 5F MPA Catheter, a 13ga Stiffening Cannula, 17ga Nitinol Needle, a 6.2F PEEK Catheter, and an 8F Cannula Sheath. The 13ga stiffening cannula with cannula sheath has a directional handle that indicates the direction of the curve. The 17ga Nitinol Needle has a directional handle that indicates the direction of the curve.

These components are used through any compatible introducer sheath to create a pathway through the liver parenchyma through which an endoprosthesis can be delivered. The SCORPION®X Needle Portal Vein Access Set is used to gain access to the hepatic vein and guide a sharp puncture tool (17ga Needle) through the parenchyma. The puncture tool (Needle) is used to make a pathway from the hepatic vein to the portal vein, and then the pathway is dilated to provide access for a larger sheath. The shunt is inserted through the sheath and deployed through the pathway. Then, all of the SCORPION®X Portal Vein Access Set components are removed.

The SCORPION®X Portal Vein Access Set is compatible with a 10 F introducer sheath no longer than 45.5cm in overall length.

Intended Purpose/ Intended Use

The SCORPION®X Portal Vein Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

Intended Patient Population

Patients may be adults of any gender, any race or ethnicity, undergoing transjugular intrahepatic portosystemic shunt (TIPS) as treatment for the following:

- Variceal bleeding, secondary prevention, and acute bleeding refractory to medical and endoscopic treatments
- Refractory ascites
- Hepatorenal syndrome (types 1 and 2)
- Budd–Chiari syndrome
- Hepatic veno-occlusive disease
- Hepatic hydrothorax
- Portal hypertensive gastropathy
- Hepatopulmonary syndrome
- Portal vein thrombosis

Indication for Use

The SCORPION®X Portal Vein Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

Intended Users

The product is intended for use by medical professionals and must only be used by qualified personnel who are familiar with the technique.

Duration/Lifetime

In use in procedures up to 4 hours

Clinical Benefits

The Portal Vein Access sets provide an indirect clinical benefit by facilitating creation of the parenchymal tract between the hepatic vein and the portal vein for delivery of a transjugular intrahepatic portosystemic shunt used to treat the complications of portal hypertension.

Risks and Side effects

- Intraperitoneal hemorrhage
- Puncture site hematoma
- Gall bladder puncture
- Cardiac arrhythmia
- Arteriovenous fistula
- Arterio-biliary fistula
- Death

Contraindications

None known

Warnings

- Do not attempt to modify the 13ga Scorpion Stiffening Cannula bend angle.
- Contents are supplied sterile and are intended for single use only. Do not re-sterilize..
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Do not use if package is open or damaged and if the expiry data has been exceeded.
- Do not continue to use if any components are damaged during the procedure.
- Do not use excessive force when advancing, retracting, or manipulating kit components and accessories.
- Use of power injection of contrast media to verify access to the portal vein may cause a loss of access.
- A patient may experience a negative reaction to angiography if the patient is allergic to contrast media or has compromised renal function.
- Risks and benefits should be carefully assessed in pregnant patients as radiation from fluoroscopic imaging may endanger the fetus.
- Do not bend the needle/catheter assembly more than 40° at any single point, as this might cause kinking of the needle/catheter.

Precautions

- Inspect package integrity before use.
- Inspect all set components prior to use
- For use by physicians trained in vascular diagnostic and interventional techniques only
- Do not overtighten hub connections.
- Use fluoroscopic guidance when advancing, retracting, or manipulating kit components and accessories in the vasculature.
- The patient should be monitored (EKG, oximetry, blood pressure, pulse) during the procedure.
- If the Portal Vein Access set is to be used for shunt delivery, the introducer sheath selected for portal vein access must also be compatible with the shunt, per the manufacturer recommendation.

Instructions for Use

1. Remove all components from packaging using sterile technique.
2. Flush components with saline (not supplied).
3. Prepare the 10F Introducer Sheath with the 10F Dilator.

4. Assemble 8F Cannula Sheath over the 13ga Stiffening Cannula, then place the safety clip.
5. Assemble the Needle PEEK Catheter over the 17ga Nitinol Needle, then place the spacer clip in between the PEEK Catheter and the 17ga Nitinol Needle.

Figure 1: Spacer Clip in Place

6. Obtain jugular access and advance a ≤ 0.035 " guidewire (not supplied) through the jugular vein, into the inferior vena cava.
7. Advance a compatible 10F Introducer Sheath and 10F Dilator (not supplied) over the guidewire into the IVC.
8. Remove the 10F Dilator.
9. Insert the 5F MPA Catheter over the guidewire into the hepatic vein. Advance the guidewire into the hepatic vein.
10. Advance the 10F Introducer Sheath the hepatic vein and maintain access without occluding the vessel.
11. Remove the 5F MPA Catheter, leaving the guidewire in place.
12. Introduce the 8F Cannula Sheath/13ga Stiffening Cannula assembly through the Introducer Sheath into the hepatic vein. Be careful to maintain hepatic access with the 13ga Stiffening Cannula while removing the guidewire.
13. Insert the 17ga Needle / Needle PEEK Catheter assembly until the black line on the Needle PEEK Catheter aligns with the end of the safety funnel on the 13ga Stiffening Cannula. The 17ga Needle tip should not be advanced past the Needle PEEK Catheter.
14. Orient the 13ga Stiffening Cannula toward the portal vein.
15. Remove the spacer clip from the 17ga Nitinol Needle / Needle PEEK Catheter assembly and retract the Needle PEEK Catheter until the hub is locked with the handle on the 17ga Nitinol Needle to expose the stylet tip.
16. Wedging the 13ga Stiffening Cannula against the vein wall and orienting the 17ga Nitinol Needle towards the portal vein, advance the 17ga Nitinol Needle forward through the parenchyma and toward the portal vein.
17. Remove 17ga Nitinol Needle from the Needle PEEK Catheter and confirm access to the portal vein.
18. Introduce a guidewire through the Needle PEEK Catheter into the portal vein.
19. If necessary, dilate the parenchymal tract to accommodate the 10F Introducer Sheath.
20. Advance the 10F Introducer Sheath across the parenchymal tract until the tip of the 10F Introducer Sheath is in the portal vein. Be careful to avoid occluding the portal vein with the Introducer Sheath.
21. Remove the SCORPION®X Needle PEEK Catheter, 13ga SCORPION®X Stiffening Cannula, the 8F Cannula Sheath, leaving the 10F Introducer Sheath in place.
22. Proceed with any indicated diagnostic or interventional procedures through the 10F Introducer Sheath.

Storage

Store at controlled room temperature.

Disposal

After use, handle and dispose in accordance with hospitals policies and procedures concerning biohazard materials and waste.

CAUTION

Federal Law (USA) restricts this device to sale by or on the order of a physician.

Note: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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© 2021 Argon Medical Devices, Inc
SCORPION®X is a registered trademark of Hatch Medical, L.L.C.