CLEANER Vac™ Thrombectomy System

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669

EMERGO EUROPE

Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

MedEnvoy Switzerland
Golthardstrasse 28
6302 Zug
Switzerland

<https://www.argonmedical.com/products/>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът за символите в електронен вариант се намира на www.argonmedical.com/symbols
Croatian	Elektronički glosar simbola nalazi se na stranici www.argonmedical.com/symbols
Czech	Vysvětlivky symbolů jsou dostupné v elektronické podobě na adrese www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
Dutch	De verklarende woordenlijst van symbolen is elektronisch te vinden op www.argonmedical.com/symbols
Estonian	Sümbolite sõnastik asub elektrooniliselt aadressil www.argonmedical.com/symbols
Finnish	Symbolien sanasto sijaitsee sähköisesti osoitteessa www.argonmedical.com/symbols
French	Le glossaire des symboles est disponible sous forme électronique à l'adresse suivante : www.argonmedical.com/symbols
German	Ein elektronisches Symbolverzeichnis finden Sie unter www.argonmedical.com/symbols
Greek	Το γλωσσάρι συμβόλων βρίσκεται στην ηλεκτρονική διεύθυνση www.argonmedical.com/symbols
Hungarian	A szimbólumok szöszedete elektronikusan található a következő hivatkozás alatt: www.argonmedical.com/symbols

Italian	Il glossario dei simboli è disponibile alla pagina web www.argonmedical.com/symbols
Latvian	Elektroniskā simbolu vārdnīca ir pieejama www.argonmedical.com/symbols
Lithuanian	Simbolių žodynėlio elektroninė versija pateikta adresu www.argonmedical.com/symbols
Polish	Słownik symboli jest dostępny w formie elektronicznej na stronie www.argonmedical.com/symbols
Portuguese (PT)	O glossário de símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Romanian	Glosarul de simboluri se găsește în varianta electronică pe www.argonmedical.com/symbols
Slovak	Slovník symbolov je umiestnený elektronicky na adrese www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols
Spanish	El glosario de símbolos se puede encontrar en línea en www.argonmedical.com/symbols
Swedish	Ordlistan med symboler finner du under www.argonmedical.com/symbols

IFU9020M

Date of Issue: 2025/0920 Rev A

CLEANER Vac™ is a trademark of Argon Medical Devices, Inc. © 2023

English

CLEANER Vac™ Thrombectomy System

MEDICAL – GENERAL MEDICAL EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY
IN ACCORDANCE WITH IEC 60601-1:2005/A1:2012, CAN/CSA-C22.2
No. 60601-1:2014 and ANSI/AAMI ES60601-1:2005/A1:2012

Device Description

The CLEANER Vac™ Thrombectomy System is comprised of:

- CLEANER Vac™ Aspiration Canister
- CLEANER Vac™ 18F Aspiration Catheter and Handpiece

The CLEANER Vac™ Aspiration Canister is a fully disposable and sterile 400cc Aspiration Canister with an integrated pump and connected aspiration tubing. The Aspiration Canister contains a power switch with adjacent vacuum indicator (A), removeable lid, and Vacuum Release Button (B). The aspiration tubing contains a Tubing Clamp (C) and connector (D) to the Handpiece. The Aspiration Canister is shown in Figure 1.


Figure 1: Aspiration Canister

The CLEANER Vac™ 18F Aspiration Catheter and Handpiece consists of a Handpiece (Figure 2) with Aspiration Control Lever (E), connection to the Aspiration Catheter (F), and connection to the aspiration tubing (G) as well as a Flushing Adapter (H).


Figure 2: Handpiece

The CLEANER Vac™ 18F Aspiration Catheter and Handpiece also includes an 18F (OD) x 115cm Aspiration Catheter and Dilator. The Aspiration Catheter (Figure 3) contains a radiopaque marker band, Side Port with 3-way stopcock (I), Aspiration Port (J), and Dilator Port (K). The Side Port allows for infusion of fluids. The Aspiration Port connects to the Handpiece, and the Dilator Port allows for guidewire and Dilator insertion. The Dilator includes a rotating collar used to secure the dilator to the catheter. There is a hemostasis valve at the Aspiration Port and Dilator Port to create a seal and prevent leakage of fluids.


Figure 3: Aspiration Catheter with Dilator

The emissions characteristics of this equipment make it suitable for use in hospitals (CISPR 11 Class A), except near active HF surgical equipment and the RF shielded room of an MRI, where EM disturbances are high.

The Essential Performance for the CLEANER Vac™ system is the ability of the physician to control aspiration using the Aspiration Control Lever on the Handpiece.

The use of the device is supported from 18 - 24 degrees Celsius, 30% to 60% relative humidity, and an altitude up to 2000m.

The device does not contain natural rubber latex.

Intended Purpose/ Intended Use

The CLEANER Vac™ Thrombectomy System is intended for the removal of fresh, soft thrombi and emboli from the vessels of the peripheral venous vasculature.

Intended Patient Population

Adults presenting with thrombi or emboli, within the peripheral venous vasculature, requiring intervention via percutaneous mechanical thrombectomy.

Indication for Use

The CLEANER Vac™ Thrombectomy System is indicated for the removal of fresh, soft thrombi and emboli from the vessels of the peripheral venous vasculature, and for the infusion of physician-specified fluids, including thrombolytics. The CLEANER Vac™ Thrombectomy System is not intended for use in the pulmonary vasculature for treating of pulmonary embolism.

Intended Users

Intended for use by licensed physicians or clinicians trained in vascular diagnostic and interventional techniques who are familiar with the procedure.

Duration/Lifetime

Transient (less than 60 minutes)

Clinical Benefits

The Cleaner Vac™ Thrombectomy System provides a direct clinical benefit with the immediate reduction of clot burden facilitating restoration of blood flow to the target vessels.

Risks and side effects

List of Complications for Peripheral Venous Thrombectomy

Thrombosis & Embolization Risks

- Acute occlusion
- Distal embolization
- Embolism (including paradoxical embolism in patients with PFO)
- Inability to completely remove thrombus
- Thrombotic events (e.g., re-thrombosis)

Vascular & Access Site Complications

- Access site hematoma
- Access site hemorrhage
- Excessive blood loss
- Vessel thrombosis
- Vessel dissection
- Vessel perforation
- Vessel stenosis (Venous)
- Arteriovenous fistula (if unintended arterial access occurs during cannulation)
- Pseudoaneurysm (at access site)
- Aneurysm

Hemodynamic & Blood-Related Risks

- Hemolysis
- Hemoglobinuria
- Heparin-induced thrombocytopenia (HIT)
- Hypotension or hemodynamic instability

Infection & Inflammatory Risks

- Infection (including cellulitis or sepsis)
- Inflammatory response
- Drug reaction/allergic reaction (contrast, thrombolytics, anticoagulants)
- Adverse reaction to device materials

Organ & Systemic Risks

- Anemia (secondary to blood loss)
- Hypoxemia (due to embolization or large clot burden)
- Renal failure (secondary to contrast-induced nephropathy or embolic events)
- Organ impairment (if embolization affects organ perfusion)
- Arrhythmias
- Dyspnea
- Ischemia
- Respiratory arrest
- Cardio-respiratory arrest
- Respiratory failure

Procedural & Device-Related Risks

- Air embolism
- Foreign body embolism
- Vascular spasm
- Spontaneous movement or migration of catheter during procedure
- Intimal disruption
- Fistulation

Neurological Risks (Paradoxical Embolism-Related)- Patients with PFO

- Stroke / Transient Ischemic Attack (TIA) (in patients with right-to-left shunting)
- Neurological deficit (if embolism affects cerebral circulation)
- Peripheral nerve damage

Non-Specific Procedural Symptoms

- Pain, tenderness, or discomfort
- Nausea / vomiting
- Fever
- Death

Contraindications

Not intended for use in the coronaries and neurovasculature.

Warnings

- Reuse or reprocessing has not been evaluated and may lead to device failure and subsequent patient illness, infection, or other injury.
- Inspect package integrity before use. Do not use if the package is open or damaged and if the expiry date has been exceeded.
- Inspect component integrity before use. Do not use any damaged component.
- Do not continue to use any component damaged during the procedure.
- Advance and retract the device under fluoroscopic guidance.
- Gently advance the Dilator and Aspiration Catheter through thrombus to avoid risk of embolization.

- Avoid over injection of fluids to minimize risk of embolization.
 - Withdrawing the Aspiration Catheter with caution is recommended when sharp radii are encountered. A withdrawal rate of 1-2 cm/second is recommended.
 - If the device needs to be reinserted, to avoid reintroduction of thrombus, flush the device with heparinized saline before proceeding with device reinsertion.
 - Do not advance or retract if resistance is met without first determining the cause of resistance under fluoroscopy and taking any necessary action. Excessive force against resistance may result in damage to the device or vasculature.
 - Prior to introduction, and anytime the Aspiration Catheter is removed from the vascular system, the Aspiration Catheter should be flushed.
 - Do not initiate aspiration unless proper device positioning is confirmed within the vasculature.
 - Do not use the aspiration system with a pump other than the CLEANER Vac™ Aspiration Canister.
 - When performing aspiration, ensure that the Aspiration Control Lever is only open for the minimum time needed to remove thrombus.
 - Do not re-infuse blood or fluid from the canister back into the patient.
 - Do not power inject through the Side Port.
 - Before injection of fluids, aspirate through the side port to remove any residual air within the catheter.
 - Minor air ingress through the Aspiration Catheter's hemostasis valve during aspiration of the Side Port tubing is possible; do not reintroduce air into the catheter.
 - Follow all contraindications, warnings, cautions, precautions, and instructions for use for all infusates, including contrast medium, as specified by their manufacturer.
 - Excessive blood loss may occur if not monitored throughout the procedure.
 - No modification of this equipment is allowed.
 - The flush adapter should not be used to inject fluids into the patient.
3. Remove all sterile components from packaging using sterile technique and inspect components for damage. Replace any damaged components; do not use until all damaged components have been replaced. Ensure that the battery pull tab on the Aspiration Canister is removed and the pin on the Handpiece is removed.
 4. Position the Aspiration Canister on a flat, stable surface capable of bearing its weight. The canister should be in a position that is viewable by the physician or staff and the aspiration tubing can reach the sterile field.
 5. Confirm that the canister lid is fully closed, and close the Tubing Clamp Then, turn the Aspiration Canister on by switching the power switch up on the base of the Aspiration Canister. Confirm the Aspiration Canister is on through audible and visual feedback (i.e., the bottom bar of the vacuum indicator should illuminate in blue).
Note: If the Aspiration Canister doesn't power on, remove the battery door, confirm proper contact, and replace the battery door.
 6. Let the Aspiration Canister sit until the vacuum indicator illuminates fully (i.e., all bars are illuminated in blue), at which point the Aspiration Canister has reached optimal conditions for aspiration. The vacuum indicator shall illuminate fully when the Aspiration Canister has reached greater than 25 inHg. If the vacuum indicator does not fully illuminate, stop and replace the canister.
 7. Flush the Dilator through the hub and flush the Aspiration Catheter through the Side Port.
 8. Insert the Dilator into the Dilator Port of the Aspiration Catheter.
 9. Flush the Handpiece prior to use. To flush the Handpiece, connect the Flushing Adapter to the back of the Handpiece. Then, connect a syringe with heparinized saline to the Flushing Adapter and infuse through the Handpiece while pressing down on the Aspiration Control Lever.
 10. Connect the aspiration tubing from the Aspiration Canister to the connector on the back of the Handpiece. Ensure that the tubing is connected to the Handpiece and open the Tubing Clamp.

Precautions

- The Aspiration Canister should be turned off when not in use.
- If needed, the device should be power cycled to reset the Aspiration Canister.
- If aspiration cannot be turned off, close the Tubing Clamp located on the aspiration tubing to stop aspiration.
- If the Aspiration Canister falls off the table, the user shall discontinue use of that device because it may be damaged and non-sterile.
- Before the canister reaches first full capacity, an audible tone will sound, at which point the canister should be turned off before the volume reaches 400cc. The audible tone that indicates the canister is full will not sound on subsequent cycles after draining canister.
- If use of a defibrillator is needed, the CLEANER Vac™ device needs to be removed from the patient before its use
- Use of this device adjacent to other equipment, including equipment that are known sources of electromagnetic disturbances such as diathermy, electrocautery, RFID equipment and Magnetic Resonance (MR) scanners, should be avoided because it could result in improper operation. If such use is necessary, the device and other equipment should be observed to verify that they are operating normally. Note that some of these RF emitters (e.g., RFID) in the intended environment of use might be concealed and the device can potentially be exposed to fields from these RF emitters without the user's awareness.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of the Cleaner™ system. Otherwise, degradation of the performance of this device could result.
- Use of this device near Magnetic Resonance (MR) equipment and scanners is not supported since the device is not MR safe.

Pre-Procedure Preparation

1. Select an appropriately sized introducer sheath to accommodate the 18F (OD) CLEANER Vac™ Aspiration Catheter.
2. Select a guidewire to accommodate passage of the 115cm CLEANER Vac™ Aspiration Catheter and CLEANER Vac™ Dilator. The CLEANER Vac™ Dilator is compatible with a guidewire ≤ 0.038".

Directions for Use

1. Under fluoroscopic guidance, gently advance the Dilator and Aspiration Catheter over the guidewire to the target site.
2. Remove the Dilator.
3. Connect the Aspiration Port of the Aspiration Catheter to the connector on the front of the Handpiece. Ensure the Aspiration Catheter is connected securely to the Handpiece before proceeding with aspiration.
4. Confirm device positioning under fluoroscopy.
5. Perform aspiration. Press and hold the Aspiration Control Lever to activate aspiration. Note: To deactivate aspiration, release the Aspiration Control Lever.
6. Monitor blood loss through the canister while aspirating.
Note: Before the canister reaches first full capacity, an audible tone will sound, at which point the canister should be turned off if possible.
Note: Once the audible tone sounds, the Aspiration Canister will continue to evacuate for 25 more seconds, and then automatically shut off to allow the canister to be emptied.
Note: An overflow valve will disable the generation of vacuum within the canister if the canister reaches full capacity.
Note: To empty the canister, turn the Aspiration Canister off by pressing down the power switch, then press the Vacuum Release Button on the canister lid to remove the vacuum within the canister. Remove the canister lid and empty contents, per facility protocol.
Note: If additional passes are needed, empty the canister contents and manually reset the float valve by pushing it down. Then, place the lid back on the canister, and confirm the lid is fully closed before proceeding. Turn the Aspiration Canister on to re-establish vacuum in the canister.
Note: The audible tone that indicates the canister is full will sound on first fill only. Any subsequent cycles must be visually monitored by the user to avoid excessive blood loss.
7. If the system clogs, deactivate aspiration, and remove the device from the patient. Then, flush the Aspiration Catheter and Handpiece to remove any thrombus.
Note: To flush the Handpiece, disconnect the Handpiece from the Aspiration Catheter, close the Tubing Clamp on the aspiration tubing, then disconnect the aspiration tubing from the Handpiece and connect the Flushing Adapter. Then, with the Aspiration Control Lever depressed, connect a

syringe with heparinized saline to the Flushing Adapter and flush the device to remove any thrombus.

Note: If clogging still exists, replace the clogged device.

8. Aspiration Catheter may be used for injection of heparinized saline, contrast, or thrombolytics through the Side Port during the procedure, per the physician's discretion.

Table 1: Dead Space Volume

CLEANER Vac™ Aspiration Catheter 18F x 115 cm	33.2 mL
--	---------

9. If additional passes of the device are needed, repeat the procedure until acceptable thrombus removal is achieved. When thrombus removal is complete, turn off the Aspiration Canister and remove the device.

Storage and Handling Conditions

- The System is intended to be stored in a controlled room temperature.
- Do not expose the system to organic solvents, ionizing radiation, or ultraviolet light.

Disposal

- After use, handle and dispose in accordance with hospitals policies and procedures concerning biohazard materials and waste
- To empty the canister, press the Vacuum Release Button on the canister lid to remove the vacuum within the canister. Then, remove the canister lid and empty contents, per facility protocol.
- Device batteries in the Aspiration Canister are lithium ion and may be removed for separate disposal if required.
- Dispose of the system in accordance with the Waste Electrical and Electronic Equipment Directive 2012 (WEEE) and according to the standard institutional procedures for medical waste including single-use, blood contacting devices.

Note: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

Disclaimer of Warranty and Limitation of Remedy

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Manufacturer or its Distributors product(s) described in this publication. Under no circumstances shall the Manufacturer or its Distributor be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind the Manufacturer or its Distributor to any representation or warranty except as specifically set forth herein. Descriptions or specifications in the manufacturer and distributor's printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Manufacturer and Distributor will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

Electromagnetic Compatibility (EMC) Safety

EMC Test	Basic EMC Standard	Test Level
RF Conducted and Radiated Emissions	CISPR 11	Class A – Group 1
Electrostatic Discharge	IEC 61000-4-2	±8kV Contact ±15kV Air
Radiated RF EM Field	IEC 61000-4-3	3V/m – 80 MHz to 2.7 GHz, and frequencies from IEC 60601-1-2:2014 Table 9
Power Frequency Magnetic Fields	IEC 61000-4-8	30 A/m

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Varumärket **CLEANER Vac™** tillhör Argon Medical Devices, Inc.© 2023