




**Argon**

Medical Devices

# Product Catalog



# General Information




# Table of Contents

<b>GENERAL INFORMATION</b>	<b>A</b>	Cournand Style	F-2
Table of Contents	A-1	Safety Arterial	F-2
<b>THROMBECTOMY PRODUCTS</b>	<b>B</b>	<b>GUIDE WIRES</b>	<b>G</b>
Cleaner Thrombectomy System	B-1	Straight	G-1
		J Tip	G-2
		Specialty	G-3
<b>HEART BIOPSY FORCEPS</b>	<b>C</b>	<b>VESSEL DILATORS</b>	<b>H</b>
Jawz EMB Forceps	C-1	Vessel Dilators	H-1
Biopsy Accessories	C-2		
<b>PERCUTANEOUS INTRODUCERS</b>	<b>D</b>	<b>MANIFOLDS</b>	<b>I</b>
AXcess Introducer	D-1	Manifolds (High Pressure)	I-1
OneStic Introducers	D-2	Stopcock (High Pressure)	I-1
MicroStic Introducer	D-3		
Valve less Introducers	D-4	<b>SYRINGES</b>	<b>J</b>
Detachable Valve Introducers	D-4, D-6	Coronary Control	J-1
1- Piece Introducers	D-5, D-6	Inflation Device	J-1
Tuohy-Borst Introducers	D-5, D-6	Medication	J-2
Multi-Purpose Curve Introducers	D-7	Labeling Set	J-2
Specialty Trays	D-8		
Accessories	D-8	<b>FLUID MANAGEMENT TUBING</b>	<b>K</b>
<b>DRAPES &amp; EQUIPMENT COVERS</b>	<b>E</b>	High Pressure Injection Lines	K-1
Angiography	E-1	Pressure Monitoring Lines	K-1
Pacing	E-2	Contrast Management Systems	K-2
Specialty	E-2	IV Spikes	K-2
Equipment Covers	E-3	Waste Management Systems	K-3
<b>ARTERIAL NEEDLES</b>	<b>F</b>	<b>PRESSURE TRANSDUCERS</b>	<b>L</b>
Thin Wall	F-1	Transducers	L-1
Percutaneous Entry	F-2	Single Pressure Kits	L-1, L-2
		Multiple Pressure Kits	L-2
		KeepSafe Closed Sampling	L-3

# Table of Contents


<b>MONITORING ACCESSORIES</b>	<b>M</b>
Pressure Monitoring Lines	M-1
Stopcocks (Low Pressure)	M-2
Extension Sets	M-2
Flush Devices	M-2
Cardiac Output Bridges	M-3
Pole Mount Manifold	M-3
Luers, Adapters and Plugs	M-4
Brackets and Holders	M-4
<b>MONITORING CATHETER KITS/TRAYS</b>	<b>N</b>
Arterial Catheter Kits	N-1
Pulmonary Artery Catheters	N-2
Thermodilution Catheters	N-2
Temporary Pacing Trays	N-3
Central Venous Catheter Kits/Trays	N-4
<b>INDEX</b>	<b>O</b>
Ordering Information/Product Return Policy	O-1
Catalog Alpha-Numeric Reference	O-1

# Thrombectomy Products




# Mechanical Thrombectomy




Product Code	Description	Catheter Length	Diameter	Sheath Included/Fr.	Units per Box
<b>Thrombectomy Systems</b>					
	<b>700009</b> Cleaner 6F Thrombectomy Device	65cm	9mm		1
	<b>700109</b> Cleaner 6F Thrombectomy Device	135cm	9mm		1
	<b>700115</b> Cleaner 6F Thrombectomy Device	135cm	15mm		1
	<b>726009</b> Cleaner 6F Thrombectomy System	65cm	9mm	2/6F	1
	<b>727009</b> Cleaner 6F Thrombectomy System	65cm	9mm	2/7F	1
	<b>716109</b> Cleaner 6F Thrombectomy Kit	135cm	9mm	1/6F	1
	<b>717109</b> Cleaner 6F Thrombectomy Kit	135cm	9mm	1/7F	1
	<b>716115</b> Cleaner 6F Thrombectomy Kit	135cm	15mm	1/6F	1
	<b>717115</b> Cleaner 6F Thrombectomy Kit	135cm	15mm	1/7F	1

# Endomyocardial Biopsy Forceps




# Endomyocardial (EMB) Biopsy Forceps




Product Code	French Size	Length (cm)	Shaft Configuration	FEP Coating	Units Per Box
<b>2.4mm Forceps</b>					
190020	7.5F	50cm	Straight		5
190030	7.5F	50cm	Pre-Curved	Yes	5
190031	7.5F	50cm	Maxi-Curved	Yes	5
190040	7.5F	105cm	Straight		5
<b>2.2mm Forceps</b>					
190080	7F	50cm	Pre-Curved	Yes	5
190081	7F	50cm	Maxi-Curved	Yes	5
190085	7F	105cm	Straight		5
190086	7F	105cm	Straight	Yes	5
<b>2.2mm Forceps (Non-Formable)</b>					
190065	6F	50cm	Straight		5
190075	6F	105cm	Straight		5
<b>1.8mm Forceps</b>					
190060	5.5F	50cm	Straight		5
190051	5.5F	50cm	Maxi-Curved	Yes	5
190070	5.5F	105cm	Straight		5
<b>1.5mm Forceps</b>					
190090	5F	50cm	Straight		5
190095	5F	105cm	Straight		5

# Biopsy Accessories

Product Code	French Size	Length (cm)	Guide Wire Diameter/Length	Hub Color	Units Per Box
<b>Access Introducers with Double ended wire</b>					
496005	5F	11cm	0.038" / 45cm	Gray	10
496006	6F	11cm	0.038" / 45cm	Green	10
496007	7F	11cm	0.038" / 45cm	Orange	10
496008	8F	11cm	0.038" / 45cm	Blue	10
<b>OneStic Introducers with Stainless Steel Mandrel Wire and 21G Needles</b>					
497830	5F	11cm	0.018" / 45cm	Gray	10
497840	6F	11cm	0.018" / 45cm	Green	10
497861	7F	11cm	0.018" / 45cm	Orange	10
497862	8F	11cm	0.018" / 45cm	Blue	10
<b>OneStic Introducers with Nitinol Mandrel Wire and 21G Needles</b>					
497831	5F	11cm	0.018" / 45cm	Gray	10
497841	6F	11cm	0.018" / 45cm	Green	10
<b>Introducer with Multi-Purpose Curve Tip</b>					
191545	5F	45cm		Gray	10
191598	5F	98cm		Gray	10
191645	6F	45cm		Green	10
191698	6F	98cm		Green	10
191745	7F	45cm		Orange	10
191798	7F	98cm		Orange	10
191845	8F	45cm		Blue	10
191898	8F	98cm		Blue	10


Product Code	Description	Drape Material	Length	Absorbent Pad	Fenestration	Units Per Box
 396303	Heart Biopsy Drape With Facial Port	SMS	77" x 100 (195 x 254cm)	32"x34" (81 x 86cm)	3.5" (9cm)	20

# Percutaneous Introducer Sheaths



# AXcess™ Percutaneous Introducer Kits

**AXcess**  
Introducer



Product Code	French Size	Length (cm)	Guide Wire Diameter/Length	Hub Color	Units per Box
<b>Pediatric Length with wire</b>					
496064	4F	5.75cm	0.021" x 40cm	Red	10
496065	5F	5.75cm	0.025" x 40cm	Gray	10
496066	6F	5.75cm	0.032" x 40cm	Green	10
496067	7F	5.75cm	0.032" x 40cm	Orange	10
<b>Standard Length with Double ended wire</b>					
496004	4F	11cm	0.035" x 45cm	Red	10
496005	5F	11cm	0.038" x 45cm	Gray	10
496006	6F	11cm	0.038" x 45cm	Green	10
496007	7F	11cm	0.038" x 45cm	Orange	10
496008	8F	11cm	0.038" x 45cm	Blue	10
496009	9F	11cm	0.038" x 45cm	Black	10
<b>Standard Length - without Wire</b>					
496035	5F	11cm	Accepts 0.035" - 0.038"	Gray	10
496036	6F	11cm	Accepts 0.035" - 0.038"	Green	10
496037	7F	11cm	Accepts 0.035" - 0.038"	Orange	10
496038	8F	11cm	Accepts 0.035" - 0.038"	Blue	10
<b>Long Sheath - without Wire</b>					
490005	5F	23cm	Accepts 0.035" - 0.038"	Gray	10
490006	6F	23cm	Accepts 0.035" - 0.038"	Green	10
490027	7F	23cm	Accepts 0.035" - 0.038"	Orange	10
490028*	8F	23cm	Accepts 0.035" - 0.038"	Blue	10
490029*	9F	23cm	Accepts 0.035" - 0.038"	Black	10

\* Dual Dilator Systems

# OneStic™ Percutaneous Introducer Kits




Product Code	French Size	Length (cm)	Guide Wire Diameter/Length	Guide Wire Core Tip	Needle Gauge/Lengths	Units per Box
<b>Stainless Steel Mandrel Wires</b>						
497850	6F	5.75cm	0.018" x 45cm	SS/SS	21G x 5.1cm	10
497854	6F	5.75cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497863	7F	5.75cm	0.018" x 45cm	SS/SS	21G x 5.1cm	10
497852	8F	5.75cm	0.018" x 45cm	SS/SS	21G x 5.1cm	10
497820	4F	11cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497845	5F	11cm	0.018" x 45cm	SS/SS	21G x 2.5cm	10
497830	5F	11cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497840	6F	11cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497861	7F	11cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497862	8F	11cm	0.018" x 45cm	SS/SS	21G x 7cm	10
<b>Nitinol Wires with Platinum Tip</b>						
497868	4F	5.75cm	0.018" x 45cm	NT/PT	21G x 2.5cm	10
497865	5F	5.75cm	0.018" x 45cm	NT/PT	21G x 2.5cm	10
497864	6F	5.75cm	0.018" x 45cm	NT/PT	21x 2.5cm	10
497853	6F	5.75cm	0.018" x 45cm	NT/PT	21G x 5.1cm	10
497849	6F	5.75cm	0.018" x 45cm	NT/PT		10
497867	7F	5.75cm	0.018" x 45cm	NT/PT	21G x 2.5cm	10
497821	4F	11cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497831	5F	11cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497841	6F	11cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497844	6F	11cm	0.018" x 45cm	NT/PT	21G x 5.1cm	10
<b>Echogenic Needles</b>						
497860	7F	5.75cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497858	4F	11cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497857	4F	11cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497859	5F	11cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497856	6F	11cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497855	6F	11cm	0.018" x 45cm	NT/PT	21G x 7cm	10

# MicroStic™ Percutaneous Introducer Kits



Product Code	French Size	Length (cm)	Guide Wire Diameter/Length	Guide Wire Core Tip	Needle Gauge/Lengths	Units per Box
<b>Stainless Steel Mandrel Wire</b>						
497801	4F	10cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497802	5F	10cm	0.018" x 45cm	SS/SS	21G x 7cm	10
<b>Stainless Steel Mandrel Wire and Echogenic Needle</b>						
497816	4F	10cm	0.018" x 45cm	SS/SS	21G x 7cm	10
497817	5F	10cm	0.018" x 45cm	SS/SS	21G x 7cm	10
<b>Nitinol Wire with Platinum Tip</b>						
497811	4F	10cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497812	5F	10cm	0.018" x 45cm	NT/PT	21G x 7cm	10
<b>Nitinol Wire with Platinum Tip and Echogenic Needle</b>						
497813	4F	10cm	0.018" x 45cm	NT/PT	21G x 7cm	10
497814	5F	10cm	0.018" x 45cm	NT/PT	21G x 7cm	10
<b>Nitinol Wire with Stainless Steel Tip and Echogenic Needle</b>						
497818	4F	10cm	0.018" x 45cm	NT/SS	21G x 7cm	10
497819	5F	10cm	0.018" x 45cm	NT/SS	21G x 7cm	10
Product Code	Diameter	Length (cm)	Guide Wire Diameter/Length	Guide Wire Core Tip	Units per Box	
<b>MicroStic Plus Stiffen Introducer with Stainless Steel Mandrel Wire</b>						
497880	4F	10cm	0.018" x 45cm	SS/SS	10	
497881	5F	10cm	0.018" x 45cm	SS/SS	10	
<b>MicroStic Plus Stiffen Introducer with Nitinol Wire with Stainless Tip and Echogenic Needle</b>						
497884	4F	10cm	0.018" x 45cm	NT/SS	10	
497885	5F	10cm	0.018" x 45cm	NTSS	10	
<b>MicroStic Plus Stiffen Introducer with Nitinol Wire with Platinum Tip</b>						
497882	4F	10cm	0.018" x 45cm	NT/PT	10	
497883	5F	10cm	0.018" x 45cm	NT/PT	10	

# Introducer Kits for Venous Applications



Product Code	French Size (cm)	Length	Wire Diameter / Needle	Hub Color	Units per Box
--------------	------------------	--------	------------------------	-----------	---------------

## PCI Kit/3 - Introducer without Hemostasis Valve

496140	4F	11cm	0.021" / 21Gx1.0"	Red	10
496141	5F	11cm	0.028" / 19Gx1.5"	Gray	10
496142	6F	11cm	0.035" / 18Gx2.75"	Green	10

### Kit Contents:

- |  |                               |
|--|-------------------------------|
| 1- Disposable Towel                            | 1- 22g x 1.5" (3.8cm) Needle  |
| 3- 4" x 4" Gauze Sponges                       | 1- 5cc Syringe                |
| 1- Vessel Dilator                              | 1- Radiopaque Sheath          |
| 1- #11 Mini Scalpel                            | 1- Arterial Needle, Clear Hub |
| 1- Guide Wire (Double ended 50cm, 3mmJ and ST) |                               |



Product Code	French Size (cm)	Length	Guide Wire Diameter/length	Hub Color	Units per Box
--------------	------------------	--------	----------------------------	-----------	---------------

## PCI Kit/2 - Introducer with Detachable Hemostasis Valve

496126	4F	5.75cm	Accepts 0.021" wire	Red	10
496127	5F	5.75cm	Accepts 0.028" wire	Gray	10
496128	6F	5.75cm	Accepts 0.038" wire	Green	10
496124	8F	11cm	Accepts 0.038" wire	Blue	10
496125	9F	11cm	Accepts 0.038" wire	Black	10

## PCI Kit/2 - Introducer with Detachable Hemostasis Valve, Feathered Tip Dilator and Wire

496308	4F	11cm	0.021" x 80cm*	Red	10
496309	5F	11cm	0.021" x 80cm*	Gray	10
496310	6F	11cm	0.021" x 80cm*	Green	10
496311	7F	11cm	0.021" x 80cm*	Orange	10
496312	8F	11cm	0.021" x 80cm*	Blue	10

\*Guide Wire (Double ended, 3mmJ and ST)

# Introducer Kits for Venous Applications

Product Code	French Size (cm)	Length	Guide Wire Diameter/length	Hub Color	Units per Box
--------------	------------------	--------	----------------------------	-----------	---------------

## PCI Kit/5 - Introducer with Detachable Hemostasis Valve



<b>496181</b>	5F	11cm	0.028" x 50cm	Gray	10
<b>496182</b>	6F	11cm	0.035" x 45cm	Green	10
<b>496183</b>	7F	11cm	0.035" x 45cm	Orange	10

### Kit Contents:

- 1- Disposable Towel
- 1- Vessel Dilator
- 1- Guide Wire (Double ended, 3mmJ and ST)
- 3- 4" x 4" Gauze Sponges
- 1- Sheath with Detachable Hemostasis Valve

Product Code	French Size (cm)	Length	Guide Wire Diameter/length	Hub Color	Units per Box
--------------	------------------	--------	----------------------------	-----------	---------------

## One Piece Introducer with Integral Hemostasis Valve



<b>496150</b>	4F	11cm	0.025" x 45cm	Red	10
<b>496151</b>	5F	11cm	0.035" x 45cm	Gray	10
<b>496152</b>	6F	11cm	0.035" x 45cm	Green	10
<b>496153</b>	7F	11cm	0.035" x 45cm	Orange	10
<b>496154</b>	8F	11cm	0.035" x 45cm	Blue	10
<b>496155</b>	8.5F	11cm	0.035" x 45cm	White	10
<b>496156</b>	9F	11cm	0.035" x 45cm	Black	10

## Tuohy-Borst Introducer with Integral Hemostasis Valve






<b>496132</b>	6F	11cm	0.035" x 45cm	Green	10
<b>496133</b>	7F	11cm	0.035" x 45cm	Orange	10
<b>496134</b>	8F	11cm	0.035" x 45cm	Blue	10
<b>496135</b>	8.5F	11cm	0.035" x 45cm	White	10
<b>496136</b>	9F	11cm	0.035" x 45cm	Black	10

### Kit Contents:


- 1- Radiopaque Sheath
- 1- Guide Wire (Double ended, 3mmJ and ST)
- 1- Vessel Dilator

# Introducer Trays



	Product Code	French Size (cm)	Length	Guide Wire Diameter/length	Hub Color	Units per Box
<b>Introducer Tray with Detachable Hemostasis Valve</b>						
	497422	6F	11cm	0.035" x 45cm	Green	10
	497423	7F	11cm	0.035" x 45cm	Orange	10
	497424	8F	11cm	0.035" x 45cm	Blue	10
	497425	8.5F	11cm	0.035" x 45cm	White	10
	497426	9F	11cm	0.035" x 45cm	Black	10
<b>One Piece Introducer Tray with Integral Hemostasis Valve</b>						
	497432	6F	11cm	0.035" x 45cm	Green	10
	497433	7F	11cm	0.035" x 45cm	Orange	10
	497434	8F	11cm	0.035" x 45cm	Blue	10
	497435	8.5F	11cm	0.035" x 45cm	White	10
	497436	9F	11cm	0.035" x 45cm	Black	10
<b>Tuohy-Borst Introducer Tray with Integral Hemostasis Valve</b>						
	497452	6F	11cm	0.035" x 45cm	Green	10
	497453	7F	11cm	0.035" x 45cm	Orange	10
	497454	8F	11cm	0.035" x 45cm	Blue	10
	497455	8.5F	11cm	0.035" x 45cm	White	10
	497456	9F	11cm	0.035" x 45cm	Black	10
<b>Tray Contents:</b>						
3- Chloraprep™ Prep Swab		1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub				
1- Disposable Towel		1- 5cc Ampule Lidocaine HCL 1%				
3- 4" x 4" Gauze Sponges		1- Guide Wire (Double ended, 3mmJ and ST)				
1- Fenestrated Drape		2- 5cc Syringe L/S				
1- Introducer Sheath		1- Vessel Dilator				
1- 3cc Syringe L/L		1- 25g x 5/8" (1.6cm) Needle				
1- 22g x 1.5" (3.8cm) Needle		1- 000 Silk Suture				
1- #11 Mini Scalpel		1- Povidone/Iodine Ointment				
1- Septishield Catheter Sleeve		1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter				



# Introducer Kits with Multi-Purpose Curve




Product Code	French Size (cm)	Length	Guide Wire Diameter/length	Hub Color	Units per Box
<b>191545</b>	5F	45cm	Accepts 0.035" - 0.038"	Gray	10
<b>191598</b>	5F	98cm	Accepts 0.035" - 0.038"	Gray	10
<b>191645</b>	6F	45cm	Accepts 0.035" - 0.038"	Green	10
<b>191698</b>	6F	98cm	Accepts 0.035" - 0.038"	Green	10
<b>191745</b>	7F	45cm	Accepts 0.035" - 0.038"	Orange	10
<b>191798</b>	7F	98cm	Accepts 0.035" - 0.038"	Orange	10
<b>191845</b>	8F	45cm	Accepts 0.035" - 0.038"	Blue	10
<b>191898</b>	8F	98cm	Accepts 0.035" - 0.038"	Blue	10

# Specialty Trays and Accessories











Product Code	Description	Units per Box		
<b>500200</b>	<b>Universal Percutaneous Introducer Tray</b>	10		
<b>Tray Contents:</b>				
3- Chloraprep™ Prep Swab	1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub			
2- Disposable Towel	2- 5cc Ampule Lidocaine HCL 1%			
6- 4" x 4" Gauze Sponges	1- .035" x 65cm FC, SS, 3mmJ & ST Guide Wire			
1- Fenestrated Drape w/Ripline	2- 5cc Syringe L/L			
1- Injection Site	1- 8F Vessel Dilator			
1- 25g x 5/8" (1.6cm) Needle	1- 22g x 1.5" (3.8cm) Needle			
1- 000 Silk Suture	1- #11 Mini Scalpel			
1- Povidone/Iodine Ointment	1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter			
Product Code	Description	Units per Box		
<b>195110</b>	<b>Multi-Purpose Biopsy Kit (without Specimen Tubes or Slides)</b>	10		
<b>195100</b>	<b>Multi-Purpose Biopsy Tray</b>	10		
<b>Tray Contents:</b>				
3- Chloraprep™ Prep Swab	1- 20g x 1.5" (3.8cm) Needle			
1- Disposable Towel	1- 5cc Ampule Lidocaine HCL 1%			
4- 4" x 4" Gauze Sponges	1- 3cc Syringe L/L			
1- Fenestrated Drape w/Ripline	1- 5cc Syringe L/S			
1- Plastic Ruler 6" (15cm)	1- 20cc Syringe L/L			
1- 25g x 5/8" (1.6cm) Needle	1- 22g x 1.5" (3.8cm) Needle			
1- Adhesive Bandage	1- #11 Long Handle Scalpel			
6- Microscope Slides	1- Medicine Cup 30cc			
2- Specimen Tube w/cap				
Product Code	Description	French Size	Units per Box	
<b>Accessories</b>				
	<b>497296</b>	Detachable Hemostasis Valve with Stopcock	4F - 6F	10
	<b>497297</b>	Detachable Hemostasis Valve with Stopcock	7F - 9F	10
	<b>497302</b>	Septishield Catheter Sleeve (80cm) <i>(Not for use with AXcess Introducers)</i>		10

# Drapes and Equipment Covers











# Angiography Drapes

Units per Box: 20 each

Product Code	Description	Drape Material	Length	Absorbent Pad	Fenestrations	Transparent Area
	<b>396201</b> ANGIO DRAPE PK/1 - Adult	SMS	77" x 125" (195 x 317cm)	32" x 44" (81 x 112cm)	5.5" (14cm)	None
	<b>396225</b> ANGIO DRAPE PK/25 - Adult	SMS	77" x 110" (195 x 279cm)	32" x 34" (81 x 86cm)	5.5" (14cm)	22" x 47" (56 x 119cm)
	<b>396208</b> ANGIO DRAPE PK/8 - Adult	Dexter	77" x 110" (195 x 279cm)	32" x 34" (81 x 86cm)	5.5" (14cm)	22" x 110" (56 x 279cm)
	<b>396319</b> ANGIO DRAPE PK/19 - Adult	SMS	84" x 125" (213 x 318cm)	32" x 50" (81 x 127cm)	5.5" (14cm)	22" x 125" (56 x 318cm)
	<b>396100</b> EPD ANGIO DRAPE - Adult	Poly	75" x 125" (191 x 318cm)	22" x 125" (56 x 318cm)	6" (15cm)	22" x 125" (56 x 318cm)
	<b>396110</b> EPD BASIC PAD DRAPE - Adult	Poly	75" x 125" (191 x 318cm)	22" x 125" (56 x 318cm)	6" (15cm)	22" x 125" (56 x 318cm)
	<b>396101</b> EPD ANGIO DRAPE - Neonatal	Poly	75" x 125" (191 x 318cm)	22" x 125" (56 x 318cm)	2" (5cm)	22" x 125" (56 x 318cm)
	<b>396102</b> EPD ANGIO DRAPE - Pediatric	Poly	75" x 125" (191 x 318cm)	22" x 125" (56 x 318cm)	3" (8cm)	22" x 125" (56 x 318cm)




# Specialty Drapes

Units per Box: 20 each

Product Code	Description	Drape Material	Length	Absorbent Pad	Fenestrations	Transparent Area
 <b>397100</b>	EPD PACING DRAPE	Poly	75" x 125" (191 x 318cm)	22" x 125" (56 x 318cm)	6" x 6.75" (15 x 17cm)	22" x 125" (56 x 318cm)
 <b>397201</b>	PACING DRAPE	SMS	77" x 125" (195 x 318cm)	32" x 44" (81 x 112cm)	6" x 6.75" (15 x 17cm)	
 <b>397101</b>	EPD PACING DRAPE - Long	Poly	75" x 150" (191 x 400cm)	22" x 150" (56 x 400cm)	6" x 6.75" (15 x 17cm)	22" x 125" (56 x 318cm)
Product Code	Description	Drape Material	Length	Absorbent Pad	Fenestrations	Units per Box
 <b>396260</b>	General Purpose Utility Drape	SMS	40" x 60" (101 x 152cm)	32" x 34" (81 x 86cm)	6.5" x 9" (16 x 23cm)	20
 <b>396204</b>	Brachial/Axillary/Cerebral Drape	Dexter	38.5" x 60" (98 x 152cm)	32" x 34" (81 x 86cm)	3.5" (9cm)	20
 <b>396302</b>	Minor Procedure Drape with Ripline	Dexter	22" x 22" (56 x 56cm)		3.5" (9cm)	10
 <b>396303</b>	Heart Biopsy Drape with Facial Port	SMS	77" x 100" (195 x 254cm)	32" x 34" (81 x 86cm)	3.5" (9cm)	20
 <b>396301</b>	Special Procedure	Poly	32" x 34" (81 x 86cm)	32" x 34" (81 x 86cm)	3.5" (9cm)	20



# Equipment Covers

Units per Box: 20 each




Product Code	Description	Dimensions Width x Height	Special Features
 <b>396223</b>	Bag Type with Elastic Band	20" x 12" (51 x 30cm)	With Corners
<b>396222</b>	Bag Type with Elastic Band	30" x 17" (76 x 43cm)	With Corners
<b>396003</b>	Bag Type with Elastic Band	30" x 36" (76 x 91cm)	With Corners
<b>396221</b>	Bag Type with Elastic Band	36" x 20" (91 x 51cm)	With Corners
<b>396215</b>	Bag Type with Elastic Band	36" x 28" (91 x 71cm)	With Corners
<b>396005</b>	Bag Type with Elastic Band	36" x 36" (91 x 91cm)	With Corners
<b>396007</b>	Bag Type with Elastic Band	42" x 38" (107 x 96cm)	With Corners
<b>396009</b>	Bag Type with Elastic Band	48" x 36" (122 x 91cm)	With Corners
<b>396011</b>	Dome Style with Elastic Band	14" Depth (35cm)	
<b>396012</b>	Dome Style with Elastic Band	18" Depth (46cm)	
<b>396013</b>	Dome Style with Elastic Band	22" Depth (56cm)	
<b>396397</b>	Dome Style with Elastic Band	27" Depth (69cm)	
 <b>396210</b>	Fluoroscope Cover	36" x 44" (91 x 111cm)	Adhesive on four sides
 <b>396348</b>	Back Stand Table Cover	44" x 76" (111 x 193cm)	

Units per Box: 40 each





# Arterial Needles




# Arterial Needles

Product Code	Description	Gauge	Length (cm)	For Guide wire sizes	Special Feature	Units per Box
						
<b>195312</b>	AMC/3 - Arterial Needle with removable Needle Adapter	18	7 cm	0.032" - 0.038"		10
<b>195315</b>	AMC/3 - Arterial Needle with removable Needle Adapter	19	3.8 cm	0.025" - 0.028"		10
<b>195311</b>	AMC/3 - Arterial Needle with removable Needle Adapter	21	2.5 cm	0.018" - 0.021"		10
						
<b>195317</b>	AMC/4 - Clear Flashback Needle	18	3.8 cm	0.035" - 0.038"		10
<b>195318</b>	AMC/4 - Clear Flashback Needle	18	7 cm	0.035" - 0.038"		10
<b>195319</b>	AMC/4 - Clear Flashback Needle	19	3.8 cm	0.025" - 0.028"		10
<b>195355</b>	AMC/4 - Clear Flashback Needle	19	7 cm	0.032" - 0.035"		10
<b>195322</b>	AMC/4 - Clear Flashback Needle	21	2.5 cm	0.018" - 0.021"		10
<b>195323</b>	AMC/4 - Clear Flashback Needle	21	5.1 cm	0.018" - 0.021"		10
<b>195310</b>	AMC/4 - Clear Flashback Needle	21	7 cm	0.018" - 0.021"	Extra Thin Wall	10
						
<b>195352</b>	AMC/4 - Winged Arterial Needle	18	7 cm	0.035" - 0.038"		10
<b>195330</b>	AMC/4 - Winged Arterial Needle	19	7 cm	0.032" - 0.035"	Extra Thin Wall	10

# Arterial Needles

Product Code	Description	Gauge	Length (cm)	For Guide wire sizes	Units per Box
					
<b>195326</b>	<b>Percutaneous Entry Needle</b> <i>Two-Part Thin Wall Tapered Cannula and Inner Arterial Beveled Stylet</i>	18	7.3 cm	0.032" - 0.038"	10
<b>195351</b>	<b>Percutaneous Entry Needle With Vented Stylet</b> <i>Two-Part Thin Wall Tapered Cannula and Vented Arterial Beveled Stylet Needle</i>	18	7.3 cm	0.032" - 0.038"	10
					
<b>195336</b>	<b>Cournand Style Needle</b> <i>Three-Part: Outer Thin Wall Tapered Cannula, Inner Arterial Beveled, and Round ended Obturator</i>	18	5.4 cm	0.032" - 0.038"	10
					
<b>195350</b>	<b>Modified Cournand Style Needle</b> <i>Two-Part: Outer Thin Wall Tapered Can-nula and Inner Arterial Beveled Needle</i>	18	5.4 cm	0.032" - 0.038"	10
					
<b>S03018-275</b>	Safety Arterial Needle	18	7.3 cm	0.032" - 0.038"	25
<b>S03018-275W</b>	Safety Arterial Needle with Wings	18	7.3 cm	0.032" - 0.038"	25

# Guide Wires



# Guide Wires - Fixed Core Straight



Product Code	Diameter	Length (cm)	Tip Curve Radius	Units per Box
<b>Stainless Steel</b>				
395451	.018 inch	40cm	Straight	10
395200	.018 inch	80cm	Straight	10
395201	.018 inch	100cm	Straight	10
395536	.021 inch	25cm	Straight	10
395513	.021 inch	40cm	Straight	10
395203	.021 inch	80cm	Straight	10
395205	.021 inch	125cm	Straight	10
395509	.021 inch	150cm	Straight	10
395206	.025 inch	40cm	Straight	10
395711	.035 inch	75cm	Straight	10
<b>PTFE Coated Stainless Steel</b>				
395225	.025 inch	150cm	Straight	10
395687	.035 inch	150cm	Straight	10
395234	.038 inch	150cm	Straight	10
<b>Exchange Wires - PTFE Coated</b>				
395367	.025 inch	260cm	Straight	10
395369	.035 inch	260cm	Straight	10
395370	.038 inch	260cm	Straight	10
<b>Mandrel Wires for Micro-Puncture Technique</b>				
395125	.018 inch	45cm	Straight	SS 10
395126	.018 inch	45cm	Straight	Nitinol 10
<b>Mullins Tip Deflector Wires - Stainless Steel</b>				
395945	.018 inch	100cm	Straight	10
395948	.021 inch	100cm	Straight	10

# Guide Wires - Fixed Core J Tip




Product Code	Diameter	Length (cm)	Tip Curve Radius	Special Feature	Units per Box
<b>PTFE Coated Stainless Steel</b>					
395420	.025 inch	150cm	1.5mm		10
395787	.035 inch	150cm	1.5mm		10
395523	.021 inch	150cm	3mm		10
395559	.025 inch	150cm	3mm		10
395600	.032 inch	150cm	3mm		10
395282	.035 inch	150cm	3mm		10
388282	.035 inch	150cm	3mm	Classic	10
393282	.035 inch	150cm	3mm	Flexible	10
393182	.035 inch	150cm	3mm	Z-Tip	10
390282	.035 inch	150cm	3mm	XTB	10
393186	.035 inch	175cm	3mm	Z-Tip	10
390186	.035 inch	175cm	3mm	XTB	10
395284	.038 inch	150cm	3mm		10
388284	.038 inch	150cm	3mm	Classic	10
393184	.038 inch	150cm	3mm	Z-Tip	10
390284	.038 inch	150cm	3mm	XTB	10
393187	.038 inch	175cm	3mm	Z-Tip	10
<b>PTFE Coated Stainless Steel Double Ended</b>					
395998	.035 inch	50cm	Str & 3mm		10
395919	.035 inch	150cm	Str & 1.5mm		10
<b>Exchange Wires - PTFE Coated</b>					
395585	.025 inch	260cm	3mm		10
395722	.035 inch	200cm	3mm		10
395794	.035 inch	260cm	3mm		10

# Guide Wires - Specialty



Product Code	Diameter	Length (cm)	Tip Curve Radius	Special Feature	Units per Box
<b>Stainless Steel</b>					
395510	.021 inch	40cm	3mm		10
395573	.025 inch	40cm	3mm		10
395271	.035 inch	40cm	3mm		10
395654	.035 inch	50cm	3mm		10
395752	.035 inch	75cm	3mm		10
<b>Stainless Steel Double Ended</b>					
395415	.025 inch	65cm	Str & 3mm		10
395539	.032 inch	65cm	Str & 3mm		10
395680	.035 inch	45cm	Str & 3mm		10
395899	.035 inch	65cm	Str & 3mm		10
395868	.035 inch	80cm	Str & 3mm		10
<b>Bentson Cerebral Wires - PTFE Coated</b>					
395930E	.035 inch	150cm			10
395170	.035 inch	180cm			10
390033	.035 inch	260cm			10
395933C	.038 inch	150cm			10
394930D	.035 inch	150cm		Heparin	10
<b>Movable Core PTFE Coated</b>					
395258	.035 inch	150cm	Straight		10
395309	.032 inch	150cm	3mm		10
395721	.035 inch	145cm	3mm	Tapered	10
395320	.035 inch	150cm	3mm		10
388300	.035 inch	150cm	3mm	Tapered	10
395312	.035 inch	150cm	3mm	Classic	10
394694	.035 inch	150cm	3mm	Heparin	10
395315	.038 inch	150cm	3mm		10

# Vessel Dilators




# Vessel Dilators



Product Code	French Size	Length (cm)	Guide Wire Diameter/length	Units Per Box
<b>Feathered Tip</b>				
392354	4F	19.8cm	Accepts 0.018" - 0.021"	10
392355	5F	19.8cm	Accepts 0.018" - 0.021"	10
392356	6F	19.8cm	Accepts 0.018" - 0.021"	10
<b>Standard Taper</b>				
392299	4F	13.2cm	Accepts 0.032" - 0.035"	10
392310	5F	19.8cm	Accepts 0.035" - 0.038"	10
392305	6F	19.8cm	Accepts 0.035" - 0.038"	10
392306	7F	19.8cm	Accepts 0.035" - 0.038"	10
392307	8F	19.8cm	Accepts 0.035" - 0.038"	10
392309	9F	19.8cm	Accepts 0.035" - 0.038"	10

# Manifolds



# Manifolds - 500psi with Rotating Adapter



**500psi**

2 Port  
Open Body

3 Port  
Open Body

4 Port  
Open Body

5 Port  
Block Body

3 Port  
Block Body

4 Port  
Block Body

Units per Box

**Right Orientation**

ON Handles	193101	193117	193133	193109	193125	193149	25
OFF Handles	193103	193119	193135	193111	193127	193151	25



**500psi**

2 Port  
Open Body

3 Port  
Open Body

4 Port  
Open Body

2 Port  
Block Body

3 Port  
Block Body

4 Port  
Block Body

Units per Box

**Left Orientation**

ON Handles							25
OFF Handles	193107	193123			193131		25

**500psi**

2 Port  
Block Body

3 Port  
Block Body

Units per Box

**Tubing Extension 6" (15cm)**



**Right Orientation**

ON Handles		193145	193141	25
OFF Handles		193146	193142	25

**1050 PSI**

Product Code

Description

Units per Box

**High Pressure 3-Way Stopcock with Rotating Adapter**




196590

OFF Handle

25

# Control and Medication Syringes



# Control Syringes



10mL      12mL



10mL Soft Touch      12mL Soft Touch



10mL Dosage Control      12mL Dosage Control

Units per Box

## Thumb Ring/Finger Ring

Rotating Adapter	193206	193210	193249	193222	25
Fixed Male Luer	193205	193209	193248	193221	25



## Thumb Ring/Finger Grip

Rotating Adapter	193214	25
Fixed Male Luer	193213	25



## Palm Pad/Finger Ring

Rotating Adapter	193212	193243	25
Fixed Male Luer	193207	193211	25



## Palm Pad/Finger Grip

Rotating Adapter	193237	193239	25	
Fixed Male Luer	193203	193236	193238	25

Product Code

Description

Units per Box



194100

### EKM Inflation Device

30 Atm inflation device with analog display and 3-way stopcock

5

194101

### EKM Inflation Device Kit

30 Atm inflation device with analog display and 3-way stopcock. Kit includes hemostasis valve, guide wire insertion tool, and torque device.

5

# Medication Syringes - 10mL

Blue      Yellow      Red      White      Light Green      Dark Green      Purple


## Unprinted Syringes



<b>382001</b>	<b>382011</b>	<b>382021</b>	<b>382031</b>	<b>382041</b>	<b>382051</b>	<b>382061</b>
---------------	---------------	---------------	---------------	---------------	---------------	---------------

Packaged 25 Units per Box. Each package includes clear adhesive labels, printed in black with common injectables. Labels: Lidocaine, Heparin/Saline, Heparin, Contrast, Nitroglycerin and blank (2ea)

Blue      Yellow      Red      White      Light Green      Dark Green      Purple



## Pre-Printed Syringes

<i>Contrast</i>	<b>382005</b>	<b>382015</b>	<b>382025</b>	<b>382035</b>	<b>382045</b>	<b>382055</b>	<b>382065</b>
<i>Lidocaine</i>	<b>382003</b>	<b>382013</b>	<b>382023</b>	<b>382033</b>	<b>382043</b>		
<i>Saline</i>	<b>382006</b>	<b>382016</b>	<b>382026</b>	<b>382036</b>	<b>382046</b>		
<i>Nitro</i>	<b>382002</b>	<b>382012</b>	<b>382022</b>	<b>382032</b>		<b>382052</b>	
<i>Heparin</i>	<b>382004</b>		<b>382024</b>	<b>382034</b>			
<i>Flush</i>				<b>382037</b>			
<i>Atropine</i>				<b>382038</b>			
<i>Arterial</i>	<b>382000</b>						
<i>Venous</i>				<b>382039</b>			

Packaged 25 Units per Box.

Product Code	Description	Units per Box
--------------	-------------	---------------

## Medication Labeling Kit





<b>197823</b>	<b>MLK Label Set with Pen</b>	<b>5</b>
---------------	-------------------------------	----------

Each set contains 40 (0.04" x 1.0") pre-printed, moisture resistant labels and a permanent marker.

# Fluid Management Tubing and Accessories










# High Pressure Injection Lines - 1200psi

	10" Length (25cm)	20" Length (50cm)	30" Length (76cm)	48" Length (122cm)	60" Length (152cm)	72" Length (183cm)	96" Length (244cm)	Units per Box
								
<b>Flexible Braided</b>								
Male to Female	193358	193359	193360	193361	193376	193377	193381	25
Rotating Adapter	193362	193363	193364	193365	193375	193378	193382	25
								
<b>Clear PVC</b>								
Male to Female	193350	193351	193352	193353	193374	193366		25
Rotating Adapter	193354	193355	193356	193357	193368	193373		25



# Pressure Monitoring Lines

	6" Length (15cm)	12" Length (30cm)	24" Length (61cm)	36" Length (91cm)	48" Length (122cm)	60" Length (152cm)	72" Length (183cm)	96" Length (244cm)
<b>Clear PVC</b>								
Male to Female	040100005A	040101005A	040102005A	040103005A	040104005A	040105005A	040106005A	040108005A
Male to Male	040200005A	040201005A	040202005A	0402030005A	040204005A	040205005A	040206005A	040208005A


# Contrast Management Systems and Spikes

Product	Description	Units per Box
 <b>197652</b>	<b>Contrast Control System</b> <i>With 72" (183cm) large bore tubing, one way check valve, 30ml vented Chamber, 3-way stopcock and fixed male luers</i>	20
 <b>197653</b>	<b>Contrast Control System</b> <i>With 60" (152cm) large bore tubing, Two 1-way check valves, 30ml Burette Chamber, 1-way stopcock, 3-way stopcock and fixed male luers</i>	20
 <b>197641</b>	<b>72" Fluid Delivery System</b> <i>72" (183cm) administration set with Two 1-way check valves and 30ml luer lock syringe</i>	20
 <b>197643</b>	<b>IV Spike 6"</b> <i>Vented Contrast management spike with large bore tubing, (For use with 197652)</i>	20
 <b>197640</b>	<b>IV Spike 6"</b> <i>Vented Contrast management spike with 1-way stopcock. (For use with 197653)</i>	20
 <b>039821599A</b>	<b>Macro IV Spike 60" (152cm)</b> <i>Non-vented chamber with 15 Micron screen filter</i>	10
 <b>039829999A</b>	<b>Contrast Media IV Spike</b> <i>72" (183cm) Contrast management spike with vented universal spike</i>	10






# Waste Management Systems

Product	Description	Units per Box
 <b>197530</b>	<b>Argo-Bagz™ with Universal Spike</b> 1500 ml PVC bag, 72" (183cm) large bore tubing with clasp and Dual Check Valve	20
 <b>197540</b>	<b>Argo-Bagz™ with Universal Drip Chamber</b> 1500 ml PVC bag, 72" (183cm) large bore tubing and Dual Check Valve	20
 <b>197560</b>	<b>Argo-Bagz™ with One Way Valve</b> 1500 ml PVC bag, 72" (183cm) large bore tubing, and One way Check Valve	20
 <b>197580</b>	<b>Argo-Bagz™ with Dual Spikes</b> 1500 ml PVC bag, 72" (183cm) large bore tubing, 72" (183cm) Vented universal spike, 72" (183cm) Contrast Spike and Dual Check Valve with 3-way Stopcock	20
 <b>197581</b>	<b>Argo-Bagz™ with Dual Spikes</b> 1500 ml PVC bag, 72" (183cm) large bore tubing, 72" (183cm) Vented universal spike, 72" (183cm) Contrast Spike and 2 Dual Check Valves	20
 <b>197590</b>	<b>Argo-Bagz™ with Male Extension</b> 1500 ml PVC bag, 72" (183cm) large bore tubing, 36" (91cm) extension with male connector and Dual Check Valve	20
 <b>197611</b>	<b>Argo-Bagz™ with Dual Spikes</b> 1500 ml PVC bag, 72" (183cm) large bore tubing, 2 each 72" (183cm) Vented universal spikes, and Dual Check Valve with 3-way Stopcock	20

# Pressure Transducers




# ArgoTrans™ Transducers

Product Code	Description	Cable Length	Stopcock	Units per Box
<b>Single Disposable Transducers</b>				
	<b>041516504A</b> ArgoTrans with Male/Female Luers	14" (35cm)		20
	<b>041517504A</b> ArgoTrans with 1-way Stopcock	42" (107cm)	1-Way	20
	<b>041587505A</b> ArgoTrans with 1-way Stopcock	60" (152cm)	1-Way	20
	<b>041589507A</b> ArgoTrans with 4-way & 1-Way Stopcocks	14" (35cm)	1-Way & 4-Way	20
	<b>041596505A</b> ArgoTrans with integral 4-way Stopcock and Flush Device	14" (35cm)	4-Way	20

Product Code	Primary Tubing Length	Arterial Extension Tubing Length	Drip Chamber	Units per Box
--------------	-----------------------	----------------------------------	--------------	---------------


**Patient Mount Disposable Transducer Kit with 3mL/hr Flush Device and 4-Way Stopcocks for Continuous Monitoring**

**041582505A** 9" (23cm) Macro-Drip 10




**Neonatal Disposable Transducer Kit with 30mL/hr Flush Device and Double 4-Way Stopcocks for Continuous Monitoring**

**041588505A** 36" (91cm) 10




# ArgoTrans™ Transducers Kits

Product Code	Primary Tubing Length	Arterial Extension Tubing Length	Drip Chamber	Units per Box
<b>Single Pressure Disposable Transducer Kits with Drip Chamber, 3mL/hr Flush Device and 4-way Stopcocks for Continuous Monitoring</b>				
041583505A	48" (122cm)	9" (23cm)	Macro-Drip	10
041584505A	48" (122cm)	9" (23cm)	Micro-Drip	10
041585505A	72" (183cm)	9" (23cm)	Macro-Drip	10
041586505A	72" (183cm)	9" (23cm)	Micro-Drip	10




## Double Pressure Disposable Transducer Kits with Bifurcated Drip Chamber, Two 3mL/hr Flush Device and 4-way Stopcocks for Continuous Monitoring

041589505A	72" (183cm)	12" (30cm)	Macro-Drip	10
------------	-------------	------------	------------	----




## Triple Pressure Disposable Transducer Kits with Trifurcated Drip Chamber, 3mL/hr Flush Device and 4-way Stopcocks for Continuous Monitoring (Color coded tubing 72" Red, 48" Blue & 48" Yellow)

041591505A	72" (183cm)	12" (30cm)	Macro-Drip	10
------------	-------------	------------	------------	----




# KeepSafe Closed Sampling Transducer Kits

Product Code	Primary Tubing Length	Arterial Extension Tubing Length	Drip Chamber	Sampling Stopcocks	Units per Box
<b>KeepSafe Closed Sampling Transducer Kit with Drip Chamber, 3mL/hr Flush Device, Covered Reservoir Syringe, Sampling Stopcock, Arterial extension and 4-way stopcocks</b>					
041592505A	48" (122cm)	9" (23cm)	Macro-Drip	1	10
041593505A	72" (183cm)	9" (23cm)	Macro-Drip	1	10
041597505A	48" (122cm)	9" (23cm)	Macro-Drip	2	10





**Triple Pressure KeepSafe Closed Sampling Kits with Trifurcated Drip Chamber, Three 3mL/hr Flush Device and 4-way Stopcocks for Continuous Monitoring**  
*(Color coded tubing 72" Red, 48" Blue & 48" Yellow)*

041594505A	72" (183cm)	12" (30cm)	Macro-Drip	1	10
------------	-------------	------------	------------	---	----



# Pressure Transducer Accessories












# Pressure Monitoring Lines







Product Code	Length	Tube Material	Units per Box
<b>Pressure Monitoring Tubing with Female/Male Luers</b>			
040100005A	6" (15cm)	PVC	25
040101005A	12" (30cm)	PVC	25
040102005A	24" (61cm)	PVC	25
040103005A	36" (91cm)	PVC	25
040104005A	48" (122cm)	PVC	25
040105005A	60" (152cm)	PVC	25
040106005A	72" (183cm)	PVC	25
040108005A	96" (244cm)	PVC	25
040100001A	6" (15cm)	Polyethylene	25
040101001A	12" (30cm)	Polyethylene	25
040102001A	24" (61cm)	Polyethylene	25
040103001A	36" (91cm)	Polyethylene	25
040104001A	48" (122cm)	Polyethylene	25
040106001A	72" (183cm)	Polyethylene	25
040108001A	96" (244cm)	Polyethylene	25
<b>Pressure Monitoring Tubing with Male//Male Luers</b>			
040200005A	6" (15cm)	PVC	25
040201005A	12" (30cm)	PVC	25
040202005A	24" (61cm)	PVC	25
040203005A	36" (91cm)	PVC	25
040204005A	48" (122cm)	PVC	25
040205005A	60" (152cm)	PVC	25
040206005A	72" (183cm)	PVC	25
040208005A	96" (244cm)	PVC	25
040203001A	36" (91cm)	Polyethylene	25
040204001A	48" (122cm)	Polyethylene	25
040206001A	72" (183cm)	Polyethylene	25
040208001A	96" (244cm)	Polyethylene	25













# Stopcocks and Flush Devices

	Product Code	Description	Units per Box
	041209001A	4-Way Stopcock with Male Luer Slip	25
	041210002A	4-Way Stopcock with Male Luer Lock	25
	041216000A	1-Way Stopcock with Male Luer Lock	25
	041220001A	Double 4-Way Stopcock with Male Luer Lock	25
	041240002A	4-Way Stopcock with Swivel Male Luer	25
	041240104A	Elcam MRVLS Sampling Stopcock with Red Cap	25
	041115002A	6" (15cm) PVC Monitoring Line with 4-Way Stopcock	25
	041112002A	12" (30cm) PVC Monitoring Line with 4-Way Stopcock	25
	041104002A	48" (122cm) PVC Monitoring Line with 4-Way Stopcock	25
	041107002A	84" (218cm) PVC Monitoring Line with 4-Way Stopcock	25
	041300003A	3mL Flush Device (Male port Out)	25
	041300013A	3mL Flush Device (Male and Female ports Out)	25


# Output Bridges and Pole Mounted Manifolds

Product Code	Description	Units per Box
<b>Cardiac Output Bridges</b>		
	<b>041400350A</b> In-Line Cardiac Output Bridge with 4-Way Stopcock and Color Coded (PA & CVP) Ports	10
	<b>041400351A</b> Vertical Cardiac Output Bridge with 4-Way Stopcock and Color Coded (PA & CVP) Ports	10
	<b>041400352A</b> In-Line Cardiac Output Bridge with 3-Way Stopcock and Color Coded (PA & CVP) Ports	10
Product Code	Description	Units per Box
<b>Manifolds - Pole Mounting</b>		
	<b>041273501A</b> 3-Port Manifold, All Female Ports and Color Coded Handles	20
	<b>041275501A</b> 5-Port Manifold, All Female Ports and Color Coded Handles	20
	<b>041278501A</b> 5-Port Manifold, Female Ports (6) and Malte Luer Lock end port	20




# Monitoring Accessories

	Product Code	Description	Units per Box
	040183000A	Double Female Luer Lock Adapter, Clear Polycarbonate	25
	040184000A	Double Male Luer Lock Adapter, Clear Polycarbonate	25
	041280000A	Rotating Adapter, Male L/L And Female L/L, Clear Polycarbonate	50
	041281000A	Rotating Adapter, Double Male L/L, Clear Polycarbonate	50
	041290000A	Zeroport Plug Male L/L With Extra Long Threads	100
	041290001A	Male L/L Plug With Extension	100
	041295000A	Male/Female L/L Plug	50
	041519000A	Manifold Holder, Plastic (Non-Sterile)	1
	041520000A	3-Position Half-Moon Bracket And Plastic Manifold Holder (Non-Sterile)	1
	041520100A	4-Position Half-Moon Bracket And Plastic Manifold Holder (Non-Sterile)	1

# Monitoring Catheter Kits & Trays



# Arterial Catheter Kits and Trays

	Product Code	Catheter	Length	Guide Wire Diameter/Length	Needle Gauge/Length	Units per Box
<b>Arterial Catheter Mini-Kits</b>						
	498107	20G	3" (7.6cm)	0.021" x 25cm	21Gx1" (2.5cm)	10
	498101	20G	6" (15cm)	0.021" x 40cm	21Gx2" (5.1cm)	10
	498111	18G	3" (7.6cm)	0.028" x 25cm	19Gx1.5" (3.8cm)	10
	498105	18G	6" (15cm)	0.028" x 40cm	19Gx2.25" (5.7cm)	10
<b>Kit Contents:</b>						
1- Radiopaque PTFE Catheter With Female Luer Hub						
1- Guide Wire (FC/SS/ST)						
1- Arterial Needle, Clear Flashback Hub						
<b>Arterial Catheter Kits</b>						
	498106	20G	3" (7.6cm)	0.021" x 25cm	21Gx1" (2.5cm)	10
	498100	20G	6" (15cm)	0.021" x 40cm	21Gx2" (5.1cm)	10
	498110	18G	3" (7.6cm)	0.028" x 25cm	19Gx1.5" (3.8cm)	10
	498104	18G	6" (15cm)	0.028" x 40cm	19Gx2.25" (5.7cm)	10
<b>Kit Contents:</b>						
1- Disposable Towel		1- Radiopaque PTFE Catheter With Female Luer Hub				
1- 5cc Syringe		1- 000 Silk Suture				
1- #11 Mini Scalpel		1- Arterial Needle, Clear Flashback Hub				
3- 4" x 4" Gauze Sponges		1- Guide Wire (FC/SS/ST)				
1- 9" (23cm) PVC Extension with Stopcock						
<b>Arterial Catheter Tray</b>						
	498411B	20G	3" (7.6cm)	0.021" x 25cm	21Gx1" (2.5cm)	10
	498410B	20G	6" (15cm)	0.021" x 40cm	21Gx2" (5.1cm)	10
	498413A	18G	3" (7.6cm)	0.028" x 25cm	19Gx1.5" (3.8cm)	10
	498412B	18G	6" (15cm)	0.028" x 40cm	19Gx2.25" (5.7cm)	10
<b>Tray Contents:</b>						
1- Fenestrated Drape 22" x 22" (56cm x 56cm) With Ripline		1- #11 Mini Scalpel		1- Radiopaque PTFE Catheter With Female Luer Hub		
3- 4" x 4" Gauze Sponges		1- Arterial Needle, Clear Flashback Hub		1- 5cc Syringe		
4- 2" x 2" Gauze Sponges		1- Guide Wire (FC/SS/3mmJ & ST)		1- 25g x 5/8" (1.6cm) Needle		
1- Povidone/Iodine Ointment		1- 5cc Ampule 1% Lidocaine HCL				
1- 000 Silk Suture						

# Thermodilution and Pulmonary Artery Catheters

Product Code	French Size	Lumen	Length	Guide Wire	Special Feature	Units per Box
--------------	-------------	-------	--------	------------	-----------------	---------------

## Pulmonary Artery Catheters (Heparin Coated with Latex Balloon)



<b>497527</b>	6F	2	110cm	Accepts .018" Wire		5
<b>497523</b>	7F	2	110cm	Accepts .025" Wire		5
<b>497540</b>	7F	2	110cm	Accepts .035" Wire		5

Product Code	French Size	Lumen	Length	Guide Wire	Special Feature	Units per Box
--------------	-------------	-------	--------	------------	-----------------	---------------

## Thermodilution Catheters - Standard (Heparin Coated with Latex Balloon)



<b>497537</b>	6F	4	110cm	Accepts .018" Wire	Safety Wedge/Sleeve	5
<b>497521</b>	7F	4	110cm	Accepts .025" Wire	Safety Wedge/Sleeve	5
<b>497539</b>	7F	4	110cm	Accepts .025" Wire	Safety Wedge	5
<b>497526</b>	7F	4	110cm	Accepts .025" Wire		5
<b>497520</b>	7.5	5	110cm	Accepts .025" Wire	Safety Wedge/Sleeve	

## Thermodilution Catheters - S-Curve Tip (Heparin Coated with Latex Balloon)



<b>497534</b>	7F	4	110cm	Accepts .025" Wire	Safety Wedge/Sleeve	5
<b>497536</b>	7F	4	110cm	Accepts .025" Wire		5

## Thermodilution Catheters - Soft Body (Heparin Coated with Latex Balloon)

<b>497531</b>	7.5F	5	110cm	Accepts .025" Wire	Safety Wedge/Sleeve	5
---------------	------	---	-------	--------------------	---------------------	---

# Temporary Pacing Electrode Trays



Product Code	French Size	Length	Special Feature	Units Per Box
<b>Balloon Flow Assisted Bipolar Electrode Tray</b>				
008566	5F	110cm		5
008560	5F	110cm	Without Lidocaine	5
<b>Special Care Bipolar Electrode Tray</b>				
008416	5F	115cm		5

**Tray Contents:**

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>1- #11 Mini Scalpel</li> <li>1- 6F Vessel Dilator</li> <li>1- Fenestrated Drape</li> <li>5- 4" x 4" Gauze Sponges</li> <li>1- 25g x 5/8" (1.6cm) Needle</li> <li>1- 21g x 1.5" (3.8cm) Needle</li> <li>2- 10cc Syringe L/L</li> <li>1- 5cc Syringe L/L</li> <li>1- TPE Catheter Kit</li> <li>1- 000 Silk Suture</li> </ul> | <ul style="list-style-type: none"> <li>2- Disposable Towels</li> <li>1- Guide Wire (0.035" x 45cm FC/SS/3mmJ &amp; Straight)</li> <li>1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub</li> <li>1- Povidone/Iodine Prep Swab</li> <li>1- 6F Tuohy Borst Introducer</li> <li>1- Povidone/Iodine Ointment</li> <li>1- 5cc Ampule Lidocaine HCL 1%</li> <li>1- Septishield Catheter Sleeve</li> <li>1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter</li> </ul> |
|---|--|

# Central Venous Catheter Kits and Trays



Product Code	Catheter	Length	Guide Wire Diameter/Length	Units per Box
<b>Blitt Basic Kit/1 Central Venous Catheter Kit</b>				
497681	16G	8" (20cm)	0.035" x 45cm	10
497682	14G	8" (20cm)	0.035" x 45cm	10
<b>Kit Contents:</b>				
1- Disposable Towel		1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter		
1- 5cc Syringe		1- Radiopaque Polyurethane CVC Catheter		
1- #11 Mini Scalpel		1- 000 Silk Suture		
3- 4" x 4" Gauze Sponges		1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub		
1- 22g x 1.5" (3.8cm) Needle		1- Guide Wire (Double ended 3mmJ & Straight)		
1- 7F Vessel Dilator				
<b>Blitt Maxi Kit/1 Central Venous Catheter Kit</b>				
497683	16G	8" (20cm)	0.035" x 45cm	10
497661	16G	6" (15cm)	0.035" x 45cm	10
497660	14G	6" (15cm)	0.035" x 45cm	10
<b>Kit Contents:</b>				
1- Fenestrated Drape		2- 5cc Syringe		
1- 25g x 5/8" (1.6cm) Needle		1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub		
1- Povidone/Iodine Ointment		3- Chloraprep™ Prep Swab		
4- 4" x 4" Gauze Sponges		1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter		
1- 22g x 1.5" (3.8cm) Needle		1- Guide Wire (Double ended 3mmJ & Straight)		
1- 7F Vessel Dilator		1- #11 Mini Scalpel		
1- 000 Silk Suture		1- Radiopaque Polyurethane CVC Catheter		

# Index

## Alpha-Numeric



# Argon Ordering Information

## Mailing Address:

Customer Service Department  
Argon Medical Devices, Inc.  
1445 Flat Creek Road  
Athens, TX 75751

Direct orders to Customer Service Department, weekdays (excluding holidays) between 7:00 am and 5:00 pm (CST)

**Toll Free: 800-927-4669**

**Local: 903-675-9321**

**Fax: 903-677-9396**

## Terms of Sale:

Net 30 days

Normal ground delivery: 5-7 business days. All orders are F.O.B. shipping point.

Handling Fee: \$4 USD will be assessed for processing each shipment.

Shipments are invoiced at the time the shipment leaves Argon's facility and ownership of all products is considered transferred at this time. Prices are subject to change without notice.

## Invoice Billings:

1. Invoice payments must be remitted per terms. Cash discounts are not allowed on any portion of the invoice. 2. Any disputed items should be explained by forwarding information with remittance or documentation mailed to Attn: A/R Department, Argon Medical Devices, 1445 Flat Creek road, Athens, TX 75751. 3. Balances remaining unpaid at the due date are subject to a finance charge of 1½% per month or the highest rate permitted by law, whichever is lower, until paid. Any discounts or credits will be applied against delinquent balances before reimbursement is made. If payments become more than sixty (60) days in arrears and it becomes necessary for Argon Medical to solicit outside help in order to collect, the customer agrees to pay all valid, reasonable attorney and/or collection fees.

## Disclaimer of Warranty & Limitation of Remedy:

THERE IS NOT EXPRESS OR IMPLIED WARRANTY, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON ARGON PRODUCTS, INCLUDING THE PRODUCTS DESCRIBED IN THIS PUBLICATION. Under no circumstances shall Argon be liable for any direct, incidental or consequential damages other than as expressly provided by specific laws. No person has the authority to bind to any representation or warranty except as specifically set forth herein.

## Product Return Policy Authorization:

Before an Argon product is returned for credit consideration, approval for the return must be obtained. Authorization of a return can be requested by contacting Argon Customer Service. Upon the approval of the return, customer service will provide an authorization number.

## Return of Products:

Once an authorization number is assigned, products should be returned PREPAID in proper protective packaging to:

## Returns Goods, QA/RA

The following information must be provided with returned products before credit can be considered:

- The Argon Return Goods Authorization number
- Customer purchase order number
- The Argon invoice number
- Lot or serial number of each item being returned
- The date the returned merchandise was shipped by Argon
- The quantity shipped
- The catalog number
- Please only return full box quantities

## Credit Policy for Unused Products:

Full Credit for merchandise shipped in error by Argon and merchandise returned at the request of Argon.

Full Credit less service charge of 20% of the billed price (Merchandise returned in re-saleable condition within 90 days after shipping date)

## Merchandise that will not be considered for credit includes:

- Merchandise returned without proper authorization documentation
- Merchandise not in the original sealed package or other Argon standard packaging
- Merchandise in un-saleable condition, as determined by Argon on receipt
- Custom Products can not be returned
- Merchandise held by the customer for more than 90 days

## Standing Purchase Orders:

Argon encourages standing purchase orders. Standing purchase orders are agreements to purchase specific quantities of product for delivery on a predetermined schedule. To establish a standing purchase order, contact customer service or your local Argon Sales Representative.

## Custom Products:

Products not listed in this catalog may be manufactured to your specifications. Contact your local Argon Sales Representative for further information.

## Intended Product Usage/Storage Usage:

All products in this catalog are supplied sterile and intended for single use (Except Reusable Transducer Cables and Brackets). Do not use any product if its sterile package is damaged. Prior to use, refer to Instructions For Use in each Box.

## Storage:

Store products in a cool, dark, dry place. Use sterile products prior to the USE BEFORE date. Do not expose to organic solvents.

## Caution:

Federal (USA) law restricts these devices to sale by or on the order of a physician.

<u>Product Code</u>	<u>Page</u>	<u>Product Code</u>	<u>Page</u>	<u>Product Code</u>	<u>Page</u>	<u>Product Code</u>	<u>Page</u>
008416	N-3	041220001A	M-2	190060	C-1	193210	J-1
008560	N-3	041240002A	M-2	190065	C-1	193211	J-1
008566	N-3	041240104A	M-2	190070	C-1	193212	J-1
039821599A	K-2	041273501A	M-3	190075	C-1	193213	J-1
039829999A	K-2	041275501A	M-3	190080	C-1	193214	J-1
040100001A	M-1	041278501A	M-3	190081	C-1	193217	J-1
040100005A	M-1	041280000A	M-4	190085	C-1	193221	J-1
040101001A	M-1	041281000A	M-4	190086	C-1	193222	J-1
040101005A	M-1	041290000A	M-4	190090	C-1	193236	J-1
040102001A	M-1	041290001A	M-4	190095	C-1	193237	J-1
040102005A	M-1	041295000A	M-4	191545	C-2, D-7	193238	J-1
040103001A	M-1	041300003A	M-2	191598	C-2, D-7	193239	J-1
040103005A	M-1	041300013A	M-2	191645	C-2, D-7	193243	J-1
040104001A	M-1	041400350A	M-3	191698	C-2, D-7	193248	J-1
040104005A	M-1	041400351A	M-3	191745	C-2, D-7	193249	J-1
040105005A	M-1	041400352A	M-3	191798	C-2, D-7	193350	K-1
040106001A	M-1	041516504A	L-1	191845	C-2, D-7	193351	K-1
040106005A	M-1	041517504A	L-1	191898	C-2, D-7	193352	K-1
040108001A	M-1	041519000A	M-4	193101	I-1	193353	K-1
040108005A	M-1	041520000A	M-4	193103	I-1	193354	K-1
040183000A	M-4	041520100A	M-4	193107	I-1	193355	K-1
040184000A	M-4	041582505A	L-1	193109	I-1	193356	K-1
040200005A	M-1	041583505A	L-2	193111	I-1	193357	K-1
040201005A	M-1	041584505A	L-2	193117	I-1	193358	K-1
040202005A	M-1	041585505A	L-2	193119	I-1	193359	K-1
040203001A	M-1	041586505A	L-2	193123	I-1	193360	K-1
040203005A	M-1	041587505A	L-1	193125	I-1	193361	K-1
040204001A	M-1	041588505A	L-1	193127	I-1	193362	K-1
040204005A	M-1	041589505A	L-2	193131	I-1	193363	K-1
040205005A	M-1	041589507A	L-1	193133	I-1	193364	K-1
040206001A	M-1	041591505A	L-2	193135	I-1	193365	K-1
040206005A	M-1	041592505A	L-3	193141	I-1	193366	K-1
040208001A	M-1	041593505A	L-3	193142	I-1	193368	K-1
040208005A	M-1	041594505A	L-3	193145	I-1	193373	K-1
041104002A	M-2	041596505A	L-1	193146	I-1	193374	K-1
041107002A	M-2	041597505A	L-3	193149	I-1	193375	K-1
041112002A	M-2	190020	C-1	193151	I-1	193376	K-1
041115002A	M-2	190030	C-1	193203	J-1	193377	K-1
041209001A	M-2	190031	C-1	193205	J-1	193378	K-1
041210002A	M-2	190040	C-1	193206	J-1	193381	K-1
041216000A	M-2	190051	C-1	193207	J-1	193382	K-1

<u>Product Code</u>	<u>Page</u>	<u>Product Code</u>	<u>Page</u>	<u>Product Code</u>	<u>Page</u>	<u>Product Code</u>	<u>Page</u>
194100	J-1	382011	J-2	392310	H-1	395559	G-2
194101	J-1	382012	J-2	392354	H-1	395573	G-3
195100	D-8	382013	J-2	392355	H-2	395585	G-2
195110	D-8	382015	J-2	392356	H-1	395600	G-2
195310	F-1	382016	J-2	393182	G-2	395654	G-3
195311	F-1	382021	J-2	393184	G-2	395680	G-3
195312	F-1	382022	J-2	393186	G-2	395687	G-1
195315	F-1	382023	J-2	393187	G-2	395711	G-1
195317	F-1	382024	J-2	393282	G-2	395721	G-3
195318	F-1	382025	J-2	394694	G-3	395722	G-2
195319	F-1	382026	J-2	394930D	G-3	395752	G-3
195322	F-1	382031	J-2	395125	G-1	395787	G-2
195323	F-1	382032	J-2	395126	G-1	395794	G-2
195326	F-2	382033	J-2	395170	G-3	395868	G-3
195330	F-1	382034	J-2	395200	G-1	395899	G-3
195336	F-2	382035	J-2	395201	G-1	395919	G-2
195350	F-2	382036	J-2	395203	G-1	395930E	G-3
195351	F-2	382037	J-2	395205	G-1	395933C	G-3
195352	F-1	382038	J-2	395206	G-1	395945	G-1
195355	F-1	382039	J-2	395225	G-1	395948	G-1
196590	I-1	382041	J-2	395234	G-1	395998	G-2
197530	K3	382043	J-2	395258	G-3	396003	E-3
197540	K-3	382045	J-2	395271	G-3	396005	E-3
197560	K-3	382046	J-2	395282	G-2	396007	E-3
197580	K-3	382051	J-2	395284	G-2	396009	E-3
197581	K-3	382052	J-2	395309	G-3	396011	E-3
197590	K-3	382055	J-2	395312	G-3	396012	E-3
197611	K-3	382061	J-2	395315	G-3	396013	E-3
197640	K-2	382065	J-2	395320	G-3	396100	E-1
197641	K-2	388282	G-2	395367	G-1	396101	E-1
197643	K-2	388284	G-2	395369	G-1	396102	E-1
197652	K-2	388300	G-3	395370	G-1	396110	E-1
197653	K-2	390033	G-3	395415	G-3	396201	E-1
197823	J-2	390186	G-2	395420	G-2	396204	E-2
382000	J-2	390282	G-2	395451	G-1	396208	E-1
382001	J-2	390284	G-2	395509	G-1	396210	E-3
382002	J-2	392299	H-1	395510	G-3	396215	E-3
382003	J-2	392305	H-1	395513	G-1	396221	E-3
382004	J-2	392306	H-1	395523	G-2	396222	E-3
382005	J-2	392307	H-1	395536	G-1	396223	E-3
382006	J-2	392309	H-1	395539	G-3	396225	E-1

**Product Code Page**

396260 E-2  
 396301 E-2  
 396302 E-2  
 396303 C-2, E-2  
 396319 E-1  
 396348 E-3  
 396397 E-3  
 397100 E-2  
 397101 E-2  
 397201 E-2  
 490005 D-1  
 490006 D-1  
 490027 D-1  
 490028 D-1  
 490029 D-1  
 496004 D-1  
 496005 C-2, D-1  
 496006 C-2, D-1  
 496007 C-2, D-1  
 496008 C-2, D-1  
 496035 D-1  
 496036 D-1  
 496037 D-1  
 496038 D-1  
 496064 D-1  
 496065 D-1  
 496066 D-1  
 496067 D-1  
 496124 D-4  
 496125 D-4  
 496126 D-4  
 496127 D-4  
 496128 D-4  
 496132 D-5  
 496133 D-5  
 496134 D-5  
 496135 D-5  
 496136 D-5  
 496140 D-4  
 496141 D-4  
 496142 D-4

**Product Code Page**

496150 D-5  
 496151 D-5  
 496152 D-5  
 496153 D-5  
 496154 D-5  
 496155 D-5  
 496156 D-5  
 496181 D-5  
 496182 D-5  
 496183 D-5  
 496308 D-4  
 496309 D-4  
 496310 D-4  
 496311 D-4  
 496312 D-4  
 497296 D-8  
 497297 D-8  
 497302 D-8  
 497422 D-6  
 497423 D-6  
 497424 D-6  
 497425 D-6  
 497426 D-6  
 497432 D-6  
 497433 D-6  
 497434 D-6  
 497435 D-6  
 497436 D-6  
 497452 D-6  
 497453 D-6  
 497454 D-6  
 497455 D-6  
 497456 D-6  
 497520 N-2  
 497521 N-2  
 497523 N-2  
 497526 N-2  
 497527 N-2  
 497531 N-2  
 497534 N-2  
 497536 N-2

**Product Code Page**

497537 N-2  
 497539 N-2  
 497540 N-2  
 497660 N-4  
 497661 N-4  
 497681 N-4  
 497682 N-4  
 497683 N-4  
 497801 D-3  
 497802 D-3  
 497811 D-3  
 497812 D-3  
 497813 D-3  
 497814 D-3  
 497816 D-3  
 497817 D-3  
 497818 D-3  
 497819 D-3  
 497820 D-2  
 497821 D-2  
 497830 C-2, D-2  
 497831 C-2, D-2  
 497840 C-2, D-2  
 497841 C-2, D-2  
 497844 D-2  
 497845 D-2  
 497849 D-2  
 497850 D-2  
 497852 D-2  
 497853 D-2  
 497854 D-2  
 497855 D-2  
 497856 D-2  
 497857 D-2  
 497858 D-2  
 497859 D-2  
 497860 D-2  
 497861 C-2, D-2  
 497862 C-2, D-2  
 497863 D-2  
 497864 D-2

**Product Code Page**

497865 D-2  
 497867 D-2  
 497868 D-2  
 497880 D-3  
 497881 D-3  
 497882 D-3  
 497883 D-3  
 497884 D-3  
 497885 D-3  
 498100 N-1  
 498101 N-1  
 498104 N-1  
 498105 N-1  
 498106 N-1  
 498107 N-1  
 498110 N-1  
 498111 N-1  
 498410B N-1  
 498411B N-1  
 498412B N-1  
 498413A N-1  
 500200 D-8  
 700009 B-2  
 700109 B-2  
 700115 B-2  
 716009 B-2  
 716109 B-2  
 716115 B-2  
 717009 B-2  
 717109 B-2  
 717115 B-2  
 S030318-275 F-2  
 S030318-275W F-2