

Ordering Information

Each Set Includes:

- (1) SKATER Catheter
- (1) Removable pigtail straightener
- (1) Metal Stiffening Cannula
- (1) Flexible Stiffening Cannula (Not available on 6F and 7F)
- (1) Locking system with snap-off tab (not available on non-locking pigtail catheters)

SKATER All Purpose and Nephrostomy Locking (Box of 5)			
CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET
756606015	6F x 15cm	.035"	Yes
756606020	6F x 20cm	.035"	Yes
756606025	6F x 25cm	.035"	Yes
756606035	6F x 35cm	.035"	No
756607015	7F x 15cm	.035"	Yes
756607020	7F x 20cm	.035"	Yes
756607025	7F x 25cm	.035"	Yes
756608015	8F x 15cm	.038"	Yes
756608020	8F x 20cm	.038"	Yes
756608025	8F x 25cm	.038"	Yes
756608030	8F x 30cm	.038"	Yes
756608035	8F x 35cm	.038"	No
756608045	8F x 45cm	.038"	No
756610015	10F x 15cm	.038"	Yes
756610020	10F x 20cm	.038"	Yes
756610025	10F x 25cm	.038"	Yes
756610030	10F x 30cm	.038"	Yes
756610035	10F x 35cm	.038"	No
756610045	10F x 45cm	.038"	No
756610060	10F x 60cm	.038"	No
756612015	12F x 15cm	.038"	Yes
756612020	12F x 20cm	.038"	Yes
756612025	12F x 25cm	.038"	Yes
756612030	12F x 30cm	.038"	Yes
756612035	12F x 35cm	.038"	No
756612045	12F x 45cm	.038"	No
756612060	12F x 60cm	.038"	No
756614015	14F x 15cm	.038"	Yes
756614020	14F x 20cm	.038"	Yes
756614025	14F x 25cm	.038"	Yes
756614030	14F x 30cm	.038"	Yes
756614035	14F x 35cm	.038"	No
756614045	14F x 45cm	.038"	No
756614060	14F x 60cm	.038"	No
756616025	16F x 25cm	.038"	Yes
756616030	16F x 30cm	.038"	Yes
756616035	16F x 35cm	.038"	No
756616045	16F x 45cm	.038"	No
756616060	16F x 60cm	.038"	No

*Flexible stiffening cannula not available on 6F and 7F.

*Additional ordering information is available on the back

SKATER All Purpose and Nephrostomy Non-Locking (Box of 5)			
CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET
756106015	6F x 15cm	.035"	Yes
756106020	6F x 20cm	.035"	Yes
756106025	6F x 25cm	.035"	Yes
756107020	7F x 20cm	.035"	Yes
756107025	7F x 25cm	.035"	Yes
756108020	8F x 20cm	.038"	Yes
756108025	8F x 25cm	.038"	Yes
756108030	8F x 30cm	.038"	Yes
756110020	10F x 20cm	.038"	Yes
756110025	10F x 25cm	.038"	Yes
756110030	10F x 30cm	.038"	Yes
756112020	12F x 20cm	.038"	Yes
756112025	12F x 25cm	.038"	Yes
756114020	14F x 20cm	.038"	Yes
756114025	14F x 25cm	.038"	Yes
756116025	16F x 25cm	.038"	Yes

*Flexible stiffening cannula not available on 6F and 7F.

SKATER Mini-Loop Drainage Catheter w/Locking Pigtail (Box of 5)			
CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET
758606015	6F x 15cm	.035"	Yes
758606020	6F x 20cm	.035"	Yes
758606025	6F x 25cm	.035"	Yes
758607015	7F x 15cm	.035"	Yes
758607025	7F x 25cm	.035"	Yes
758608015	8F x 15cm	.035"	Yes
758608025	8F x 25cm	.035"	Yes
758610015	10F x 15cm	.038"	Yes
758610025	10F x 25cm	.038"	Yes
758612015	12F x 15cm	.038"	Yes
758612025	12F x 25cm	.038"	Yes
758614025	14F x 25cm	.038"	Yes

SKATER Biliary Drainage Catheter w/Locking Pigtail (Box of 5)		
CATALOG #	FRENCH x LENGTH	MAX GUIDEWIRE SIZE
755308040	8F x 40cm	.038"
755310040	10F x 40cm	.038"
755312040	12F x 40cm	.038"

2600 Dallas Parkway
Suite 440
Frisco, TX 75034
Customer Service: 800.927.4669
Fax: 903.677.9396

© 2021 Argon Medical Devices, Inc. All Rights Reserved.
No reproduction of any portion of this work shall be made without prior permission of the author.

Complete Drainage System

SKATER™




ARGON
MEDICAL DEVICES
96-6300-14/0521G

ARGON
MEDICAL DEVICES

The SKATER™ Drainage System provides a complete portfolio of catheters and accessories for all-purpose, nephrostomy, biliary and centesis procedures.

Formulated to promote maximum drainage flow



All-Purpose Drainage

- Catheters for drainage of abscesses and fluid collections
- Choice Lock™ trocar stylet allows for single-handed abscess drainage placement
 - Available in locking and non-locking configurations
 - Expansive size offering of 6F - 16F
 - Can be used safely with ethyl alcohol and doxycycline

Nephrostomy Drainage

- Catheters for drainage of the urinary tract
- Designed for Seldinger technique
 - Expansive size offering of 6F - 14F
 - Can be used safely with ethyl alcohol and doxycycline

Biliary Drainage

- Catheters for drainage of the biliary tract
- Radiopaque marker band provides clear view of the most proximal drainage hole for accurate placement in the biliary tract
 - Additional drainage holes for biliary duct to duodenum drainage
 - Can be used safely with ethyl alcohol and doxycycline

Also Available:

- SKATER Centesis
- SKATER Introducer Sets
- Introducer Sheaths/Needles (ISN)
- Hawkins™ Blunt Needle Access System
- Drainage Bags
- SKATER Fix Catheter Fixation Device

