



**JAWZ™**

The cutting edge in  
endomyocardial biopsies




# A Clean-Cut Approach to High-Quality Specimens<sup>1</sup>




**Perform precise cuts with clean margins**  
Dual cutting head with a proprietary alloy designed to maintain sharpness.

**Gain greater flexibility and ease of use**  
The single-use device has a dependable response and rotating thumb ring for comfort and control.



**Tailored to Your Technique**  
JAWZ configurations include:  
• Formable or Non-Formable  
• Stainless Steel or FEP Coated  
• Straight, Pre-Curved, or Maxi-Curved



**Ergonomic handle**  
maintains desired thumb position while providing comfort and control

**Rigid and flexible metal coil**  
to advance into the tissue and safely navigate to the target



**Advanced dual cutting head**  
designed for sharper, cleaner cuts

# Expect Clean-Cut Samples



## Ordering Information

### Box of 1

| Catalog # | Jaw size | Jaw OD (French) | Jaw volume          | Length | Curve       | Coating | Formability  | Fenestration    |
|-----------|----------|-----------------|---------------------|--------|-------------|---------|--------------|-----------------|
| 190090    | 1.5mm    | 5F              | 1.80mm <sup>3</sup> | 50cm   | Straight    | None    | Formable     | Fenestrated     |
| 190095    | 1.5mm    | 5F              | 1.80mm <sup>3</sup> | 105cm  | Straight    | None    | Formable     | Fenestrated     |
| 190051    | 1.8mm    | 6F              | 2.46mm <sup>3</sup> | 50cm   | Maxi-Curved | FEP     | Formable     | Non-Fenestrated |
| 190060    | 1.8mm    | 6F              | 2.46mm <sup>3</sup> | 50cm   | Straight    | None    | Formable     | Non-Fenestrated |
| 190065    | 1.8mm    | 6F              | 2.46mm <sup>3</sup> | 50cm   | Straight    | None    | Non-Formable | Non-Fenestrated |
| 190066    | 1.8mm    | 6F              | 2.46mm <sup>3</sup> | 70cm   | Straight    | FEP     | Formable     | Non-Fenestrated |
| 190070    | 1.8mm    | 6F              | 2.46mm <sup>3</sup> | 105cm  | Straight    | None    | Formable     | Non-Fenestrated |
| 190075    | 1.8mm    | 6F              | 2.46mm <sup>3</sup> | 105cm  | Straight    | None    | Non-Formable | Non-Fenestrated |
| 190080    | 2.2mm    | 7F              | 4.97mm <sup>3</sup> | 50cm   | Pre-Curved  | FEP     | Formable     | Non-Fenestrated |
| 190081    | 2.2mm    | 7F              | 4.97mm <sup>3</sup> | 50cm   | Maxi-Curved | FEP     | Formable     | Non-Fenestrated |
| 190085    | 2.2mm    | 7F              | 4.97mm <sup>3</sup> | 105cm  | Straight    | None    | Formable     | Non-Fenestrated |
| 190086    | 2.2mm    | 7F              | 4.97mm <sup>3</sup> | 105cm  | Straight    | FEP     | Formable     | Non-Fenestrated |
| 190020    | 2.4mm    | 7.5F            | 5.43mm <sup>3</sup> | 50cm   | Straight    | None    | Formable     | Non-Fenestrated |
| 190030    | 2.4mm    | 7.5F            | 5.43mm <sup>3</sup> | 50cm   | Pre-Curved  | FEP     | Formable     | Non-Fenestrated |
| 190031    | 2.4mm    | 7.5F            | 5.43mm <sup>3</sup> | 50cm   | Maxi-Curved | FEP     | Formable     | Non-Fenestrated |
| 190040    | 2.4mm    | 7.5F            | 5.43mm <sup>3</sup> | 105cm  | Straight    | None    | Formable     | Non-Fenestrated |

1. Jawz Endomyocardial Biopsy Forceps from Argon Medical: A Review of Safety and Efficacy Based on Literature Review and Bench Model Comparison With the Cordis Bioptome. P. Stibbs, J. Reinault, M. Parsons (2023). DOI: 10.14293/S2199-1006.1.SOR-.PPEZMK7.v1



**Retrieve  
High-Quality  
Samples Today**

For full Indication and Risk information, refer to **ArgonMedical.com**

© 2025 Argon Medical Devices, Inc. All rights reserved. Argon Medical Devices, the stylized logo, and the product and service marks mentioned herein are trademarks or registered trademarks of Argon Medical Devices, Inc. in the United States and other countries.

7800 Dallas Parkway, Suite 200, Plano, TX 75024 P: 800.927.4669 F: 903.677.9396

96-6000-07/0125C **Rx Only**

