


# V•Stick™ Vascular Access Sets

## Access with Ease

V•Stick Vascular Access Sets are designed to facilitate a smooth and easy introduction into the vessel for guidewires up to .038" and minimize trauma to the surrounding tissue.





## Ordering Information

V-Stick Vascular Access Sets w/Standard Co-axial Introducer (Box of 10)							
SHEATH		NEEDLE		DILATOR		GUIDEWIRE	
CATALOG#	SIZE	GAUGE X LENGTH	ECHOGENIC	STIFFNESS	LOCKING NUT COLOR	DIAMETER X LENGTH	MATERIAL
391504500E	4F	21G x 7cm	NO	Standard	Red	.018" x 40cm	SS
391505000E	5F	21G x 7cm	YES	Standard	Red	.018" x 40cm	SS
391505500E	5F	21G x 7cm	NO	Standard	Red	.018" x 40cm	SS

V-Stick Vascular Access Sets w/Stiff Co-axial Introducer (Box of 10)							
SHEATH		NEEDLE		DILATOR		GUIDEWIRE	
CATALOG#	SIZE	GAUGE X LENGTH	ECHOGENIC	STIFFNESS	LOCKING NUT COLOR	DIAMETER X LENGTH	MATERIAL
391604000E	4F	21G x 7cm	YES	Stiff	White	.018" x 40cm	SS
391604500E	4F	21G x 7cm	NO	Stiff	White	.018" x 40cm	SS
391605000E	5F	21G x 7cm	YES	Stiff	White	.018" x 40cm	SS

SS = Stainless Steel