

The Future of Portal Vein Access

SCORPION[®]
PORTAL VEIN ACCESS SERIES

A blue silhouette of a scorpion is positioned behind the letter 'S' in the word 'SCORPION'.


ARGON
MEDICAL DEVICES

Revolutionizing Portal Vein Access Procedures

Go beyond the ordinary with Argon Medical's innovative Scorpion® Portal Vein Access Series.

Greater Visibility.

- **Echogenic tips** enable visibility of the needle or stylet on various imaging modalities for easy detection while piercing through liver to portal vein.
- **Alignment marker** helps identify when the needle/stylet is flush with the end of the stiffening cannula.
- **Radiopaque band** on the stiffening cannula provides amplified visibility under fluoroscopy to assist with locating and orienting each component.



Alignment Marker

Improved Durability.

- **Needle/stylet catheters** are constructed of PEEK material instead of traditional PTFE to reduce the risk of fishmouthing and buckling when advancing through cirrhotic parenchyma.



PEEK Needle/Stylet Catheter

Unparalleled Steerability.

- **Independently steer the curved needle/stylet and the curved stiffening cannula** to expertly maneuver through the vasculature and parenchyma and to change trajectory if necessary.
- **Directional handles** allow precise directional control of the stiffening cannula and needle/stylet to enter the portal vein and change direction when advancing through the parenchyma.



Directional Needle/Stylet

SCORPION® PORTAL VEIN ACCESS SERIES

Revolutionize your portal vein access procedures through *enhanced visibility* and *improved durability* combined with novel bi-directional steerability with our Scorpion Portal Vein Access Series.

SCORPION®
PORTAL VEIN ACCESS SET



0.040" Stylet

SCORPION^X®
PORTAL VEIN ACCESS SET



17Ga Needle

360° of bi-directional steerability
with curved stiffening cannula
and curved needle/stylet

Access Set Components

Scorpion® Access Set

Catalog # **TPS005**

14ga x 52cm Stiffening Cannula
7F x 52cm Cannula Sheath
0.040" x 73cm Curved Nitinol Stylet
5F x 71cm PEEK Stylet Catheter
5F x 80cm MPA Catheter
Spacer Clip

Scorpion® X Access Set

Catalog # **TPS006**

13ga x 52cm Stiffening Cannula
8F x 52cm Cannula Sheath
17ga x 73cm Curved Nitinol Needle
6.2F x 71cm PEEK Needle Catheter
5F x 80cm MPA Catheter
Spacer Clip

How does the Scorpion® Portal Vein Access Series work?

Scan QR below to view



Indication

The Scorpion® Portal Vein Access Sets are intended for transjugular liver access in diagnostic and interventional procedures.

Contact your local Argon Medical representative or customer service at customer.service@argonmedical to learn more about Scorpion®.

© 2021 Argon Medical Devices, Inc. All rights reserved. Argon Medical Devices, the stylized logo, and the product and service marks mentioned herein are trademarks or registered trademarks of Argon Medical Devices, Inc. in the United States and other countries.

Scorpion is a trademark of Hatch Medical, L.L.C. US pat. 8287481, 8628491, 9072596, EP2294177A1, All other trademarks are the property of their respective owners.

2600 Dallas Parkway, Suite 440, Frisco, TX 75034
P: 800.927.4669 F: 903.677.9396 argonmedical.com

MADE IN THE USA



96-6500-27/1021A