

Access with Ease


V-STICK™

VASCULAR ACCESS SETS


Now with Talwire for smooth access in peripheral vessels

V-Stick Vascular Access Sets are designed to facilitate smooth and easy introduction into the vessel for guidewires up to .038".

Access Needle




Coaxial Introducer



Talwire Guidewire



Talwire is the only articulating tip access wire designed to reduce the number of failed 1st attempts in radial, distal radial, pedal, and vein access procedures.^{1,2}



Backed by Clinical Data

Study #1

97

distal radial
procedures

65

radial access
procedures

88%

wires placed at
first attempt

Study #2

91

pedal access
procedures
attempted

97%

1st attempt
success rate
(88/91)

100%

all-attempt
success was
achieved

V-Stick Vascular Access Sets w/Coaxial Introducer (Box of 10)

Catalog #	Sheath	Needle			Dilator		Guidewire		
	Size	Gauge	Length	Echogenic	Stiffness	Locking Nut Color	Diameter	Length	Material
391504300E	4F	21G	7cm	Yes	Standard	Red	.018"	40cm	NT
391504400E	4F	21G	7cm	No	Standard	Red	.018"	40cm	NT
391505300E	5F	21G	7cm	Yes	Standard	Red	.018"	40cm	NT
391505400E	5F	21G	7cm	No	Standard	Red	.018"	40cm	NT
391604300E	4F	21G	7cm	Yes	Stiff	White	.018"	40cm	NT
391604400E	4F	21G	7cm	No	Stiff	White	.018"	40cm	NT
391605300E	5F	21G	7cm	Yes	Stiff	White	.018"	40cm	NT
391605400E	5F	21G	7cm	No	Stiff	White	.018"	40cm	NT

¹Tal, M. G., Covey, A., Qaqish, S., Livne, R., & Klass, D. (2023). Prospective evaluation of efficacy and safety of distal radial and radial artery access using a novel articulating-tip guidewire. *The Journal of Vascular Access*, 1–6. <https://doi.org/10.1177/11297298231212227>

²Nair, S., & Sastry, A. (2025). Retrograde tibial-pedal access can be safely and effectively obtained using a wire with an articulating tip. Sunrise Vascular Center.

Learn About the

Talwire

Difference



For Full Indication and Risk Information, refer to **ArgonMedical.com**

© 2026 Argon Medical Devices, Inc. All rights reserved. Argon Medical Devices, the stylized logo, and the product and service marks mentioned herein are trademarks or registered trademarks of Argon Medical Devices, Inc. in the United States and other countries.

7800 Dallas Parkway, Suite 200, Plano, TX 75024 P: 800.927.4669 F: 903.677.9396

96-6500-47/03042026 REV A **Rx Only**

